



UniCel[®] DxC
Synchron[®] Clinical Systems
Host Interface Specifications

This manual is intended for use with

UniCel[®] DxC 600

UniCel[®] DxC 800

UniCel[®] DxC 600i

WARRANTY AND DISCLAIMER NOTICE

BY ACCEPTING OR USING THIS DOCUMENT, THE USER EXPRESSLY AGREES TO ALL TERMS AND CONDITIONS CONTAINED IN THIS NOTICE. If the User rejects any or all of this notice, User shall return this document, unused, to Beckman Coulter for refund or credit of any amount invoiced User for this document. Beckman Coulter makes no warranty whatsoever with respect to the accuracy or completeness of this document and reserves the right to modify the interface covered by this document ("Interface") and/or this document at any time. User assumes full responsibility for verifying the correct operation of any computer system or device connected to the Interface and the accuracy and/or completeness of any data to or from the Interface. User assumes full liability for all costs associated with and all consequences of any kind resulting from errors or omissions in this document or changes which Beckman Coulter may make to this document and/or the Interface, including—but not limited to—reprogramming costs. BECKMAN COULTER EXPRESSLY DISCLAIMS ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, WHETHER EXPRESS OR IMPLIED, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, OR WHICH MIGHT OTHERWISE ARISE FROM COURSE OF DEALING BETWEEN THE PARTIES OR USAGE OR TRADE.

Table Of Contents

CHAPTER 1	Host Interface Specifications	1-1
	Overview	1-2
	Hardware Interface	1-6
	Data Transmission	1-7
	Unidirectional Protocol	1-8
	Bidirectional	1-11
CHAPTER 2	LX20/DxC Interface	2-1
	LX20/DxC Messages	2-2
	Stream 800 - Special Functions	2-4
	Stream 801 - Sample/Cup Program	2-6
	Stream 802 - Results	2-15
	Stream 803 - Instrument Status	2-29
	Stream 804 - Setup Status	2-34
CHAPTER 3	CX7 Compatible Interface	3-1
	CX7 Compatible Messages	3-2
	Stream 700 - Special Functions	3-4
	Stream 701 - Sample/Cup Program	3-6
	Stream 702 - Results	3-15
	Stream 703 - Instrument Status	3-29
	Stream 704 - Setup Status	3-34
CHAPTER 4	Tables/Codes	4-1
	Test Type Codes	4-2
	Sample Type Codes	4-3
	Chemistry Codes	4-4
	Unit Codes	4-10
	Instrument Codes	4-12
	Result Error Codes	4-18
	Special Calculation Status Codes	4-24
	Instrument Status Codes	4-25
	Access Codes Converted to DxC Codes for DxC 600i Reporting	4-26
	Key Code Conversion for Local Languages	4-27

INDEX

CHAPTER 1 Host Interface Specifications

Table of Contents

Host Interface Specifications	1-2
Overview	1-2
Introduction	1-2
UniCel LX20/DxC/CX7 Compatible	1-2
Definitions	1-4
Interface Protocol Options	1-5
Hardware Interface	1-6
Introduction	1-6
Connector Pin Assignments	1-6
Cable Length	1-6
Data Transmission	1-7
Introduction	1-7
Transmission Format	1-7
Unidirectional Protocol	1-8
Introduction	1-8
Transmitted Information	1-8
Unidirectional Software Control	1-9
Protocol Control Characters	1-9
Unidirectional Hardware Control	1-9
Unidirectional Transmission Example 1	1-10
Unidirectional Transmission Example 2	1-10
Bidirectional	1-11
Introduction	1-11
Bidirectional Protocol Control Characters	1-11
Line Bidding	1-11
Data Transfer	1-12
Successful Data Transfer	1-13
Unsolicited Messages	1-15
Unsolicited Message Transmission Example	1-17
Solicited Messages	1-17
Solicited Message Transmission Example	1-18
Downloading Sample Programming from Host to UniCel DxC	1-18
Host Downloading Sample Program Example	1-19

Host Interface Specifications

Overview

Introduction

These instructions to the Host Interface Specification provide the necessary information to interface the UniCel DxC System to a Laboratory Information System (LIS). Use this computer communications specification with software version 1.0 or higher. This manual has four chapters:

1. Host Interface Specifications
2. LX20/DxC Interface
3. CX7 Compatible Interface
4. Tables/Codes

Host Interface Specifications and Tables/Codes are general sections. *LX20/DxC Interface* and *CX7 Compatible Interface* are specific sections; refer to the section that describes the interface of your system.

UniCel LX20/DxC/CX7 Compatible

The following table describes the differences between the LX20/DxC Interface and the CX7 Compatible Interface on a UniCel DxC system. It also lists the results if you use the CX7 Compatible Interface.

Table 1.1 LX20/DxC/CX7 Compatible

Feature	LX20/DxC Interface Specification	CX7 Compatible Interface Specification
Patient ID	15 Alphanumeric	12 Alphanumeric <i>The sample result is not transmitted if the Patient ID is more than 12 characters.</i>
Sex	M, F, B <i>Converts M/F to B</i>	M, F <i>Converts M/F to M</i>
Patient Comment	45 Alphanumeric	25 Alphanumeric <i>If more than 25 characters are used, only the first 25 are sent.</i>
Sample ID	15 Alphanumeric	11 Alphanumeric <i>The Sample result is not transmitted if the Sample ID exceeds 11 characters.</i>

(1 of 2)

Table 1.1 LX20/DxC/CX7 Compatible, continued

Feature	LX20/DxC Interface Specification	CX7 Compatible Interface Specification
Control Lot Numbers	12 Characters	Not Supported
Sector/Rack Number	4 digits <i>0000 to 0999 can be accepted from the host (LIS). Up to 9999 can be sent from the instrument.</i>	2 digits <i>If the 3 digit check box is selected, 1 character field reserved for the cup data is used as a place to store the third digit of the rack number. E.g. rack 123, position 4 is represented as 23 in rack field and 14 in cup field.</i>
Cuvette Number	3 Numeric	2 Hexadecimal <i>The cuvette number is reported in hexadecimal notation. For example, cuvette 100 is reported as 64; cuvette 109 is reported as 6D.</i>
Sample Comment Code 1	25 Alphanumeric	25 Alphanumeric
Sample Comment Code 2	Not Supported	25 Alphanumeric <i>Not Used</i>
Dates	8 digits <i>Long date. Format is LDATE (ddmmyyyy)</i>	6 digits <i>Format is DATE (ddmmyy)</i>
Clear Rack/ Sample IDs	Stream 801 - Function 3 <i>Rack number field supports any rack to be cleared.</i>	Stream 701 - Function 3 <i>Samples in racks above rack 99 must be cleared with Sample ID, or manually cleared.</i>

(2 of 2)

Definitions

The following terms and definitions are useful in understanding this interface specification.

Table 1.2 Interface Specification Terms

Term	Definition
UniCel Dx C	UniCel Dx C 600 or UniCel Dx C 800
Dx C 600i or Dx C 800i	UniCel Dx C 600i or UniCel Dx C 800i
Dx C	UniCel Dx C 600 or UniCel Dx C 800 analyzer
Host	Clinical laboratory computer system
ALPHA	Alphanumeric
DATE	Date. The format is ddmm yy.
Download ^a	The process of the host sending sample/rack programming to the UniCel Dx C System.
Handshake ^a	The process of the UniCel Dx C System and the host communicating with each other. Permission to transmit is requested, permission is granted, and receipt of transmission is acknowledged (ACK – NAK protocol).
LDATE	Long Date. The format is ddmm yyyy.
Message	A collection of related information, in the form of records, on a single patient, sent at one time.
Null Component	A message component that contains no data. Null components are represented by two component delimiters with no character in between.
Null Field	A message field that contains no data. Null fields are represented by two field delimiters with no character in between.
Text	The set of displayable characters, from ASCII characters decimal 32 to decimal 126.

^a Bidirectional only

Interface Protocol Options

The UniCel DxC System provides four interface protocol options which are established using the Host Communications option in *Setup*. These protocol options are defined in the table below.

Table 1.3 Interface Protocol Options

Protocol	Description
Unidirectional	Information is transmitted from the UniCel DxC System to the host computer. The host uses flow control to control incoming data. The three flow control options are: <ul style="list-style-type: none">• Software Flow Control (XON/XOFF)• Hardware Flow Control (RTS/CTS) or• No Flow
Bidirectional	Information is transmitted from the UniCel DxC System to the host computer and from the host computer to the UniCel DxC System. The two systems use a software handshake to communicate with each other to request permission to send information and to acknowledge receipt of information (ACK-NAK protocol).
Bidirectional with Query	Follows the same protocol as bidirectional interface, with the added capability to request sample programs from the host. "Double host query" is not supported. (That is, the DxC queries the Datalink and the Datalink queries the LIS.) "Requery" allows the DxC to Query the Host when a sample is loaded and has a "Complete" status at the DxC.
Off	No information is sent or received by the UniCel DxC System.

Hardware Interface

Introduction

The hardware interface to the UniCel Dx C System external communication port uses a standard RS-232-C 9-pin or 25-pin D-connector (male). Communication is done with asynchronous byte transmission.

Connector Pin Assignments

Connector pin assignments are shown below.

Table 1.4 Connector Pin Assignments

25-Pin	9-Pin	Signal	Direction
1	-	Protected ground – used on one side only	
2	3	Transmit data	Out
3	2	Receive data	In
4 ^a	8 ^a	Clear to send (CTS)	In
5 ^b	7 ^b	Request to send (RTS)	Out
6 ^c	4 ^c	Data terminal ready (DTR)	Out
7	5	Signal ground	
20 ^c	6 ^c	Data set ready (DSR)	In
-	1	Data Carrier Detect (DCD)	In
-	9	Ring Indicator (RI)	

^a CTS is required before data will be transmitted.

^b RTS is turned on when operational.

^c Pins 6 and 20 and 4 and 6 are only needed when a modem is used.

Note: Pins 4 (CTS) and 5 (RTS) and 7 (RTS) and 8 (CTS) are enabled only in unidirectional mode.

th

Use a cable length of 40 feet (12 m) or less with the RS-232-C interface.

Data Transmission

Introduction

The UniCel DxC System uses 7-bit ASCII (provides the 128 character ASCII set when 7 data bits are selected). For local languages, 8-bit ASCII must be used (provides the 256-character ASCII set).

Transmission Format

The serial transmission format is configured from the UniCel DxC System console using the Host Communications option as outlined in CHAPTER 4, *System Setup* of the UniCel DxC Synchron Clinical Systems *Instructions for Use* manual and detailed in the UniCel DxC Synchron Clinical Systems *Reference Manual*. The following options are available.

Table 1.5 Serial Transmission Format

Parameter	Options	Default
Data Transmission Mode	Unidirectional, Bidirectional, Bidirectional with Query, Off	Unidirectional
Data Bits	7, 8 ^a	8
Stop Bits	1, 2	1
Parity	Even, Odd, None	None
Flow control	Software (XON/XOFF), Hardware (RTS/CTS), None	Software (XON/XOFF)
Device ID	00-99	00
Baud Rates	300, 1200, 2400, 4800, 9600, 19200	9600
Re-query	Selected, Not Selected	Not Selected

^a For use with non-English keyboards, select 8-bit ASCII for the host interface in order to send and receive special characters (for key code conversion information, refer to Figure 4.1 Key Code Conversion Chart).

Unidirectional Protocol

Introduction

In unidirectional transmission, information is transmitted in one direction only from the UniCel DxC System to the host computer. This information is transmitted as soon as it is available.

Transmitted Information

The information transmitted includes:

NOTICE

The 800 stream below is referenced to the LX20/DxC Interface. The CX7 compatible interface is similar in data streams and should be interpreted as a 700 stream where an 800 stream is displayed.

Table 1.6 Transmitted Information

Information	Description
Cup Header	The Cup Header information for each cup in the rack is transmitted when a cup has finished and results are ready to send to the host. (Stream 802 - Function 1)
Test Results	A Test Result message is transmitted as each test is completed. (Stream 802 - Functions 3, 11, 13, 81)
End of Cup	An End-of-Cup message is transmitted when all tests programmed for a sample are completed. (Stream 802 - Function 5)
End of Run	An End-of-Run message is transmitted when all tests programmed have been completed and the DxC has gone into the <i>Standby mode</i> . (Stream 803 - Function 17)
Power Up	The Power Up message is transmitted when the DxC is booted. (Stream 803 - Function 1)
Host Setup	The Host Setup message is transmitted when the DxC is booted and when there are changes in the Host Communications in System Setup. (Stream 800 - Function 2)

Output field width is dependent on the units and decimal precision selected. This should not exceed the fixed field width specified in this document.

Unidirectional Software Control

The flow of information from the UniCel DxC System can be controlled by the host through an XON/XOFF protocol. The sequence of events is:

- The host transmits the character XOFF prior to its buffer overflowing.
- The DxC immediately suspends transmission.
- When the host's buffer is ready to accept information again, the character XON is transmitted.
- The DxC then resumes transmission.
- The host should transmit the XON character after power up to indicate it is ready to receive data.

Protocol Control Characters

The software protocol control characters follow:

Table 1.7 Protocol Control Characters

Character	Name	ASCII	HEX
XON	Resume transmission	DL1	11
XOFF	Suspend transmission	DL3	13

Unidirectional Hardware Control

The flow of information from the UniCel DxC System can be controlled by the host through the CTS hardware signal. When the CTS signal is brought negative, the DxC will suspend transmission. The DxC will resume transmission when the CTS signal is brought positive. Examples of Unidirectional Transmission are cited below.

NOTICE

Do not select the Hardware option for flow control. With the Hardware option, ESD interference can affect the performance.

Unidirectional Transmission Example 1

Unidirectional mode is selected and the host is not using software control of transmission. The accession number is 168, the rack is 12, the position is 1, and the chemistries are Sodium (01A), Potassium (01B), and Chloride (04A). An operator-defined special calculation (SPC_CALC) is transmitted.

[illegible]

(- Indicates space holder)

Unidirectional Transmission Example 2

The same information is transmitted as in *Example 1*, but software flow control has been implemented by the host, selected by the DxC operator, and used by the Host where indicated.

UniCel DxC

HOST

```

===== XON

[-0,802,01,25091998,080534,--168,RE,--12,-1,RO,#####,SE,121-----,
-----,Harris-----,
Maria-----,J,100A10B100-----,Bond-----,25091998,0100,
ER-----,-47,5,26051951,B,Lipemic-----,
#####,####,####,#####,--3,01A-,01B-,04A-]C2<CR><LF>

<===== XOFF

<===== XON

[-0,802,03,25091998,080812,--168,-----116,--12,-1,121-----,01A-,###,
#####,###,-1,----104.7,#####,2,0,04,LO,NR,0,NA,104.65540,
-----,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,1.0000,
#####]0D<CR><LF>
[-0,802,03,25091998,080811,--168,-----116,--12,-1,121-----,01B-,###,
#####,###,-1,----2.45,#####,2,0,04,LO,NR,0,NA,2.4467065,-----,
NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,1.0000,
#####]13<CR><LF>
[-0,802,03,25091998,080812,--168,-----116,--12,-1,121-----,04A-,###,
#####,###,-1,----77.8,#####,2,0,04,LO,NR,0,NA,77.755951,-----,
NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,1.0000,
#####]FB<CR><LF>
[-0,802,11,25091998,080534,--168,--12,-1,121-----,-1,
SPC_CALC-----,OK,19.017143,-----]B9<CR><LF>
[-0,802,05,25091998,082242,--168,121-----,-12,-1]D5<CR><LF>

```

(- Indicates space holder)

Bidirectional

Introduction

In bidirectional transmission, information is transmitted from the UniCel DxC System to the host, and from the host to the DxC System. The two systems communicate with each other using a software handshake. Permission to send information is requested, permission is granted, and receipt of information is acknowledged (ACK - NAK protocol). Table 1.8 describes the bidirectional protocol control characters that are referred to throughout this section.

Bidirectional Protocol Control Characters

Bidirectional protocol control characters are listed below.

Table 1.8 Bidirectional Protocol Control Characters

Acronym	Name	ASCII	HEX
LB	Line Bid	SOH	01
LBO	Line Bid Override Request	STX	02
LF	Line Feed	LF	0A
ENQ	Enquiry (used to re-establish communications)	ENQ	05
EOT	End of Transmission	EOT	04
ACK-0	Even Acknowledgment	ACK	06
ACK-1	Odd Acknowledgment	ETX	03
NAK	Negative Acknowledgment	NAK	15

Line Bidding

Before either the UniCel DxC System or the host transmits a message, they must bid for the communication line. Line bidding uses the LB, LBO, EOT, ACK-0, and NAK characters.

Stage	Description
1	To bid for the line, <EOT> <LB> is transmitted when the line is idle.
2	The EOT clears the line.
3	The receiving system responds by transmitting ACK-0 to acknowledge the line bid and allow message transmission, or by transmitting NAK to deny the line bid.
4	If the line bid is denied, the originator waits a short period of time (e.g., one second) then bids for the line again.

- **LINE BID** – When the DxC is very busy, the response to a LINE BID may take up to 15 seconds.
1. When the DxC has some messages in the output queue, it does a LINE BID. When the LINE BID is accepted by the host, the instrument sends one set of messages at a time. It then gives up the line and waits for 2 seconds before continuing with the other messages.
 2. If the host has messages to transmit, it can do a LINE BID (EOT LB) even if the instrument is busy transmitting messages (after sending a message, the DxC looks for a NAK, ACK, or LB). When the DxC sends a message, the host can respond with a line bid <EOT><SOH>. The message in progress will be saved by the DxC. The line will be granted to the host by sending an <ACK>.

NOTICE

Do not activate the LINE BID while the instrument is transmitting messages because it would stop the collation of results.

- **LINE BID TIME OUT** occurs if the receiving system does not respond to the EOT LB within 15 seconds, or responds with something other than an ACK-0 or LBO. After 7 consecutive time outs, or unrecognizable responses (i.e., not ACK-0 or NAK), the originator considers the line DOWN. The originator waits 20 seconds and tries the line bid again.
- **CONTENTION** occurs when both systems bid for the line at the same time. The DxC will be considered the master and the host should respond with ACK-0. However, the host may override the DxC line bid by transmitting LBO in response. The DxC will respond by transmitting ACK-0 to acknowledge the line bid override and allow message transmission, or it will ignore the message if DxC was not bidding for the line when the message was sent.

Data Transfer

After successfully bidding for the line, the originator or sending system transmits its message. If the sending device does not transmit a message within 20 seconds, the receiving device times out. Further communication will require another line bid.

Successful Data Transfer

A successful data transfer consists of: [<MESSAGE>] <CS><CR><LF>.

- CS – Check sum. It is two hexadecimal digits.
256 - ((Sum of ASCII value of all characters in <MESSAGE>, including spaces, commas, ' ' and ' ') MODULE 256) = CS
- CR – Carriage return
- LF – Line feed

If the checksum is incorrect or if any element is missing, the data transfer is unsuccessful.

After data transmission, the receiving system acknowledges transmission as follows:

Stage	Description
1	If data transfer was successful, the receiving system alternately returns ACK-1 and ACK-0 after each message. Since the receiving system responded ACK-0 to the line bid, ACK-1 is the correct acknowledgment to the first message, ACK-0 to the second message, ACK-1 to the next message, and so forth. When data transfer is complete, the sending system transmits EOT.
2	If data transfer was unsuccessful due to a bad checksum or other problem, the receiving system responds NAK. The sending system re-transmits the message.
3	<p>If the sending system does not receive an acknowledgment within 15 seconds after data transfer, or if an incorrect acknowledgment is received, it transmits ENQ. The receiving system retransmits its last acknowledgment (e.g., ACK-0, ACK-1, or NAK). The sending system responds to that acknowledgment in one of the following manners:</p> <ul style="list-style-type: none"> • If the correct acknowledgment is transmitted, the sending system transmits the next message or EOT. • If an incorrect acknowledgment or NAK is transmitted, the sending system retransmits the last message.

(1 of 2)

Stage	Description
4	<p>TIME OUT occurs:</p> <ul style="list-style-type: none"> • If the sending system does not transmit data in 20 seconds after successfully bidding for the line, the receiving system times out and returns to idle. • If the sending system does not receive an acknowledgment within 15 seconds after data transfer, or if an incorrect acknowledgment is received, it then transmits ENQ. If an acknowledgment is not received after seven consecutive ENQ transmissions, the Dx C waits 20 seconds before it attempts a line bid and continues to do so until an acknowledgment is received.

(2 of 2)

Examples of data transfer (These examples are also correct if the labels UniCell Dx C and HOST are reversed):

NAK Example

UniCel Dx C (Sending System)

HOST (Receiving System)

EOT SOH	----->	
	<-----	ACK
Message 1	----->	
	<-----	ETX
Message 2	----->	
	<-----	ACK
Message 3	----->	
	<-----	NAK
Message 3	----->	
	<-----	ETX
EOT	----->	

No Response

UniCel Dx C (Sending System)

HOST (Receiving System)

EOT SOH	----->	
	<=====	ACK
Message 1	=====>	
	<=====	ETX
Message 2	=====>	
		No Response (15 secs)
ENQ	=====>	
	<=====	ETX
Message 2	=====>	
	<=====	ACK
EOT	=====>	

No Response Return to Idle

UniCel Dx C (Sending System)

HOST (Receiving System)

ACK	<=====	EOT SOH
	=====>	
IN IDLE STATE		No Response (>20 secs)
MESSAGE 1 IGNORED	<=====	MESSAGE 1
	<=====	
ACK	<=====	EOT SOH
	=====>	
	<=====	MESSAGE 1
	<=====	
ETX	<=====	EOT
	=====>	
	<=====	

Request for Idle State

UniCel Dx C (Sending System)

HOST (Receiving System)

EOT SOH	=====	
	<=====	ACK
Message 1	=====	
	<=====	ETX
Message 2	=====	
	<=====	ACK
EOT	=====	
IN IDLE STATE		IN IDLE STATE

Unsolicited Messages

Unsolicited messages are those messages which are automatically transmitted by the UniCel Dx C System as the information becomes available. The following unsolicited messages are transmitted in the bidirectional mode.

NOTICE

The 800 stream below is referenced to the LX20/Dx C Interface. The CX7 compatible interface is similar in data streams and should be interpreted as a 700 stream where an 800 stream is displayed.

Table 1.9 Unsolicited Messages

Message	Description
Cup Header	The Cup Header information for each cup in the rack is transmitted when a rack is loaded and accepted. (Stream 802 - Function 1)
Test Results	A Test Result message is transmitted as each test is completed. (Stream 802 - Functions 3, 11, 13, 81)
End of Cup	An End-of-Cup message is transmitted when all tests programmed for a sample are completed. (Stream 802 - Function 5)

(1 of 2)

Table 1.9 Unsolicited Messages, continued

Message	Description
End of Run	An End-of-Run message is transmitted when all tests programmed have been completed and the DxC has gone into the <i>Standby</i> mode. (Stream 803 - Function 17)
Power Up	The Power Up message is transmitted when the DxC is booted. (Stream 803 - Function 1) Note: Undefined characters may be transmitted during system boot-up.
Bidirectional Startup	The Bidirectional message is transmitted when the bidirectional or bidirectional with query interface option is enabled. (Stream 803 - Function 2)
Host Setup	The Host Setup message is transmitted when the DxC is booted and when the operator enters and exits Host Communications in System Setup. This is transmitted immediately after the interface is changed in the System Setup. (Stream 800 - Function 2)

(2 of 2)

Unsolicited Message Transmission Example

The following message is an example of an unsolicited message transmission.

UniCel DxC

HOST

```

EOT SOH  =====>
<===== ACK
[-0,802,01,25091998,080534,--168,RE,--12,-1,RO,#####,SE,121-----,
-----,-----,-----,-----,-----,-----,-----,-----,
J,100A10B100-----,Bond-----,25091998,0100,ER-----,-47,5,26051951,
B,Lipemic-----,#####,###,###,#####,-3,01A-,01B-,
04A-]C2<CR><LF>
<===== ETX
[-0,802,03,25091998,080812,--168,-----116,--12,-1,121-----,01A-,###,#####,###,
-1,-----104.7,#####2,0,04,LO,NR,0,NA,104.65540,-----,NO,NO,NO,NO,NO,NO,NO,NO,NO,
NO,NO,NO,NO,NO,NO,NO,NO,1.0000,#####]0D<CR><LF>
<===== ACK
[-0,802,03,25091998,080811,--168,-----116,--12,-1,121-----,01B-,###,#####,###,
-1,-----2.45,#####2,0,04,LO,NR,0,NA,2.4467065,-----,NO,NO,NO,NO,NO,NO,NO,NO,NO,
NO,NO,NO,NO,NO,NO,NO,NO,1.0000,#####]15<CR><LF>
<===== ETX
[-0,802,03,25091998,080812,--168,-----116,--12,-1,121-----,04A-,###,#####,###,
-1,-----77.8,#####2,0,04,LO,NR,0,NA,77.755951,-----,NO,NO,NO,NO,NO,NO,NO,NO,NO,
NO,NO,NO,NO,NO,NO,NO,NO,1.0000,#####]FB<CR><LF>
<===== ACK
[-0,802,11,25091998,080534,--168,--12,-1,121-----,-1,SPC_CALC-----,OK,
19.017143,-----]B9<CR><LF>
[-0,802,05,25091998,082242,--168,121-----,--12,-1]D5<CR><LF>
<===== ETX
[ 0,803,17,06101998,100137]4C<CR><LF>
<===== ACK
EOT =====>

```

(- Indicates space holder)

Solicited Messages

Solicited messages are transmitted by the host to request information from the UniCel DxC System. Solicited messages are used by the host to:

- Confirm that a bidirectional system is attached to the host serial port. (Stream 800 - Function 1)
- Obtain DxC instrument status information. (Stream 803 - Function 3)
- Obtain a list of the DxC installed chemistries. (Stream 804 - Functions 7)

Solicited Message Transmission Example

The following message is an example of a solicited message transmission.

UniCel DxC

HOST

	<=====	EOT SOH
ACK	=====>	
	<=====	[00,803,03]<CS><CR><LF>
ETX	=====>	
	<=====	EOT
EOT SOH	=====>	
	<=====	ACK
	[-0,803,04,07101998,111246,DxC-,02,15,#####]<CS><CR><LF>=>	
	<=====	ETX
EOT	=====>	

(<CS> Indicates checksum)

(- Indicates space holder)

Downloading Sample Programming from Host to UniCel DxC

- Sample programming may be transmitted directly from the host to the UniCel DxC Systems. The information that can be programmed includes:
 - Rack and cup number
 - Sample ID
 - Test type
 - Fluid type
 - Demographics
 - Chemistry requests (Stream 801 - Function 1)
- The host, after sending a single sample/cup must release the line by sending an EOT and then wait for the DxC to send the sample/cup return status. (Stream 801 - Function 2) This response by the DxC is used to notify the host whether or not the program was accepted (e.g., a BUSY response may indicate programming is being done at the DxC console). There is an interlock prohibiting simultaneous programming of a single sample from the host and the UniCel DxC System console.
- After the DxC releases the line (sends an EOT), the host may then repeat the process for as many sample/cups as necessary for completion of the download. If this protocol is not followed, improper chemistry requests may be downloaded to the DxC.
- In addition, the host can clear previous rack programming before transmitting the new sample programming. (Stream 801 - Function 3)

- In response to the clear rack message transmitted by the host, the DxC transmits an OK, BUSY, or SYNTAX ERROR message. (Stream 801 - Function 4) This response by the DxC is used to notify the host as to whether the rack was cleared (e.g., A BUSY response may indicate programming is being done at the DxC console or the rack is on the sample carousel.)
- There is an interlock prohibiting simultaneous clearing from the host and clearing or programming from the DxC System console.

Host Downloading Sample Program Example

The following is an example of a host downloading a sample program.

UniCel DxC	HOST
	EOT SOH
ACK	
[00,801,01,0401,01,0,RO,TU,RSG0001-----,-----,-----,Sample-has-been-hemolyzed,Rosenthel-----,Robert-----,R,PID456789012345,Dr.-H.-Schwartz---,08041998,1134,Humana-Hospital-Brea,-35,5,12041963,M,Chronic-asmatic-activity-27293133353739414345,111.1--,22.2,33.3,03,4.4444,003,01A-,0,01B-,0,02A-,0]<CS><CR><LF>	
ETX	
	EOT
EOT SOH	
	ACK
[-0,801,02,-0,--306,-401,-1,RSG0001-----]<CS><CR><LF>=>	
	ETX
EOT	

(<CS> Indicates checksum)
(- Indicates space holder)

Bidirectional

Host Downloading Sample Program Example

CHAPTER 2 LX20/DxC Interface

Table of Contents

LX20/DxC Interface	2-2
LX20/DxC Messages	2-2
Introduction	2-2
Message Format	2-2
Message Streams	2-3
Stream 800 - Special Functions	2-4
Stream 800 - Function 1: Are You There?	2-4
Stream 800 - Function 2: Host Setup	2-4
Stream 801 - Sample/Cup Program	2-6
Stream 801 - Function 1: Sample/Cup Program	2-6
Stream 801 - Function 2: Sample/Cup Return Status	2-10
Stream 801 - Function 3: Clear Rack/Sample IDs	2-12
Stream 801 - Function 4: Clear Rack/Sample IDs Status	2-13
Stream 801 - Function 6: Host Query Rack/Sample IDs	2-14
Stream 802 - Results	2-15
Stream 802 - Function 1: Cup Header	2-15
Stream 802 - Function 3: Test Results	2-17
Stream 802 - Function 5: End of Cup	2-20
Stream 802 - Function 11: Special Calculations Results	2-21
Stream 802 - Function 13: Timed Urine Results	2-22
Stream 802 - Function 81: Expanded Results	2-23
Stream 803 - Instrument Status	2-28
Stream 803 - Function 1: Power Up	2-28
Stream 803 - Function 2: Bidirectional On	2-30
Stream 803 - Function 3: Request Instrument State	2-31
Stream 803 - Function 4: Instrument State	2-31
Stream 803 - Function 17: End of Run	2-32
Stream 804 - Setup Status	2-33
Stream 804 - Function 7: Request Installed Chemistries	2-33
Stream 804 - Function 8: Installed Chemistries	2-33

LX20/DxC Interface

LX20/DxC Messages

Introduction

The data record format for all messages transmitted to and from the UniCel DxC System is:

' [' <MESSAGE> '] ' <CS><CR><LF>

This section describes the information contained in the <MESSAGE> portion.

Message Format

The format of the <MESSAGE> is:

<DEVICEID> , <STREAM> , <FUNCTION> , <FIELD1> , ..., <FIELDN> where:

- <DEVICE ID> is a number between 0 and 99. The default UniCel DxC ID is 0 but can be changed using the Host Communications Setup procedure. All messages which have a different Device ID than the Host Communications Screen will not be processed.
- <STREAM> is a number between 800 and 899.
- <FUNCTION> is a number between 1 and 99.
- <FIELD1>, ..., <FIELDN> are data fields associated with the <STREAM> and <FUNCTION>.
- All data fields <FIELDN> are fixed lengths.
- Numeric fields, excluding function numbers, are right justified and blank filled to the maximum length.
- Character and string fields are left justified and right blank filled to the maximum length.
- If numeric data exceeds the maximum field length, the field will be filled with asterisks (*). If a field does not apply in a record, it will be filled with pound signs (#).
- Alphanumeric fields received from the host cannot contain commas (.). Commas received in an alphanumeric field will generate an HCP error. A comma entered in an alphanumeric field at the instrument is transmitted as a semicolon (;). A semicolon received from the host is converted to a comma.
- Valid and Invalid characters for entries are listed in the table below:

Table 2.1 ASCII Codes

Type of Entry	THESE ARE VALID ASCII Character Codes...	EXCEPT FOR THESE INVALID Characters and Character Codes...
Sample ID	33 thru 126	A <i>space</i> and \$, ; * ? [] \ ^ &
Other Entries	32 thru 175	ASCII Character Codes 44 and 127

- **Alpha:** Any printable ASCII string (Commas are translated into semicolons when sent to the host. Semicolons are translated into commas when received from the host.)
- **LDATE:** Long Date. Format is ddmmyyyy.

Message Streams

Messages sent to and from the UniCel DxC System are divided into streams.

- Each stream corresponds to one group of related operations within the interface. Within each stream are one or more functions. For each stream only one function can be active at a time.
- Multiple streams can be active at the same time.

Stream transmission options are outlined in CHAPTER 4, *System Setup* of the UniCel DxC Synchron Clinical Systems *Instructions for Use* manual.

Stream 800 - Special Functions

Stream 800 - Function 1: Are You There?

Sent by: Host

Purpose: To make sure that a bidirectional system is attached to the host port. Upon receipt of the message 800-01, the UniCel DxC responds with the message 800-02, *Host Setup*.

Table 2.2 Are You There

Field	Length	Format	Notes
Device ID	2	0-99	Operator assigned
Stream	3	800	
Function	2	01	

Example:

[00,800,01]97<CR><LF>

Stream 800 - Function 2: Host Setup

Sent by: UniCel DxC

Purpose: Sent by the UniCel DxC in response to message 800-01. Also sent unsolicited upon power up and whenever there are changes to the Host Communications in *System Setup*.

Table 2.3 Host Setup

Field	Length	Format	Notes
Device ID	2	0-99	Operator assigned
Stream	3	800	
Function	2	02	
Date	8	LDate	
Time	6	Time	
Software Revision	10	ALPHA	Rev NNN.N
Driver Mode	1	B or U	B = Bidirectional U = Unidirectional
Stream 800 Function 2	1	0 or 1	0 = OFF 1 = ON
Stream 801 Function 6	1	0 or 1	0 = OFF 1 = ON
Stream 802 Function 3	1	0 or 1	0 = OFF 1 = ON

(1 of 2)

Table 2.3 Host Setup, continued

Field	Length	Format	Notes
Stream 802 Function 7	1	0 or 1	0 = OFF 1 = ON
Stream 802 Function 9	1	0 or 1	0 = OFF 1 = ON
Stream 802 Function 11	1	0 or 1	0 = OFF 1 = ON
Stream 802 Function 23	1	0 or 1	0 = OFF 1 = ON
Stream 802 Function 25	1	0 or 1	0 = OFF 1 = ON
Stream 802 Function 81	1	0 or 1	0 = OFF 1 = ON
Stream 802 Function 83	1	0 or 1	0 = OFF 1 = ON
Stream 803 Function 1	1	0 or 1	0 = OFF 1 = ON
Stream 803 Function 2	1	0 or 1	0 = OFF 1 = ON
Stream 803 Function 5	1	0 or 1	0 = OFF 1 = ON
Stream 803 Function 7	1	0 or 1	0 = OFF 1 = ON
Stream 803 Function 13	1	0 or 1	0 = OFF 1 = ON
Stream 803 Function 17	1	0 or 1	0 = OFF 1 = ON
Stream 802 Function 13	1	0 or 1	0 = OFF 1 = ON

(2 of 2)

Example:

```
[ -0,800,02,07101998,084523,1.05.48X--,U,1,0,1,0,0,1,0,0,1,0,1,1,0,0,0,1,1,
0,0,0,0,0,0,0]2A<CR><LF>
```

(- Indicates space holder)

Stream 801 - Sample/Cup Program

Stream 801 - Function 1: Sample/Cup Program

Sent by: Host

Purpose: To inform UniCel DxC what tests to perform for a given sample. The DxC responds with an 801-02, *Sample/Cup Return Status*, which shows whether or not the program was accepted.

Table 2.4 Sample/Cup Program

Field	Length	Format	Notes
Device ID	2	00-99	Operator assigned
Stream	3	801	
Function	2	01	
Rack Number	4	0000-0999	0 = No rack number assigned (Cup Number must also be 0). Rack numbers greater than 999 will result in a Rack Number Too Large return code. Refer to Note 3 after this table.
Cup Number	2	0-4	0 = No cup number assigned (Rack Number must also be 0). A cup number greater than 4 will result in a BAD MESSAGE return code. Refer to Note 3 after this table.
Update Flag	1	FLAG	0 = Replace all programming and patient demographics 1 = Add to existing program Refer to Note 2 and Note 4 after this table.
Test Type	2	TEST CODE	Refer to Table 4.1 Test Type Codes. Refer to Note 4 after this table.
Sample Type	2	FLUID CODE	Refer to Table 4.2 Sample Type Codes.

(1 of 3)

Table 2.4 Sample/Cup Program, continued

Field	Length	Format	Notes
Sample ID	15	ALPHA	Must be present in all modes. All lower case letters are converted to upper case. Refer to Note 4 after this table.
Control Name	20	ALPHA	Refer to Note 1 after this table.
QC Lot Number	12	ALPHA	Not used.
Sample Comment Code 1	25	ALPHA	
Name Last	18	ALPHA	
Name First	15	ALPHA	
Name Middle Initial	1	ALPHA	
Patient ID	15	ALPHA	
Doctor	18	ALPHA	
Draw Date	8	LDATE	Refer to Note 7 and Note 9 after this table.
Draw Time	4	TIME	
Location	20	ALPHA	
Age	3	NUMERIC	Default = 0
Age Units	1	NUMERIC	Default = 5 Refer to Table 4.3 Age Unit Codes.
Birth Date	8	LDATE	Refer to Note 9 after this table.
Sex	1	ALPHA	M, F, B
Patient Comments	45	ALPHA	
Timed Urine Volume	7	NUMERIC	0-99999.0
Timed Urine Period	4	NUMERIC	Collection period 0-99.0
Timed Urine Creatinine	4	NUMERIC	Serum creatinine 0-99.0

(2 of 3)

Table 2.4 Sample/Cup Program, continued

Field	Length	Format	Notes
Timed Urine Creatinine Units	2	UNIT CODE	Blank is only valid if the Timed Urine Creatinine field is zero or blank. Refer to Note 8 after this table.
Timed Urine Area	6	NUMERIC	Surface area 0 – 9.9900 Blank is OK.
Number Tests	3	NUMERIC	1-999
Test-1 Chem	4	CHEM CODE	Refer to Note 5 and Note 6 after this table.
Test-1 ORDAC	1	FLAG	1 = ORDAC
Test-2 Chem	4	CHEM CODE	Refer to Note 5 and Note 6 after this table.
Test-2 ORDAC	1	FLAG	1 = ORDAC
Test-N ^a Chem	4	CHEM CODE	Refer to Note 5 and Note 6 after this table.
Test-N ORDAC	1	FLAG	1 = ORDAC

(3 of 3)

^a N = Value from Number Tests Field

Note 1: The Control Name uniquely identifies the control of this result. Only the selected fluid type for the control can be used. Any other fluid type is a syntax error. All lower case letters are converted to upper case. If a control name is included in the message, the system will address the sample program as a control sample program.

Note 2: A new program will be entered if no program exists for the specified sample/cup. Duplicate chemistries will not be added to an existing program.

Updates consist of chemistry additions only and are allowed only under the following conditions:

- Sample Type and Program Type of program to be updated match Fluid Type and Test Type (respectively) of sample program sent by the host.
- Sample Status of the program to be updated is *Sample Required*, *In-Progress*, *Removed*, *Manual Assign*, or *Rerun*; and Control Name of the sample program sent by the host is blank for patient samples, or matches the Control Name of the existing control sample program.

Overwrites are allowed under the following conditions:

- Sample Status of the program to be overwritten is *Sample Required*, *Incomplete*, *Complete*, *Manual Assign*, or *Rerun*.

Note 3: Rack numbers greater than 999 will result in a Rack Number too Large (15) return code in the 801:2 message.

A cup number greater than 4 will result in a BAD MESSAGE return code in the 801:2 message.

If the Rack Number is 0, the Cup Number must also be 0.

If the Cup Number is 0, the Rack Number must also be 0.

Note 4: If the Test Type is ST (STAT) or RO (Routine), the sample program will be rejected if the Sample ID matches a defined control ID on the system.

Note 5: The system will accept control sample programs even if they contain chemistries not defined for the control on the instrument. It is assumed that the host will handle the extra chemistry results.

Note 6: If the host stream contains at least one of the Chem Codes used for serum index chemistries, the system will select to run all serum indices tests. The system will treat a serum index test selection as an exception to the requirement that chemistry selections in the host stream have to be configured on the instrument.

Note 7: If the host stream contains a Draw Time but does not contain a Draw Date, the system will clear the Draw Time to zero and log an event, indicating it has done so.

Note 8: In order for the Timed Urine Results Message to be sent back to the host, unit codes that are sent by the host must translate into one of the selectable CREM or CREA system (DxC) setup units.

The Timed Urine Results Message is *not* sent back to the host if Timed Urine Creatinine units sent by the host

- are blank
- do not translate into one of the selectable CREM or CREA system setup units or
- are not recognized by the instrument.

For this case, an event is logged which indicates that the data in the Timed Urine Creatinine Field has been cleared and that the Timed Urine Results will not be run.

Note 9: The DxC will reject a 2-digit year field. If the DxC receives a 2-digit field for the year, the DxC will log the event, indicating a problem was detected. The field will be cleared and process the sample with a blank date.

Example:

```
[00,801,01,0401,01,0,RO,TU,RSG0001-----,-----,
-----,Sample-has-been-hemolyzed,Rosenthel-----,
Robert-----,R,PID456789012345,Dr.-H.-Schwartz---,08041998,
1134,Humana-Hospital-Brea,-35,5,12041963,M,Chronic-asmatic-activity-
27293133353739414345,111.1--,22.2,33.3,03,4.4444,003,01A-,0,01B-,0,
02A-,0]<CS><CR><LF>
```

(- Indicates space holder)

Stream 801 - Function 2: Sample/Cup Return Status

Sent by: UniCel DxC

Purpose: Sent in response to an 801-01 message from the host. 801-02 indicates if the sample/cup program was accepted, or why it was not accepted.

Table 2.5 Sample/Cup Return Status

Field	Length	Format	Notes
Device ID	2	0-99	Operator assigned
Stream	3	801	
Function	2	02	
Return Code	2	00-99	0 = OK 1 = BAD MESSAGE – There is something wrong with the message syntax. 2 = BUSY – The cup is running or sample is being programmed. 3 = NOT CONFIGURED – Message contains at least one non –configured chemistry. 4 = NON ORDAC – Message contains at least one non-ORDAC chemistry, required for ORDAC. 5 = DILUTION ERROR – Message contains chemistries with different dilution factors. 6 = CONTROL NOT CONFIGURED 7 = CALIBRATOR Rack ONLY

(1 of 2)

Table 2.5 Sample/Cup Return Status, continued

Field	Length	Format	Notes
Return Code, continued	2	0-99	8 = Not used 9 = DxC ERROR – An internal logical error has occurred. 10 = Completed sample (Program update only) 11 = Incompatible fluid types. Used for program update or control sample program. 12 = Incompatible test types. (Program update only) 13 = Incompatible patient name. (Program update only) 14 = Sample ID matches existing Control ID. 15 = Rack Number too large. Rack Number must be between 0 and 1 and 0999. 16 = Too many tests programmed for sample. 17-20 = Future use
Accession Number	5	NUMERIC	1-65535, supplied by DxC
Rack Number	4	0000-0999	0 = No rack number assigned
Cup Number	2	00-4	
Sample ID	15	ALPHA	

(2 of 2)

Example:

[-0,801,02,-0,--306,-401,-1,RSG0001-----]<CS><CR><LF>

(- Indicates space holder)

Stream 801 - Function 3: Clear Rack/Sample IDs

Sent by: Host

Purpose: Sent by the host prior to programming samples to clear the program for a rack or a group of samples. It is not necessary for the host to send this message as each sample is cleared if the update flag is not set in the 801-01, Sample/Cup Program message.

- Receipt of this message causes the DxC to clear the program(s) for the rack or Sample IDs specified in the message.
- 801-03 messages that include both rack number and Sample IDs are rejected by sending an 801-04 message with a BAD MESSAGE status.
- DxC responds to the 801-03 message with an 801-04, Clear Rack/Sample IDs Status message indicating whether the clear command was accepted.
- If the DxC aspirated reagent for a rack or Sample ID, an 801-03 request to clear that rack or Sample ID receives an 801-04 message indicating the rack or Sample ID is busy.

Table 2.6 Clear Rack/Sample IDs

Field	Length	Format	Notes
Device ID	2	0-99	Operator assigned
Stream	3	801	
Function	2	03	
Rack Number	4	0000-0999	0 = No rack number assigned. Refer to Note 1 and Note 2 after this table.
Sample ID 1	15	ALPHA	Sample ID for a cup
Sample ID 2	15	ALPHA	Sample ID for a cup
Sample ID 3	15	ALPHA	Sample ID for a cup
Sample ID 4	15	ALPHA	Sample ID for a cup

Note 1: If one or more Sample IDs are specified, rack number must be zero. Specification of a non-zero rack number with non-blank Sample IDs will result in a 801-04 being returned with a rack return code of BAD MESSAGE.

Note 2: Rack numbers greater than 999 result in a rejected Clear Rack/Sample ID Message and cause the Clear Rack/Sample ID Status rack return code to be set to 5 - Rack Number too large.

Example:

[00,801,03,0000,SAMP1-----,SAMP2-----,SAMP3-----,SAMP4-----]6A<CR><LF>

(- Indicates space holder)

Stream 801 - Function 4: Clear Rack/Sample IDs Status**Sent by:** UniCel DxC**Purpose:** Sent in response to the 801-03 message. The host is informed if the clear command was accepted.

Table 2.7 Clear Rack/Sample IDs Status

Field	Length	Format	Notes
Device ID	2	0-99	Operator assigned
Stream	3	801	
Function	2	04	
Rack Number	4	0000-0999	0 = No rack number assigned
Rack Return Code	2	0-99	0 = Complete rack or all Sample IDs cleared. 1 = BAD MESSAGE 2 = BUSY The rack or one of the Sample IDs was running. 3 = DxC error 4 = Non-existent error 5= Rack Number too large
Sample ID 1 Return Code	2	0-99	0 = Cleared. 2 = BUSY - Not cleared 3 = DxC Error 4 = Non-existent error
Sample ID 2 Return Code	2	0-99	0 = Cleared 2 = BUSY - Not cleared 3 = DxC Error 4 = Non-existent error
Sample ID 3 Return Code	2	0-99	0 = Cleared 2 = BUSY - Not cleared 3 = DxC Error 4 = Non-existent error
Sample ID 4 Return Code	2	0-99	0 = Cleared 2 = BUSY - Not cleared 3 = DxC Error 4 = Non-existent error

Example:

[-0,801,04,-0,-0,-4,-4,-4,-4]6B<CR><LF>

(- Indicates space holder)

Stream 801 - Function 6: Host Query Rack/Sample IDs

Sent by: UniCel Dx-C

Purpose: To request the program for the Sample IDs, specified in the message from the host. This message is used in the bidirectional with query mode and finds samples for which it has no program. The host has 2.5 minutes to respond with first sample program.

Table 2.8 Host Query Rack/Sample IDs

Field	Length	Format	Notes
Device ID	2	0-99	Operator assigned
Stream	3	801	
Function	2	06	
Sample ID 1	15	ALPHA	Refer to Note 1 after this table.
Sample ID 2	15	ALPHA	Refer to Note 1 after this table.
Sample ID 3	15	ALPHA	Refer to Note 1 after this table.
Sample ID 4	15	ALPHA	Refer to Note 1 after this table.

Note 1: The message stream maintains the sample cup position of a rack. For example, if "samp 4" is in cup position 4, it will remain in the 7th field of the 801, 06 stream.

Example:

[-0,801,06,samp1-----,samp2-----,samp3-----,samp4 -----]EB<CR><LF>

(- Indicates space holder)

Stream 802 - Results

Stream 802 - Function 1: Cup Header

Sent by: UniCel DxC

Purpose: Sent before sending any test results because the header contains identification and demographic information about a sample.

Table 2.9 Cup Header

Field	Length	Format	Notes
Device ID	2	00-99	Operator assigned
Stream	3	802	
Function	2	01	
Start Date	8	LDATE	
Start Time	6	TIME	
Accession Number	5	NUMERIC	1-65535, supplied by DxC
Print Type	2	CODE	RG = Regular RE = Recall IN = Incomplete
Rack Number	4	1-0999	The rack number in which the sample was run.
Cup Number	2	1-4	The cup number in which the sample was run.
Test Type	2	TEST CODE	Refer to Table 4.1 Test Type Codes.
Future Use Space	9	ALPHA	
Sample Type	2	FLUID CODE	Refer to Table 4.2 Sample Type Codes.
Sample ID	15	ALPHA	
Control Name	20	ALPHA	Refer to Note 1 after this table.
Sample Comment Code 1	25	ALPHA	
Name Last	18	ALPHA	
Name First	15	ALPHA	
Name Middle Initial	1	ALPHA	
Patient ID	15	ALPHA	

(1 of 2)

Table 2.9 Cup Header, continued

Field	Length	Format	Notes
Doctor	18	ALPHA	
Draw Date	8	LDATE	
Draw Time	4	TIME	
Location	20	ALPHA	
Age	3	NUMERIC	Default = 0
Age Units	1	NUMERIC	Default = 5
Birth Date	8	LDATE	
Sex	1	ALPHA	M, F, B
Patient Comments	45	ALPHA	
Timed Urine Volume	7	REAL	0.0 - 99999.0 Blank is OK. Refer to Note 2 after this table.
Timed Urine Period	4	REAL	Collection period 0-99.0 Refer to Note 2 after this table.
Timed Urine Creatinine	4	REAL	Serum creatinine 0-99.0 Refer to Note 2 after this table.
Timed Urine Area	6	REAL	Surface area 1-9.9900 Refer to Note 2 after this table.
Number Chems Programmed	3	NUMERIC 000-999	Does not include replicates
Com Chem ID 1	4	CHEM CODE	
Com Chem ID 2	4	CHEM CODE	
Com Chem ID N	4	CHEM CODE	$N \leq 999$

(2 of 2)

Note 1: The Control Name uniquely identifies the control of this result.

Note 2: The number of decimal places is adjusted so the number can fit in this field.

Example:

```
[-0,802,01,25091998,080534,--168,RE,--12,-1,RO,#####,SE,121-----,
-----,Harris-----,
Maria-----,J,100A10B100-----,Bond-----,25091998,0100,
ER-----,-47,5,26051951,B,
Lipemic-----,#####,####,####,#####,
--3,01A-,01B-,04A-]C2<CR><LF>
```

(- Indicates space holder)

Stream 802 - Function 3: Test Results**Sent by:** UniCel DxC**Purpose:** Sent when all chemistry results for a cup have been completed and calculated. 802-03 can be turned off in the host setup option.

Table 2.10 Test Results

Field	Length	Format	Notes
Device ID	2	00-99	Operator assigned
Stream	3	802	
Function	2	03	
Completion Date	8	LDATE	
Completion Time	6	TIME	
Accession Number	5	NUMERIC	1-65535, supplied by DxC
Result Record Number	9	NUMERIC	1-10000
Rack Number	4	1-0999	
Cup Number	2	1-4	
Sample ID	15	ALPHA	
Com Chem ID	4	CHEM CODE	Refer to Table 4.4 and Table 4.5 for Chemistry (Chem) Codes.
Reagent Serial Number	3	ALPHA	
Reagent Lot Number	6	ALPHA	
Cuvette Number	3	NUMERIC	
Replicate Number	2	NUMERIC	
Results in Selected Units	9	ALPHA NUMERIC	Refer to Note 2 and Note 3 after this table.
Calibration Rate	9	REAL	Refer to Note 2 and Note 3 after this table.

(1 of 3)

Table 2.10 Test Results, continued

Field	Length	Format	Notes
Positive or Negative	1	CODE	0 = Negative 1 = Positive 2 = Not applicable 3 = Equivocal (DxC 600i or 800i only) 4 = Non-reactive (DxC 600i or 800i only) 5 = Reactive (DxC 600i or 800i only) 6 = Not confirmed (DxC 600i or 800i only) 7 = Confirmed (DxC 600i or 800i only) 8 = Gray zone flag (DxC 600i or 800i only) 9 = Access [®] 2 LOW flag (DxC 600i or 800i only) Refer to Note 4 after this table.
Suppress Value	1	FLAG	1 = Suppress value
Units	2	CODE	Refer to Table 4.6 Unit Codes.
Normal Range	2	CODE	Refer to range codes after this table.
Critical Range	2	CODE	Refer to range codes after this table.
ORDAC	1	FLAG	1 = ORDAC used
Control Range	2	FLAG	NA = Not applicable NR = Within 2 SD of mean H2 = 2 to 3 SD above mean H3 = Greater than 3 SD above mean H4 = Greater than 4 SD above mean L2 = 2 to 3 SD below mean L3 = Greater than 3 SD below mean L4 = Greater than 4 SD below mean IC = Incomplete CD = Control or chemistry deleted
Calculated Result	9	ALPHA NUMERIC	Refer to Note 1 and Note 3 after this table. Note: Result field reporting the results for TDM'S may contain a < (less than) or a > (greater than) sign.

(2 of 3)

Table 2.10 Test Results, continued

Field	Length	Format	Notes
Instrument Codes	9	ALPHA	Refer to Table 4.7 Instrument Codes with Remarks.
Result Error 1	2	ALPHA	Refer to Table 4.8 Result Error Codes.
Result Error 2	2	ALPHA	Refer to Table 4.8 Result Error Codes.
Result Error 3	2	ALPHA	Refer to Table 4.8 Result Error Codes.
Result Error 4	2	ALPHA	Refer to Table 4.8 Result Error Codes.
Result Error 5	2	ALPHA	Refer to Table 4.8 Result Error Codes.
Result Error 6	2	ALPHA	Refer to Table 4.8 Result Error Codes.
Result Error 7	2	ALPHA	Refer to Table 4.8 Result Error Codes.
Result Error 8	2	ALPHA	Refer to Table 4.8 Result Error Codes.
Result Error 9	2	ALPHA	Refer to Table 4.8 Result Error Codes.
Result Error 10	2	ALPHA	Refer to Table 4.8 Result Error Codes.
Result Error 11	2	ALPHA	Refer to Table 4.8 Result Error Codes.
Result Error 12	2	ALPHA	Refer to Table 4.8 Result Error Codes.
Result Error 13	2	ALPHA	Refer to Table 4.8 Result Error Codes.
Result Error 14	2	ALPHA	Refer to Table 4.8 Result Error Codes.
Result Error 15	2	ALPHA	Refer to Table 4.8 Result Error Codes.
Result Error 16	2	ALPHA	Refer to Table 4.8 Result Error Codes.
Dilution Factor	6	REAL	0.0 – 1000.0
Spare	25	ALPHA	Future use

(3 of 3)

RANGE CODES:

HI = above normal range

LO = below normal range

NR = within normal range

OR = out of range

NA = not applicable

Note 1: The number of decimal places is adjusted so the number can fit into this field.

Note 2: The number of decimal places in this field is the number in the selected precision for the selected units field in the chemistry database.

Note 3: If Suppress Value = 1, this field will be filled with #.

Note 4: For DxC chemistries, 0 and 1 are used only for drugs of abuse (DAT). On DxC 600i or DxC 800i, values from the range 0, 1, 3-9, are used by some Access 2 results, depending on the specific chemistry. Refer to the *Access Immunoassay System Manual* for additional information. All other chemistries use 2.

Example:

```
[ -0,802,03,25091998,080812,--168,-----116,--12,-1,121-----,  
01A-,###,#####,-1,----104.7,#####2,0,04,LO,NR,0,NA,104.65540,  
-----,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,1.0000,  
#####]0D<CR><LF>
```

(- Indicates space holder)

Stream 802 - Function 5: End of Cup**Sent by:** UniCel Dx-C**Purpose:** To indicate the end of a stream of test result messages for a specific cup.

Table 2.11 End of Cup

Field	Length	Format	Notes
Device ID	2	00-99	Operator assigned
Stream	3	802	
Function	2	05	
Date	8	LDATE	
Time	6	TIME	
Accession Number	5	NUMERIC	1-65535, supplied by Dx-C
Sample ID	15	ALPHA	
Rack number	4	1-0999	
Cup number	2	1-4	

Example:

[-0,802,05,25091998,082242,--168,121-----,--12,-1]D5<CR><LF>

(- Indicates space holder)

Stream 802 - Function 11: Special Calculations Results**Sent by:** UniCel Dx-C**Purpose:** Contains the result of a special calculation. This message contains accession number, rack number, cup number, and Sample ID. This message can be turned OFF in the Host Communication parameters screen in system setup.

Table 2.12 Special Calculations Results

Field	Length	Format	Notes
Device ID	2	00-99	Operator assigned
Stream	3	802	
Function	2	11	
Date	8	LDATE	
Time	6	TIME	
Accession Number	5	NUMERIC	1-65535, supplied by Dx-C
Rack Number	4	1-0999	
Cup Number	2	1-4	
Sample ID	15	ALPHA	

(1 of 2)

Table 2.12 Special Calculations Results, continued

Field	Length	Format	Notes
Replicate Number	2	NUMERIC	
Special Calc Name	20	ALPHA	Operator assigned
Special Calc Status	2	CODE	OK = Result is valid AB = One of the involved chemistries was programmed but not run. ZD = Denominator of a Ratio is zero
Special Calc Results	9	REAL	Refer to Note 1 after this table.
Special Unit String	8	ALPHA	

(2 of 2)

Note 1: The number of decimal places is adjusted to fit the number into this field.

Example:

```
[ -0,802,11,25091998,080534,--168,--12,-1,121-----,-1,
SPC_CALC-----,OK,50.367081,mmo1----]34<CR><LF>
```

(- Indicates space holder)

Stream 802 - Function 13: Timed Urine Results**Sent by:** UniCel DxC**Purpose:** Contains the result of a timed urine test. This message can be turned off on the host setup option.

Table 2.13 Timed Urine Results

Field	Length	Format	Notes
Device ID	2	0-99	Operator assigned
Stream	3	802	
Function	2	13	
Date	8	LDATE	
Time	6	TIME	
Accession Number	5	NUMERIC	1-65535, supplied by DxC
Rack number	4	1-0999	
Cup Number	2	1-4	
Sample ID	15	ALPHA	
Replicate Number	2	NUMERIC	
TU Calc Name	20	ALPHA	Operator assigned
TU Calc Status	2	CODE	OK = Result is valid AB = One of the involved chemistries was programmed but not run. ZD = Denominator of a Ratio is zero
TU Calc results	9	REAL	Refer to Note 1 after this table.
TU Unit String	8	ALPHA	

Note 1: The number of decimal places is adjusted to fit the number into this field.

Example:

```
[ -0,802,13,10041998,135429,---88,22,-1,S221-----,-1,NA-----,  
OK,244.00033,mmol/24.]6B<CR><LF>
```

(- Indicates space holder)

Stream 802 - Function 81: Expanded Results

Sent by: UniCel DxC

Purpose: Contains the expanded result of a test. This message can be turned off in the host setup option. Do not run the system with expanded results turned on unless directed by a Beckman Coulter representative.

Table 2.14 Expanded Results

Field	Length	Format	Notes
Device ID	2	0-99	Operator assigned
Stream	3	802	
Function	2	81	
Date	8	LDATE	
Time	6	TIME	
Accession Number	5	NUMERIC	1-65535, supplied by DxC
Sample Record Number	9	NUMERIC	
Result Record Number	9	NUMERIC	
Sample ID	15	ALPHA	
Rack Number	4	1-0999	
Cup Number	2	1-4	
Reagent Record Number	9	NUMERIC	
Replicate Number	2	NUMERIC	
Com Chem ID	4	CHEM CODE	Refer to Table 4.4 Chemistry (Chem) Codes.
Cuvette Number	3	NUMERIC	
Cal Level	1	0-6	Not used
Suppressed Answer	1	FLAG	1 = Answer suppressed
ORDAC Result	1	FLAG	1 = ORDAC
Noise Maximum	5	NUMERIC	
Outlier Maximum	5	NUMERIC	
Noise Threshold Total	3	NUMERIC	
Outlier Threshold Total	3	NUMERIC	

(1 of 4)

Table 2.14 Expanded Results, continued

Field	Length	Format	Notes
Initial Absorbance	9	REAL	Refer to Note 1 after this table.
Final Absorbance	9	REAL	Refer to Note 1 after this table.
Water Blank-Rate	9	REAL	Refer to Note 1 after this table.
Water Blank-Abs	9	REAL	Refer to Note 1 after this table.
Water Blank-Mean Deviation	9	REAL	Refer to Note 1 after this table.
Water Blank-Maximum Deviation	9	REAL	Refer to Note 1 after this table.
Blank-Rate	9	REAL	Refer to Note 1 after this table.
Blank-Abs	9	REAL	Refer to Note 1 after this table.
Blank-Mean Deviation	9	REAL	Refer to Note 1 after this table.
Blank-Maximum Deviation	9	REAL	Refer to Note 1 after this table.
Reaction-Rate	9	REAL	Refer to Note 1 after this table.
Reaction-Abs	9	REAL	Refer to Note 1 after this table.
Reaction-Mean Deviation	9	REAL	Refer to Note 1 after this table.
Reaction-Maximum Deviation	9	REAL	Refer to Note 1 after this table.
Raw Result	9	ALPHA	Refer to Note 1 after this table. Note: Result field reporting the results for TDM's may contain a < (less than) or > (greater than) sign.
Calculated Result	9	ALPHA	Refer to Note 1 after this table. Note: Result field reporting the results for TDM's may contain a < (less than) or > (greater than) sign.
Calibration Rate	9	REAL	Refer to Note 1 after this table.

(2 of 4)

Table 2.14 Expanded Results, continued

Field	Length	Format	Notes
Positive or Negative	1	FLAG	0 = Negative 1 = Positive 2 = Not applicable 3 = Equivocal (DxC 600i or DxC 800i only) 4 = Non-reactive (DxC 600i or DxC 800i only) 5 = Reactive (DxC 600i or DxC 800i only) 6 = Not confirmed (DxC 600i or DxC 800i only) 7 = Confirmed (DxC 600i or DxC 800i only) 8 = Gray zone flag (DxC 600i or DxC 800i only) 9 = Access 2 LOW flag (DxC 600i or DxC 800i only) Refer to Note 2 after this table.
Units	2	CODE	Refer to Table 4.6 Unit Codes.
Instrument Codes	9	ALPHA	Refer to Table 4.7 Instrument Codes with Remarks.
Result Error 1	2	CODE	Refer to Table 4.8 Result Error Codes.
Result Error 2	2	CODE	Refer to Table 4.8 Result Error Codes.
Result Error 3	2	CODE	Refer to Table 4.8 Result Error Codes.
Result Error 4	2	CODE	Refer to Table 4.8 Result Error Codes.
Result Error 5	2	CODE	Refer to Table 4.8 Result Error Codes.
Result Error 6	2	CODE	Refer to Table 4.8 Result Error Codes.
Result Error 7	2	CODE	Refer to Table 4.8 Result Error Codes.
Result Error 8	2	CODE	Refer to Table 4.8 Result Error Codes.

(3 of 4)

Table 2.14 Expanded Results, continued

Field	Length	Format	Notes
Result Error 9	2	CODE	Refer to Table 4.8 Result Error Codes.
Result Error 10	2	CODE	Refer to Table 4.8 Result Error Codes.
Result Error 11	2	CODE	Refer to Table 4.8 Result Error Codes.
Result Error 12	2	CODE	Refer to Table 4.8 Result Error Codes.
Result Error 13	2	CODE	Refer to Table 4.8 Result Error Codes.
Result Error 14	2	CODE	Refer to Table 4.8 Result Error Codes.
Result Error 15	2	CODE	Refer to Table 4.8 Result Error Codes.
Result Error 16	2	CODE	Refer to Table 4.8 Result Error Codes.
Intermediate ADC 1	5	NUMERIC	
Intermediate ADC 2	5	NUMERIC	
Intermediate ADC 3	5	NUMERIC	
Intermediate ADC 4	5	NUMERIC	
Intermediate ADC 5	5	NUMERIC	
Dilution Factor	6	REAL	0.0-1000.0
Future Use Space	23	ALPHA	

(4 of 4)

Note 1: The number of decimal places is adjusted to fit the number into this field.

Note 2: For DxC chemistries, 0 and 1 are used only for drugs of abuse (DAT). On UniCel DxC 600i or UniCel DxC 800i, values from the range 0, 1, 3-9, are used by some Access 2 results, depending on the specific chemistry. Refer to the *Access Immunoassay System Manual* for additional information. All other chemistries use 2.

Example:

[-0,802,81,25091998,080812,--168,#####,-----116,121-----,--12,
-1,385875973,-1,01A-,###,#,0,0,#####,#####,###,###,#####,#####,
#####,#####,#####,#####,#####,#####,#####,
#####,#####,#####,#####,#####,#####,-3882.800,104.65540,
#####,2,04,-----,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,
--538,-3344,-3882,---6,---3,1.0000,#####]71<CR><LF>

(- Indicates space holder)

Stream 803 - Instrument Status

Stream 803 - Function 1: Power Up

Sent by: UniCel DxC

Purpose: Sent upon DxC power up or reset. This message can be turned off in the host setup option.

Table 2.15 Power Up

Field	Length	Format	Notes
Device ID	2	0-99	Operator assigned
Stream	3	803	
Function	2	01	
Date	8	LDATE	
Time	6	TIME	
Software Set	10	ALPHA	Console Master Program Version Number
CV ICC BOOT	10	ALPHA	
CV ICC APP	10	ALPHA	
CV MSMC1 BOOT 1	10	ALPHA	
CV MSMC1 BOOT 2	10	ALPHA	
CV MSMC1 APP	10	ALPHA	
CV MSMC2 BOOT 1	10	ALPHA	
CV MSMC2 BOOT 2	10	ALPHA	
CV MSMC2 APP	10	ALPHA	
CV SM BOOT	10	ALPHA	
CV SM APP 21	10	ALPHA	
CV SM APP 11	10	ALPHA	
CV SM APP 12	10	ALPHA	
CV SM APP 13	10	ALPHA	
CV SM APP 14	10	ALPHA	
CV SM APP 31	10	ALPHA	
CV SM APP 41	10	ALPHA	

(1 of 3)

Table 2.15 Power Up, continued

Field	Length	Format	Notes
CV SM APP 42	10	ALPHA	
CV SM APP 43	10	ALPHA	
CV SM APP 51	10	ALPHA	
CV SM APP 52	10	ALPHA	
CV SM APP 61	10	ALPHA	
CV SM APP 62	10	ALPHA	
CV SM APP 63	10	ALPHA	
CV SM APP 71	10	ALPHA	
CV SM APP 81	10	ALPHA	
CV SM APP 82	10	ALPHA	
CV SM APP 91	10	ALPHA	
CV SM APP A1	10	ALPHA	
CV SM APP A2	10	ALPHA	
CV SM APP B1	10	ALPHA	
CV SM APP B2	10	ALPHA	
CV SM APP B3	10	ALPHA	
CV SM APP C1	10	ALPHA	
CV SM APP C2	10	ALPHA	
CV SM APP C3	10	ALPHA	
SPARE	10	ALPHA	All "#"s
SPARE	10	ALPHA	All "#"s
SPARE	10	ALPHA	All "#"s
SPARE	10	ALPHA	All "#"s
SPARE	10	ALPHA	All "#"s
SPARE	10	ALPHA	All "#"s
SPARE	10	ALPHA	All "#"s
SPARE	10	ALPHA	All "#"s
SPARE	10	ALPHA	All "#"s
SPARE	10	ALPHA	All "#"s
SPARE	10	ALPHA	All "#"s

(2 of 3)

Table 2.15 Power Up, continued

Field	Length	Format	Notes
SPARE	10	ALPHA	All "#"s

(3 of 3)

Example:

```
[ -0,803,01,10041998[-0,703,01,100497,164320, v01.00----,v00.00----,
v00.00----,v00.00----,v00.00----,v00.00----,v00.00----,v00.00----,
v00.00----,v01.06----,v00.13----,v04.54----,v04.54----,v04.54----,
v04.54----,v01.30----,v00.89----,v00.89----,v00.89----,v00.89----,
v00.89----,v04.54----,v04.54----,v04.54----,v00.89----,v00.1B----,
#####,v00.30----,v00.89----,v00.89----,v00.89----,v00.89----,
v00.89----,v00.89----,v00.89----,v00.89----,#####,#####,
#####,#####,#####,#####,#####,#####,#####,
#####,#####,#####,#####] 5E<CR><LF>
```

(- Indicates space holder)

Stream 803 - Function 2: Bidirectional On**Sent by:** UniCel Dx C**Purpose:** This message is sent to the host when bidirectional communication is enabled. This message can be disabled in the host setup option.

Table 2.16 Bidirectional On

Field	Length	Format	Notes
Device ID	2	0-99	Operator assigned
Stream	3	803	
Function	2	02	
Date	8	LDATE	
Time	6	TIME	
Software Set	10	ALPHA	DxC Software Version

Example:

```
[ -0,803,02,10041998,203332,A4.4.54--- ] F0<CR><LF>
```

(- Indicates space holder)

Stream 803 - Function 3: Request Instrument State

Sent by: Host

Purpose: Sent to obtain the state of the UniCel Dx C.

Table 2.17 Request Instrument State

Field	Length	Format	Notes
Device ID	2	0-99	Operator assigned
Stream	3	803	
Function	2	03	

Example:

[-0 , 803 , 03] 92 <CR> <LF>

(- Indicates space holder)

Stream 803 - Function 4: Instrument State

Sent by: UniCel Dx C

Purpose: Sent in response to a Request Instrument State, 803-03 message from the host.

Table 2.18 Instrument State

Field	Length	Format	Notes
Device ID	2	0-99	Operator assigned
Stream	3	803	
Function	2	04	
Date	8	LDATE	
Time	6	TIME	
Configuration	5	ALPHA	Dx C
MC State	2	CODE	Refer to Table 4.11 Instrument Status Codes.
CC State	2	CODE	Refer to Table 4.11 Instrument Status Codes.
Future Use	5	ALPHA	

Example:

[-0 , 803 , 04 , 11041998 , 105150 , Dx C - , 02 , 15 , ## , -1] 85 <CR> <LF>

(- Indicates space holder)

Stream 803 - Function 17: End of Run

Sent by: UniCel DxC

Purpose: Sent when the DxC completes all programmed tests. This message can be turned off in the host setup option.

Table 2.19 End of Run

Field	Length	Format	Notes
Device ID	2	0-99	Operator assigned
Stream	3	803	
Function	2	17	
Date Queued	8	LDATE	
Time Queued	6	TIME	

Example:

[-0,803,17,06101998,100137]4C<CR><LF>

(- Indicates space holder)

Stream 804 - Setup Status

Stream 804 - Function 7: Request Installed Chemistries

Sent by: Host

Purpose: Sent to obtain the list of installed chemistries from the DxC database.

Table 2.20 Request Installed Chemistries

Field	Length	Format	Notes
Device ID	2	0-99	Operator assigned
Stream	3	804	
Function	2	07	

Example:

[00,804,07]8F<CR><LF>

Stream 804 - Function 8: Installed Chemistries

Sent by: UniCel DxC

Purpose: This message is sent in response to a Request Installed Chemistries, 804-07 message from the host.

Table 2.21 Installed Chemistries

Field	Length	Format	Notes
Device ID	2	0-99	Operator assigned
Stream	3	804	
Function	2	08	
Date	8	LDATE	
Time	6	TIME	
Com Chem ID 1	4	CHEM CODE	Refer to Table 4.4 and Table 4.5 for Chemistry (Chem) Codes.
Com Chem ID N	4	CHEM CODE	Refer to Table 4.4 and Table 4.5 for Chemistry (Chem) Codes. N is incremented for each Com Chem ID that is configured.

Note: Only chemistry codes of chemistries that are configured on the DxC will be included in the message. The last chemistry code transmitted will be 0000 to indicate that no more chemistry codes follow. If no chemistries are configured on the instrument, only 0000 chemistry code will be listed.

Example:

```
[ -0, 804, 08, 11041998, 105518, 01A-, 01B-, 04A-, 02A-, 09D-, 08D-, 05D-, 03E-,  
06D-, 43D-, 07D-, 14B-, 56B-, 35A-, 31A-, 31B-, 54A-, 84A-, 10A-, 24C-, 24B-, 93B-,  
30A-, 30B-, 85A-, 86A-, 05A-, 98A-, 59A-, 44A-, 32A-, 32B-, 40A-, 87A-, 89A-, 12A-,  
62A-, 46B-, 72C-, 36A-, 06A-, 14A-, 83A-, 79B-, 52B-, 51B-, 53B-, 55A-, 33A-, 50A-,  
50B-, 07B-, 64A-, 63A-, 48A-, 92A-, 66A-, 67C-, 68C-, 65A-, 93C-, 94A-, 91A-, 11A-,  
42B-, 42C-, 88A-, 88B-, 41A-, 69C-, 70C-, 71B-, 0000] BE<CR><LF>
```

(- Indicates space holder)

Stream 804 - Setup Status

Stream 804 - Function 8: Installed Chemistries

CHAPTER 3 CX7 Compatible Interface

Table of Contents

CX7 Compatible Interface	3-2
CX7 Compatible Messages	3-2
Introduction	3-2
Message Format	3-2
Message Streams	3-3
Stream 700 - Special Functions	3-4
Stream 700 - Function 1: Are You There?	3-4
Stream 700 - Function 2: Host Setup	3-4
Stream 701 - Sample/Cup Program	3-6
Stream 701 - Function 1: Sample/Cup Program	3-6
Stream 701 - Function 2: Sample/Cup Return Status	3-9
Stream 701 - Function 3: Clear Rack/Sample IDs	3-11
Stream 701 - Function 4: Clear Rack/Sample IDs Status	3-13
Stream 701 - Function 6: Host Query Rack/Sample IDs	3-14
Stream 702 - Results	3-15
Stream 702 - Function 1: Cup Header	3-15
Stream 702 - Function 3: Test Results	3-17
Stream 702 - Function 5: End of Cup	3-21
Stream 702 - Function 11: Special Calculations Results	3-22
Stream 702 - Function 13: Timed Urine Results	3-23
Stream 702 - Function 81: Expanded Results	3-24
Stream 703 - Instrument Status	3-29
Stream 703 - Function 1: Power Up	3-29
Stream 703 - Function 2: Bidirectional On	3-31
Stream 703 - Function 3: Request Instrument State	3-31
Stream 703 - Function 4: Instrument State	3-32
Stream 703 - Function 17: End of Run	3-33
Stream 704 - Setup Status	3-34
Stream 704 - Function 7: Request Installed Chemistries	3-34
Stream 704 - Function 8: Installed Chemistries	3-34

CX7 Compatible Interface

CX7 Compatible Messages

Introduction

The data record format for all messages transmitted to and from the UniCel DxC System is:

['<MESSAGE>']<CS><CR><LF>

This section describes the information contained in the <MESSAGE> portion.

Message Format

The format of the <MESSAGE> is:

<DEVICEID>,<STREAM>,<FUNCTION>,<FIELD1>,...,<FIELDN> where:

- <DEVICE ID> is a number between 0 and 99. The default UniCel DxC ID is 0 but can be changed using the Host Communications Setup procedure. All messages which have a different Device ID than the Host Communications Screen will not be processed.
- <STREAM> is a number between 700 and 799.
- <FUNCTION> is a number between 1 and 99.
- <FIELD1>,...,<FIELDN> are data fields associated with the <STREAM> and <FUNCTION>.
- All data fields <FIELDN> are fixed lengths.
- Numeric fields, excluding function numbers, are right justified and blank filled to the maximum length.
- Character and string fields are left justified and right blank filled to the maximum length.
- If numeric data exceeds the maximum field length, the field will be filled with asterisks (*). If a field does not apply in a record, it will be filled with pound signs (#).
- Alphanumeric fields received from the host cannot contain commas (,). Commas received in an alphanumeric field will generate an HCP error. A comma entered in an alphanumeric field at the instrument is transmitted as a semicolon (;). A semicolon received from the host is converted to a comma.
- Valid and Invalid characters for entries are listed in the table below:

Table 3.1 ASCII Codes

Type of Entry	THESE ARE VALID ASCII Character Codes...	EXCEPT FOR THESE INVALID Characters and Character Codes...
Sample ID	33 thru 126	A <i>space</i> and \$, ; * ? [] \ ^ &
Other Entries	32 thru 175	ASCII Character Codes 44 and 127

- **Alpha:** Any printable ASCII string (Commas are translated into semicolons when sent to the host. Semicolons are translated into commas when received from the host.)

Message Streams

Messages sent to and from the UniCel Dx C System are divided into streams.

- Each stream corresponds to one group of related operations within the interface. Within each stream are one or more functions. For each stream only one function can be active at a time.
- Multiple streams can be active at the same time.

Stream transmission options are outlined in CHAPTER 4, *System Setup*, of the UniCel Dx C Synchron Clinical Systems *Instructions for Use* manual.

Stream 700 - Special Functions

Stream 700 - Function 1: Are You There?

Sent by: Host

Purpose: To make sure that a bidirectional system is attached to the host port. Upon receipt of the message 700-01, the UniCel DxC responds with the message 700-02, *Host Setup*.

Table 3.2 Are You There

Field	Length	Format	Notes
Device ID	2	0-99	Operator assigned
Stream	3	800	
Function	2	01	

Example:

[00,700,01]98<CR><LF>

Stream 700 - Function 2: Host Setup

Sent by: UniCel DxC

Purpose: Sent by the UniCel DxC in response to message 700-01. Also sent unsolicited upon power up and whenever there are changes to the Host Communications in *System Setup*.

Table 3.3 Host Setup

Field	Length	Format	Notes
Device ID	2	0-99	Operator assigned
Stream	3	700	
Function	2	02	
Date	6	Date	
Time	6	Time	
Software Revision	10	ALPHA	Rev NNN.N
Driver Mode	1	B or U	B = Bidirectional U = Unidirectional
Stream 700 Function 2	1	0 or 1	0 = OFF 1 = ON
Stream 701 Function 6	1	0 or 1	0 = OFF 1 = ON
Stream 702 Function 3	1	0 or 1	0 = OFF 1 = ON

(1 of 2)

Field	Length	Format	Notes
Stream 702 Function 7	1	0 or 1	0 = OFF 1 = ON
Stream 702 Function 9	1	0 or 1	0 = OFF 1 = ON
Stream 702 Function 11	1	0 or 1	0 = OFF 1 = ON
Stream 702 Function 23	1	0 or 1	0 = OFF 1 = ON
Stream 702 Function 25	1	0 or 1	0 = OFF 1 = ON
Stream 702 Function 81	1	0 or 1	0 = OFF 1 = ON
Stream 702 Function 83	1	0 or 1	0 = OFF 1 = ON
Stream 703 Function 1	1	0 or 1	0 = OFF 1 = ON
Stream 703 Function 2	1	0 or 1	0 = OFF 1 = ON
Stream 703 Function 5	1	0 or 1	0 = OFF 1 = ON
Stream 703 Function 7	1	0 or 1	0 = OFF 1 = ON
Stream 703 Function 13	1	0 or 1	0 = OFF 1 = ON
Stream 703 Function 17	1	0 or 1	0 = OFF 1 = ON
Stream 702 Function 13	1	0 or 1	0 = OFF 1 = ON

(- Indicates space holder)

Stream 701 - Sample/Cup Program

Stream 701 - Function 1: Sample/Cup Program

Sent by: Host

Purpose: To inform UniCel DxC what tests to perform for a given sample. The DxC responds with a 701-02, *Sample/Cup Return Status*, that shows whether or not the program was accepted.

Table 3.4 Sample/Cup Program

Field	Length	Format	Notes
Device ID	2	00-99	Operator assigned.
Stream	3	701	
Function	2	01	
Rack Number	2	0000-0999	0 = No rack number in bar code mode. (Cup Number must also be 0) Refer to Note 3 after this table.
Cup Number	2	0-4	0 = No cup number in bar code mode. A cup number greater than 4 will result in a BAD MESSAGE return code in the 701-02 message. Refer to Note 3 after this table.
Update Flag	1	FLAG	0 = Replace all programming and patient demographics. 1 = Add to existing program. Refer to Note 2 and Note 4 after this table.
Test Type	2	TEST CODE	Refer to Table 4.1 Test Type Codes. Refer to Note 4 after this table.
Sample Type	2	FLUID CODE	Refer to Table 4.2 Sample Type Codes.

(1 of 3)

Table 3.4 Sample/Cup Program, continued

Field	Length	Format	Notes
Sample ID	11	ALPHA	Must be present in all modes. All lower case letters are converted to upper case. Refer to Note 4 after this table.
Control Name	20	ALPHA	Refer to Note 1 after this table.
Sample Comment Code 1	25	ALPHA	
Sample Comment Code 2	25	ALPHA	Not Used by Dx C
Name Last	18	ALPHA	
Name First	15	ALPHA	
Name Middle Initial	1	ALPHA	
Patient ID	12	ALPHA	
Doctor	18	ALPHA	
Draw Date	6	DATE	Refer to Note 7 after this table.
Draw Time	4	TIME	
Location	20	ALPHA	
Age	3	NUMERIC	Default = 0
Age Units	1	NUMERIC	Default = 5 Refer to Table 4.3 Age Unit Codes.
Birth Date	6	DATE	
Sex	1	ALPHA	M, F, B
Patient Comments	25	ALPHA	
Timed Urine Volume	7	NUMERIC	0-99999.0
Timed Urine Period	4	NUMERIC	Collection period 0-99.0
Timed Urine Creatinine	4	NUMERIC	Serum creatinine 0-99.0
Timed Urine Area	6	NUMERIC	Surface area 1 – 9.9900

(2 of 3)

Table 3.4 Sample/Cup Program, continued

Field	Length	Format	Notes
Number Tests	3	NUMERIC	1-999
Test-1 Chem	4	CHEM CODE	Refer to Note 5 and Note 6 after this table.
Test-1 ORDAC	1	FLAG	1 = ORDAC
Test-2 Chem	4	CHEM CODE	Refer to Note 5 and Note 6 after this table.
Test-2 ORDAC	1	FLAG	1 = ORDAC
Test-N ^a Chem	4	CHEM CODE	Refer to Note 5 and Note 6 after this table.
Test-N ORDAC	1	FLAG	1 = ORDAC

(3 of 3)

^a N = Value from Number Tests Field

Note 1: The Control Name uniquely identifies the control of this result. Only the selected fluid type for the control can be used. Any other fluid type is a syntax error. All lower case letters are converted to upper case. If a control name is included in the message, the system will address the sample program as a control sample program.

Note 2: A new program will be entered if no program exists for the specified sample/cup. Duplicate chemistries will not be added to an existing program.

Updates consist of chemistry additions only.

Updates are allowed only under the following conditions:

- Sample Type and Program Type of program to be updated match Fluid Type and Test Type (respectively) of sample program sent by the host.
- Sample Status of the program to be updated is *Sample Required*, *In-Progress*, *Removed*, *Manual Assign*, or *Rerun*; and Control Name of the sample program sent by the host is blank for patient samples, or matches the Control Name of the existing control sample program.

Overwrites are allowed under the following conditions:

- Sample Status of the program to be overwritten is *Sample Required*, *Incomplete*, *Complete*, *Manual Assign*, or *Rerun*.

Note 3: The most significant character of the range of the Rack Number field will be located in the most significant character position of the Cup Number field.

Example: Rack 345, Cup 1 would be 45, 31.

Note 4: If the Test Type is ST (STAT) or RO(Routine), the sample program will be rejected if the Sample ID matches a defined control ID on the system.

Note 5: The system will accept control sample programs even if they contain chemistries not defined for the control on the instrument. It is assumed that the host will handle the extra chemistry results.

Note 6: If the host stream contains at least one of the Chem Codes used for serum index chemistries, the system will select to run all serum indices tests. The system will treat a serum index test selection as an exception to the requirement that chemistry selections in the host stream have to be configured on the instrument.

Note 7 : If the host stream contains a Draw Time but does not contain a Draw Date, the system will clear the Draw Time to zero and log an event, indicating it has done so.

Example:

```
[00,701,01,01,03,1,CO,SE,samp3-----,CONTROL_NAME_____,  
-----,SAMPLE_COMMENT_TWO_____,-----,  
PATIENT_F_NAME_,M,-----,-----,-----,  
-----,000,4,-----,M,-----,  
----,----,-----,001,01B-,0]2F<CR><LF>
```

(- Indicates space holder)

Stream 701 - Function 2: Sample/Cup Return Status

Sent by: UniCel Dx C

Purpose: Sent in response to a 701-01 message from the host. 701-02 indicates if the sample/cup program was accepted, or why it was not accepted.

Table 3.5 Sample/Cup Return Status

Field	Length	Format	Notes
Device ID	2	0-99	Operator assigned
Stream7	3	701	
Function	2	02	

(1 of 3)

Table 3.5 Sample/Cup Return Status, continued

Field	Length	Format	Notes
Return Code	2	00-99	<p>0 = OK</p> <p>1 = BAD MESSAGE – There is something wrong with the message syntax.</p> <p>2 = BUSY – The cup is running or sample is being programmed.</p> <p>3 = NOT CONFIGURED – Message contains at least one non – configured chemistry.</p> <p>4 = NON ORDAC – Message contains at least one non-ORDAC chemistry, required for ORDAC.</p> <p>5 = DILUTION ERROR – Message contains chemistries with different dilution factors.</p> <p>6 = CONTROL NOT CONFIGURED</p> <p>7 = CALIBRATOR Rack ONLY</p> <p>8 = Not used</p> <p>9 = DxC ERROR – An internal logical error has occurred.</p> <p>10 = Completed sample (Program update only)</p> <p>11 = Incompatible fluid types. Used for program update or control sample program.</p> <p>12 = Incompatible test types. (Program update only)</p> <p>13 = Incompatible patient name. (Program update only)</p> <p>14 = Patient Sample Program Sample ID matches existing DxC Control ID.</p> <p>15 = Not used by CX7 I/F messages.</p> <p>16 = Too many tests programmed for sample.</p> <p>17-20 = Future use.</p> <p>Refer to Note 1 and Note 2 after this table.</p>
Accession Number	5	NUMERIC	1-65535, supplied by DxC
Rack Number	2	0000-0999	<p>0 = No rack number assigned.</p> <p>Refer to Note 1 and Note 2 after this table.</p>

(2 of 3)

Table 3.5 Sample/Cup Return Status, continued

Field	Length	Format	Notes
Cup Number	2	0-4	Refer to Note 1 and Note 2 after this table.
Sample ID	11	ALPHA	

(3 of 3)

Note 1: The cup number and rack number must be filled in when the instrument is in non-bar code mode.

Note 2: The most significant character of the Rack Number field's range will be located in the most significant character position of the Cup Number field.

Example: Rack 345, Cup 1 would be 45, 31.

Example:

```
[-0,701,02,-0,--230,-1,-3,samp3-----]5D<CR><LF>
```

(- Indicates space holder)

Stream 701 - Function 3: Clear Rack/Sample IDs

Sent by: Host

Purpose: Sent by the host prior to programming samples to clear the program for a rack or a group of samples. It is not necessary for the host to send this message as each sample is cleared if the update flag is not set in the 701-01, Sample/Cup Program message.

- Receipt of this message causes the DxC to clear the program(s) for the rack or Sample IDs specified in the message.
- 701-03 messages that include both rack number and Sample IDs are rejected by sending a 701-04 message with a BAD MESSAGE status.
- DxC responds to the 701-03 message with a 701-04, Clear Rack/Sample IDs Status message indicating whether the clear command was accepted.
- If the DxC aspirated reagent for a rack or Sample ID, a 701-03 request to clear that rack or Sample ID receives a 701-04 message indicating the rack or Sample ID is busy.

Table 3.6 Clear Rack/Sample IDs

Field	Length	Format	Notes
Device ID	2	0-99	Operator assigned
Stream	3	701	
Function	2	03	

(1 of 2)

Table 3.6 Clear Rack/Sample IDs, continued

Field	Length	Format	Notes
Rack Number	2	0000-0999	0 = No rack number in bar code mode Refer to Note 1 and Note 3 after this table.
Sample ID 1	11	ALPHA	Sample ID for a cup Refer to Note 2 after this table.
Sample ID 2	11	ALPHA	Sample ID for a cup Refer to Note 2 after this table.
Sample ID N	11	ALPHA	Sample ID for a cup. N = 3 thru 7.

(2 of 2)

Note 1: The most significant character of the Rack Number field's range will be located in the most significant character position of the Cup Number field.

Example: Rack 345, Cup 1 would be 45, 31.

Example:

[-0, 701, 02, -0, --230, -1, -3, samp3-----] 5D<CR><LF>

(- Indicates space holder)

Note 2: If one or more Sample IDs are specified, rack number must be zero. Specification of a non-zero rack number with non-blank Sample IDs will result in a 701-04 being returned with a rack return code of BAD MESSAGE.

Note 3: Only rack numbers less than or equal to 99 may be cleared by rack number. All other racks must be cleared by Sample ID.

Example:

[00, 701, 03, 00, SAMP1-----, SAMP2-----, SAMP3-----, SAMP4-----, SAMP5-----, SAMP6-----, SAMP7-----] D2<CR><LF>

(- Indicates space holder)

Stream 701 - Function 4: Clear Rack/Sample IDs Status**Sent by:** UniCel DxC**Purpose:** Sent in response to the 701-03 message. The host is informed if the clear command was accepted.

Table 3.7 Clear Rack/Sample IDs Status

Field	Length	Format	Notes
Device ID	2	0-99	Operator assigned
Stream	3	701	
Function	2	04	
Rack Number	2	0000-0999	0 = No rack number in bar code mode Refer to Note 1 after this table.
Rack Return Code	2	0-99	0 = Complete rack or all Sample IDs cleared. 1 = BAD MESSAGE 2 = BUSY The rack or one of the Sample IDs was running. 3 = DxC error 4 = Non-existent error
Sample ID 1 Return Code	2	0-99	0 = Cleared. 2 = BUSY - Not cleared 3 = DxC Error 4 = Non-existent error
Sample ID 2 Return Code	2	0-99	0 = Cleared 2 = BUSY - Not cleared 3 = DxC Error 4 = Non-existent error
Sample ID N Return Code	2	0-99	0 = Cleared 2 = BUSY - Not cleared 3 = DxC Error 4 = Non-existent error N = 3 thru 7.

Note 1: The most significant character of the Rack Number field's range will be located in the most significant character position of the Cup Number field.

Example: Rack 345, Cup 1 would be 45, 31.

Example:

```
[ -0, 701, 02, -0, --230, -1, -3, samp3----- ] 5D<CR><LF>
```

Example:

```
[ -0, 701, 04, -0, -4, -4, -4, -4, -4, -4, -4, -4 ] 28 <CR> <LF>
```

(- Indicates space holder)

Stream 701 - Function 6: Host Query Rack/Sample IDs

Sent by: UniCel Dx-C

Purpose: To request the program for the Sample IDs, specified in the message from the host. This message is used in the bidirectional with query mode and finds samples for which it has no program. The host has 2.5 minutes to respond with first sample program.

Table 3.8 Host Query Rack/Sample IDs

Field	Length	Format	Notes
Device ID	2	0-99	Operator assigned
Stream	3	701	
Function	2	06	
Sample ID 1	11	ALPHA	Refer to Note 1 after this table.
Sample ID 2	11	ALPHA	Refer to Note 1 after this table.
Sample ID 3	11	ALPHA	Refer to Note 1 after this table.
Sample ID 4	11	ALPHA	Refer to Note 1 after this table.

Note 1: The message stream maintains the sample cup position of a rack. For example, if "samp 4" is in cup position 4, it will remain in the 7th field of the 701,06 stream.

Example:

```
[ -0, 701, 06, samp1-----, samp2-----, samp3-----, samp4-----, samp5-----,
samp6-----, samp7----- ] EB <CR> <LF>
```

(- Indicates space holder)

Stream 702 - Results

Stream 702 - Function 1: Cup Header

Sent by: UniCel DxC

Purpose: Sent before sending any test results because the header contains identification and demographic information about a sample.

Table 3.9 Cup Header

Field	Length	Format	Notes
Device ID	2	0-99	Operator assigned
Stream	3	702	
Function	2	01	
Start Date	6	DATE	
Start Time	6	TIME	
Accession Number	5	NUMERIC	1-65535, supplied by DxC
Print Type	2	CODE	RG = Regular RE = Recall IN = Incomplete
Rack Number	2	1-0999	The rack number in which the sample was run Refer to Note 1 after this table.
Cup Number	2	1-4	The cup number in which the sample was run Refer to Note 1 after this table.
Test Type	2	TEST CODE	Refer to Table 4.1 Test Type Codes.
Future Use Space	9	ALPHA	
Sample Type	2	FLUID CODE	Refer to Table 4.2 Sample Type Codes.
Sample ID	11	ALPHA	
Control Name	20	ALPHA	
Sample Comment Code 1	25	ALPHA	
Sample Comment Code 2	25	ALPHA	
Name Last	18	ALPHA	

(1 of 2)

Table 3.9 Cup Header, continued

Field	Length	Format	Notes
Name First	15	ALPHA	
Name Middle Initial	1	ALPHA	
Patient ID	12	ALPHA	
Doctor	18	ALPHA	
Draw Date	6	DATE	
Draw Time	4	TIME	
Location	20	ALPHA	
Age	3	NUMERIC	Default = 0
Age Units	1	NUMERIC	Default = 5
Birth Date	6	DATE	
Sex	1	ALPHA	M, F
Patient Comments	25	ALPHA	
Timed Urine Volume	7	REAL	0.0 - 99999.0 mL Refer to Note 2 after this table.
Timed Urine Period	4	REAL	Collection period 0-99.0 Refer to Note 2 after this table.
Timed Urine Creatinine	4	REAL	Serum creatinine 0-99.0 Refer to Note 2 after this table.
Timed Urine Area	6	REAL	Surface area 1-9.9900 Refer to Note 2 after this table.
Number Chems Programmed	3	NUMERIC 001-999	1-999 Does not include replicates
Com Chem ID 1	4	CHEM CODE	
Com Chem ID 2	4	CHEM CODE	
Com Chem ID N	4	CHEM CODE	N ≤ 999

(2 of 2)

Note 1: The most significant character of the Rack Number field will be located in the most significant character position of the Cup Number field.

Example: Rack 345, Cup 1 would be 45,31.

Note 2: The number of decimal places is adjusted so the number can fit in this field.

Example:

```
[ -0,702,01,100497,135429,---88,RG,22,-1,RO,#####,TU,S221-----,
-----,-----,-----,
White-----,Allen-----,B,S221-----,Nelson-----,
100497,1500,-----,-17,5,190280,M,-----,
1400.00,24.0,2.30,1.7300,--3,01A-,01B-,04A-]3A<CR><LF>
```

(- Indicates space holder)

Stream 702 - Function 3: Test Results

Sent by: UniCel Dx C

Purpose: Sent when all chemistry results for a cup have been completed and calculated. 702-03 can be turned off in the host setup option.

Table 3.10 Test Results

Field	Length	Format	Notes
Device ID	2	00-99	Operator assigned
Stream	3	702	
Function	2	03	
Completion Date	6	DATE	
Completion Time	6	TIME	
Accession Number	5	NUMERIC	1-65535, supplied by Dx C
Result Record Number	9	NUMERIC	1-10000
Rack Number	2	1-0999	Refer to Note 5 after this table.
Cup Number	2	1-4	Refer to Note 5 after this table.
Sample ID	11	ALPHA	
Com Chem ID	4	CHEM CODE	Refer to Table 4.4 and Table 4.5 for Chemistry (Chem) Codes.
Reagent Serial Number	3	ALPHA	
Reagent Lot Number	6	ALPHA	
Cuvette Number	2	HEX	

(1 of 4)

Table 3.10 Test Results, continued

Field	Length	Format	Notes
Replicate Number	2	NUMERIC	
Results in Selected Units	9	ALPHA NUMERIC	Refer to Note 2 and Note 3 after this table. Note: Result field reporting the results for TDMs may contain a < (less than) or a > (greater than) sign.
Calibration Rate	9	REAL	Refer to Note 2 and Note 3 after this table.
Positive or Negative	1	CODE	0 = Negative 1 = Positive 2 = Not applicable 3 = Equivocal (DxC 600i or DxC 800i) 4 = Non-reactive (DxC 600i or DxC 800i) 5 = Reactive (DxC 600i or DxC 800i) 6 = Not confirmed (DxC 600i or DxC 800i) 7 = Confirmed (DxC 600i or DxC 800i) 8 = Gray zone flag (DxC 600i or DxC 800i) 9 = Access 2 LOW flag (DxC 600i or DxC 800i) Refer to Note 4 after this table.
Suppress Value	1	FLAG	1 = Suppress value
Units	2	CODE	Refer to Table 4.6 Unit Codes.
Normal Range	2	CODE	Refer to range codes after this table.
Instrument Range	2	CODE	Refer to range codes after this table.
Critical Range	2	CODE	Refer to range codes after this table.
ORDAC	1	FLAG	1 = ORDAC used

(2 of 4)

Table 3.10 Test Results, continued

Field	Length	Format	Notes
Control Range	2	FLAG	NA = Not applicable NR = Within 2 SD of mean H2 = 2 to 3 SD above mean H3 = Greater than 3 SD above mean H4 = Greater than 4 SD above mean L2 = 2 to 3 SD below mean L3 = Greater than 3 SD below mean L4 = Greater than 4 SD below mean IC = Incomplete CD = Control or chemistry deleted
Calculated Result	9	ALPHA NUMERIC	Refer to Note 1 and Note 3 after this table. Note: Result field reporting the results for TDM'S may contain a < (less than) or a > (greater than) sign.
Instrument Codes	9	ALPHA	Refer to Table 4.7 Instrument Codes with Remarks.
Result Error 1	2	ALPHA	Refer to Table 4.8 Result Error Codes.
Result Error 2	2	ALPHA	Refer to Table 4.8 Result Error Codes.
Result Error 3	2	ALPHA	Refer to Table 4.8 Result Error Codes.
Result Error 4	2	ALPHA	Refer to Table 4.8 Result Error Codes.
Result Error 5	2	ALPHA	Refer to Table 4.8 Result Error Codes.
Result Error 6	2	ALPHA	Refer to Table 4.8 Result Error Codes.
Result Error 7	2	ALPHA	Refer to Table 4.8 Result Error Codes.
Result Error 8	2	ALPHA	Refer to Table 4.8 Result Error Codes.
Result Error 9	2	ALPHA	Refer to Table 4.8 Result Error Codes.
Result Error 10	2	ALPHA	Refer to Table 4.8 Result Error Codes.
Result Error 11	2	ALPHA	Refer to Table 4.8 Result Error Codes.
Result Error 12	2	ALPHA	Refer to Table 4.8 Result Error Codes.
Result Error 13	2	ALPHA	Refer to Table 4.8 Result Error Codes.
Result Error 14	2	ALPHA	Refer to Table 4.8 Result Error Codes.
Result Error 15	2	ALPHA	Refer to Table 4.8 Result Error Codes.
Result Error 16	2	ALPHA	Refer to Table 4.8 Result Error Codes.
Dilution Factor	6	REAL	0.0 – 1000.0

(3 of 4)

Table 3.10 Test Results, continued

Field	Length	Format	Notes
Reportable Range	2	CODE	Refer to range codes after this table.
Spare	20	ALPHA	Future use

(4 of 4)

RANGE CODES:

HI = above normal range

LO = below normal range

NR = within normal range

OR = out of range

NA = not applicable

Note 1: The number of decimal places is adjusted so the number can fit into this field.

Note 2: The number of decimal places in this field is the number in the selected precision for the selected units field in the chemistry database.

Note 3: If Suppress Value = 1, this field will be filled with #.

Note 4: For DxC chemistries, 0 and 1 are used only for drugs of abuse (DAT). On UniCel DxC 600i or DxC 800i, values from the range 0, 1, 3-9, are used by some Access 2 results, depending on the specific chemistry. Refer to the *Access Immunoassay System Manual* for additional information. All other chemistries use 2.

Note 5: The most significant character of the range of the Rack Number field is located in the most significant character position of the Cup Number field.

Example: Rack 345, Cup 1 would be 45, 31.

Example:

[illegible]

(- Indicates space holder)

Stream 702 - Function 5: End of Cup**Sent by:** UniCel DxC**Purpose:** To indicate the end of a stream of test result messages for a specific cup.

Table 3.11 End of Cup

Field	Length	Format	Notes
Device ID	2	00-99	Operator assigned
Stream	3	702	
Function	2	05	
Date	6	DATE	
Time	6	TIME	
Accession Number	5	NUMERIC	1-65535, supplied by DxC
Sample ID	11	ALPHA	
Rack number	2	1-0999	Refer to Note 1 after this table.
Cup number	2	1-4	Refer to Note 1 after this table.

Note 1: The most significant character of the Rack Number field's range will be located in the most significant character position of the Cup Number field.

Example: Rack 345, Cup 1 would be 45, 31.

Example:

[-0,702,05,100497,140240,---88,S221-----,22,-1]ED<CR><LF>

(- Indicates space holder)

Stream 702 - Function 11: Special Calculations Results

Sent by: UniCel Dx-C

Purpose: Contains the result of a special calculation. This message contains accession number, rack number, cup number, and Sample ID. This message can be turned OFF in the Host Communication parameters screen in system setup.

Table 3.12 Special Calculations Results

Field	Length	Format	Notes
Device ID	2	00-99	Operator assigned
Stream	3	702	
Function	2	11	
Date	6	DATE	
Time	6	TIME	
Accession Number	5	NUMERIC	1-65535, supplied by Dx-C
Rack Number	2	1-0999	Refer to Note 2 after this table.
Cup Number	2	1-4	Refer to Note 2 after this table.
Sample ID	11	ALPHA	
Replicate Number	2	NUMERIC	
Special Calc Name	20	ALPHA	Operator assigned
Special Calc Status	2	CODE	OK = Result is valid AB = One of the involved chemistries was programmed but not run. ZD = Denominator of a Ratio is zero
Special Calc Results	9	REAL	Refer to Note 1 after this table.
Special Unit String	8	ALPHA	

Note 1: The number of decimal places is adjusted to fit the number into this field.

Note 2: The most significant character of the range of the Rack Number field will be located in the most significant character position of the Cup Number field.

Example: Rack 345, Cup 1 would be 45, 31.

Example:

```
[ -0, 702, 11, 100497, 135429, ---88, 22, -1, S221-----, -1, SPC_CALC-----,
OK, 50.367081, mmol----] B5<CR><LF>
```

(- Indicates space holder)

Stream 702 - Function 13: Timed Urine Results**Sent by:** UniCel DxC**Purpose:** Contains the result of a timed urine test. This message can be turned off on the host setup option.

Table 3.13 Timed Urine Results

Field	Length	Format	Notes
Device ID	2	00-99	Operator assigned
Stream	3	702	
Function	2	13	
Date	6	DATE	
Time	6	TIME	
Accession Number	5	NUMERIC	1-65535, supplied by DxC
Rack number	2	1-0999	Refer to Note 2 after this table.
Cup Number	2	1-4	Refer to Note 2 after this table.
Sample ID	11	ALPHA	
Replicate Number	2	NUMERIC	
TU Calc Name	20	ALPHA	Operator assigned
TU Calc Status	2	CODE	OK = Result is valid AB = One of the involved chemistries was programmed but not run. ZD = Denominator of a Ratio is zero Refer to Table 4.10.
TU Calc results	9	REAL	Refer to Note 1 after this table.
TU Unit String	8	ALPHA	

Note 1: The number of decimal places is adjusted to fit the number into this field.

Note 2: The most significant character of the range of the Rack Number field will be located in the most significant character position of the Cup Number field.

Example: Rack 345, Cup 1 would be 45, 31

Example:

```
[ -0,702,13,100497,135429,---88,22,-1,S221-----,-1,NA-----,
OK,244.00033,mmol/24.]0B<CR><LF>
```

(- Indicates space holder)

Stream 702 - Function 81: Expanded Results

Sent by: UniCel Dx-C

Purpose: Contains the expanded result of a test. This message can be turned off in the host setup option. Do not run the system with expanded results turned on unless directed by a Beckman Coulter representative.

Table 3.14 Expanded Results

Field	Length	Format	Notes
Device ID	2	0-99	Operator assigned
Stream	3	702	
Function	2	81	
Date	6	DATE	
Time	6	TIME	
Accession Number	5	NUMERIC	1-65535, supplied by Dx-C
Sample Record Number	9	NUMERIC	
Result Record Number	9	NUMERIC	
Sample ID	11	ALPHA	
Rack Number	2	1-0999	Refer to Note 2 after this table.
Cup Number	2	1-4	Refer to Note 2 after this table.
Reagent Record Number	9	NUMERIC	
Replicate Number	2	NUMERIC	
Com Chem ID	4	CHEM CODE	Refer to Table 4.4 and Table 4.5 for Chemistry (Chem) Codes.
Cuvette Number	2	HEX	
Cal Level	1	0-6	Not used
Suppressed Answer	1	FLAG	1 = Answer suppressed
ORDAC Result	1	FLAG	1 = ORDAC
Noise Maximum	5	NUMERIC	
Outlier Maximum	5	NUMERIC	
Noise Threshold Total	3	NUMERIC	
Outlier Threshold Total	3	NUMERIC	
Initial Absorbance	9	REAL	Refer to Note 1 after this table.

(1 of 4)

Table 3.14 Expanded Results, continued

Field	Length	Format	Notes
Final Absorbance	9	REAL	Refer to Note 1 after this table.
Water Blank-Rate	9	REAL	Refer to Note 1 after this table.
Water Blank-Abs	9	REAL	Refer to Note 1 after this table.
Water Blank-Mean Deviation	9	REAL	Refer to Note 1 after this table.
Water Blank-Maximum Deviation	9	REAL	Refer to Note 1 after this table.
Blank-Rate	9	REAL	Refer to Note 1 after this table.
Blank-Abs	9	REAL	Refer to Note 1 after this table.
Blank-Mean Deviation	9	REAL	Refer to Note 1 after this table.
Blank-Maximum Deviation	9	REAL	Refer to Note 1 after this table.
Reaction-Rate	9	REAL	Refer to Note 1 after this table.
Reaction-Abs	9	REAL	Refer to Note 1 after this table.
Reaction-Mean Deviation	9	REAL	Refer to Note 1 after this table.
Reaction-Maximum Deviation	9	REAL	Refer to Note 1 after this table.
Raw Result	9	ALPHA	Refer to Note 1 after this table. Note: Result field reporting the results for TDM's may contain a < (less than) or > (greater than) sign.
Calculated Result	9	ALPHA	Refer to Note 1 after this table. Note: Result field reporting the results for TDM's may contain a < (less than) or > (greater than) sign.
Calibration Rate	9	REAL	Refer to Note 1 after this table.

(2 of 4)

Table 3.14 Expanded Results, continued

Field	Length	Format	Notes
Positive or Negative	1	FLAG	0 = Negative 1 = Positive 2 = Not applicable 3 = Equivocal (DxC 600i or DxC 800i only) 4 = Non-reactive (DxC 600i or DxC 800i only) 5 = Reactive (DxC 600i or DxC 800i only) 6 = Not confirmed (DxC 600i or DxC 800i only) 7 = Confirmed (DxC 600i or DxC 800i only) 8 = Gray zone flag (DxC 600i or DxC 800i only) 9 = Access 2 LOW flag (DxC 600i or DxC 800i only) Refer to Note 3 after this table.
Units	2	CODE	Refer to Table 4.6 Unit Codes.
Instrument Codes	9	ALPHA	Refer to Table 4.7 Instrument Codes with Remarks.
Result Error 1	2	CODE	Refer to Table 4.8 Result Error Codes.
Result Error 2	2	CODE	Refer to Table 4.8 Result Error Codes.
Result Error 3	2	CODE	Refer to Table 4.8 Result Error Codes.
Result Error 4	2	CODE	Refer to Table 4.8 Result Error Codes.
Result Error 5	2	CODE	Refer to Table 4.8 Result Error Codes.
Result Error 6	2	CODE	Refer to Table 4.8 Result Error Codes.
Result Error 7	2	CODE	Refer to Table 4.8 Result Error Codes.
Result Error 8	2	CODE	Refer to Table 4.8 Result Error Codes.

(3 of 4)

Table 3.14 Expanded Results, continued

Field	Length	Format	Notes
Result Error 9	2	CODE	Refer to Table 4.8 Result Error Codes.
Result Error 10	2	CODE	Refer to Table 4.8 Result Error Codes.
Result Error 11	2	CODE	Refer to Table 4.8 Result Error Codes.
Result Error 12	2	CODE	Refer to Table 4.8 Result Error Codes.
Result Error 13	2	CODE	Refer to Table 4.8 Result Error Codes.
Result Error 14	2	CODE	Refer to Table 4.8 Result Error Codes.
Result Error 15	2	CODE	Refer to Table 4.8 Result Error Codes.
Result Error 16	2	CODE	Refer to Table 4.8 Result Error Codes.
Intermediate ADC 1	5	NUMERIC	
Intermediate ADC 2	5	NUMERIC	
Intermediate ADC 3	5	NUMERIC	
Intermediate ADC 4	5	NUMERIC	
Intermediate ADC 5	5	NUMERIC	
Dilution Factor	6	REAL	0.0-1000.0
Future Use Space	23	ALPHA	

(4 of 4)

Note 1: The number of decimal places is adjusted to fit the number into this field.

Note 2: The most significant character of the range of the Rack Number field will be located in the most significant character position of the Cup Number field.

Example: Rack 345, Cup 1 would be 45, 31.

Note 3: For DxC chemistries, 0 and 1 are used only for drugs of abuse (DAT). On DxC 600i or DxC 800i, values from the range 0, 1, 3-9, are used by some Access 2 results, depending on the specific chemistry. Refer to the *Access Immunoassay System Manual* for additional information. All other chemistries use 2.

Stream 702 - Results

Stream 702 - Function 81: Expanded Results

Example:

```
[-0,702,81,100497,140128,---88,-50331669,---916097,S221-----,22,-1,
385875993,-1,05A-,-E,#,0,0,####,####,###,###,0.5623355,0.5538675,
#####,-0.127546,#####,-0.003311,0.5732892,
0.0002016,0.0003895,-0.012754,0.5577027,0.0001692,0.0002893,-0.012754,
46.466198,#####2,-0,-----,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,
NO,NO,NO,NO,#####,#####,#####,#####,#####,1.0000,#####]AE<CR><LF>
```

(- Indicates space holder)

Stream 703 - Instrument Status

Stream 703 - Function 1: Power Up

Sent by: UniCel DxC

Purpose: Sent upon DxC power up or reset. This message can be turned off in the host setup option.

Table 3.15 Power Up

Field	Length	Format	Notes
Device ID	2	0-99	Operator assigned
Stream	3	703	
Function	2	01	
Date	6	DATE	
Time	6	TIME	
Software Set	10	ALPHA	Console Master Program Version Number
CV ICC BOOT	10	ALPHA	
CV ICC APP	10	ALPHA	
CV MSMC1 BOOT 1	10	ALPHA	
CV MSMC1 BOOT 2	10	ALPHA	
CV MSMC1 APP	10	ALPHA	
CV MSMC2 BOOT 1	10	ALPHA	
CV MSMC2 BOOT 2	10	ALPHA	
CV MSMC2 APP	10	ALPHA	
CV SM BOOT	10	ALPHA	
CV SM APP 21	10	ALPHA	
CV SM APP 11	10	ALPHA	
CV SM APP 12	10	ALPHA	
CV SM APP 13	10	ALPHA	
CV SM APP 14	10	ALPHA	
CV SM APP 31	10	ALPHA	
CV SM APP 41	10	ALPHA	
CV SM APP 42	10	ALPHA	
CV SM APP 43	10	ALPHA	

(1 of 2)

Table 3.15 Power Up, continued

Field	Length	Format	Notes
CV SM APP 51	10	ALPHA	
CV SM APP 52	10	ALPHA	
CV SM APP 61	10	ALPHA	
CV SM APP 62	10	ALPHA	
CV SM APP 63	10	ALPHA	
CV SM APP 71	10	ALPHA	
CV SM APP 81	10	ALPHA	
CV SM APP 82	10	ALPHA	
CV SM APP 91	10	ALPHA	
CV SM APP A1	10	ALPHA	
CV SM APP A2	10	ALPHA	
CV SM APP B1	10	ALPHA	
CV SM APP B2	10	ALPHA	
CV SM APP B3	10	ALPHA	
CV SM APP C1	10	ALPHA	
CV SM APP C2	10	ALPHA	
CV SM APP C3	10	ALPHA	
SPARE	10	ALPHA	All "#"s
SPARE	10	ALPHA	All "#"s
SPARE	10	ALPHA	All "#"s
SPARE	10	ALPHA	All "#"s
SPARE	10	ALPHA	All "#"s
SPARE	10	ALPHA	All "#"s
SPARE	10	ALPHA	All "#"s
SPARE	10	ALPHA	All "#"s
SPARE	10	ALPHA	All "#"s
SPARE	10	ALPHA	All "#"s
SPARE	10	ALPHA	All "#"s
SPARE	10	ALPHA	All "#"s

(2 of 2)

Example:

```
[ -0,703,01,100497,164320, v01.00----,v00.00----,v00.00----,v00.00----,
v00.00----,v00.00----,v00.00----,v00.00----,v00.00----,v01.06----,
v00.13----,v04.54----,v04.54----,v04.54----,v04.54----,v01.30----,
v00.89----,v00.89----,v00.89----,v00.89----,v00.89----,v04.54----,
v04.54----,v04.54----,v00.89----,v00.1B----,#####,v00.30----,
v00.89----,v00.89----,v00.89----,v00.89----,v00.89----,v00.89----,
v00.89----,v00.89----,#####,#####,#####,#####,
#####,#####,#####,#####,#####,#####,#####,
#####,#####]5E<CR><LF>
```

(- Indicates space holder)

Stream 703 - Function 2: Bidirectional On**Sent by:** UniCel DxC**Purpose:** This message is sent to the host when bidirectional communication is enabled. This message can be disabled in the host setup option.

Table 3.16 Bidirectional On

Field	Length	Format	Notes
Device ID	2	0-99	Operator assigned
Stream	3	703	
Function	2	02	
Date	6	DATE	
Time	6	TIME	
Software Set	10	ALPHA	DxC Software Version

Example:

```
[ -0,703,02,100497,203332,A4.4.54---]F0<CR><LF>
```

(- Indicates space holder)

Stream 703 - Function 3: Request Instrument State**Sent by:** Host**Purpose:** Sent to obtain the state of the UniCel DxC.

Table 3.17 Request Instrument State

Field	Length	Format	Notes
Device ID	2	00-99	Operator assigned
Stream	3	703	
Function	2	03	

Example:

```
[ -0,703,03]93<CR><LF>
```

(- Indicates space holder)

Stream 703 - Function 4: Instrument State

Sent by: UniCel Dx-C

Purpose: Sent in response to a Request Instrument State, 703-03 message from the host.

Table 3.18 Instrument State

Field	Length	Format	Notes
Device ID	2	00-99	Operator assigned
Stream	3	703	
Function	2	04	
Date	6	DATE	
Time	6	TIME	
Configuration	5	ALPHA	DxC
MC State	2	CODE	Refer to Table 4.11 Instrument Status Codes.
CC State	2	CODE	Refer to Table 4.11 Instrument Status Codes.
Future Use	2	ALPHA	
Instrument Mode	2	CODE	0 = Rack/cup priority 1 = Bar code priority.

Example:

`[-0,703,04,110497,105150,DxC-,02,15,##,-1]87<CR><LF>`

(- Indicates space holder)

Stream 703 - Function 17: End of Run

Sent by: UniCel DxC

Purpose: Sent when the DxC completes all programmed tests. This message can be turned off in the host setup option.

Table 3.19 End of Run

Field	Length	Format	Notes
Device ID	2	00-99	Operator assigned
Stream	3	703	
Function	2	17	
Date Queued	6	DATE	
Time Queued	6	TIME	

Example:

[-0,703,17,100497,140245]1F<CR><LF>

(- Indicates space holder)

Stream 704 - Setup Status

Stream 704 - Function 7: Request Installed Chemistries

Sent by: Host

Purpose: Sent to obtain the list of installed chemistries from the DxC database.

Table 3.20 Request Installed Chemistries

Field	Length	Format	Notes
Device ID	2	00-99	Operator assigned
Stream	3	704	
Function	2	07	

Example:

```
[00,704,07]8E<CR><LF>
```

Stream 704 - Function 8: Installed Chemistries

Sent by: UniCel DxC

Purpose: This message is sent in response to a Request Installed Chemistries, 704-07 message from the host.

Table 3.21 Installed Chemistries

Field	Length	Format	Notes
Device ID	2	00-99	Operator assigned
Stream	3	704	
Function	2	08	
Date	6	DATE	
Time	6	TIME	
Com Chem ID 1	4	CHEM CODE	Refer to Table 4.4 and Table 4.5 for Chemistry (Chem) Codes.
Com Chem ID N	4	CHEM CODE	Refer to Table 4.4 and Table 4.5 for Chemistry (Chem) Codes.

Note: Only chemistry codes of chemistries that are configured on the DxC will be included in the message. The last chemistry code transmitted will be "0000" to indicate that no more chemistry codes follow. If no chemistries are configured on the instrument, only "0000" chemistry code will be listed.

Example:

```
[ -0, 704, 08, 110497, 105518, 01A-, 01B-, 04A-, 02A-, 09D-, 08D-, 05D-, 03E-, 06D-,
43D-, 07D-, 14B-, 56B-, 35A-, 31A-, 31B-, 54A-, 84A-, 10A-, 24C-, 24B-, 93B-, 30A-,
30B-, 85A-, 86A-, 05A-, 98A-, 59A-, 44A-, 32A-, 32B-, 40A-, 87A-, 89A-, 12A-, 62A-,
46B-, 72C-, 36A-, 06A-, 14A-, 83A-, 79B-, 52B-, 51B-, 53B-, 55A-, 33A-, 50A-, 50B-,
07B-, 64A-, 63A-, 48A-, 92A-, 66A-, 67C-, 68C-, 65A-, 93C-, 94A-, 91A-, 11A-, 42B-,
42C-, 88A-, 88B-, 41A-, 69C-, 70C-, 71B-, 0000] BE<CR><LF>
```

(- Indicates space holder)

Stream 704 - Setup Status

Stream 704 - Function 8: Installed Chemistries

CHAPTER 4 Tables/Codes

Table of Contents

Tables/Codes	4-2
Test Type Codes	4-2
Test Type Codes	4-2
Sample Type Codes	4-3
Sample Type Codes	4-3
Age Unit Codes	4-3
Chemistry Codes	4-4
Chemistry Codes	4-4
Access Chemistry Codes	4-8
Unit Codes	4-10
Unit Codes	4-10
Instrument Codes	4-12
Instrument Codes	4-12
Result Error Codes	4-18
Result Error Codes	4-18
Result Error Code Descriptions	4-20
Special Calculation Status Codes	4-24
Special Calculation Status Codes	4-24
Instrument Status Codes	4-25
Instrument Status Codes	4-25
Access Codes Converted to DxC Codes for DxC 600i Reporting	4-26
Access Codes Converted to DxC Codes (DxC 600i only)	4-26
Key Code Conversion for Local Languages	4-27
Key Code Conversion	4-27

Tables/Codes

Test Type Codes

Test Type Codes

Table 4.1 Test Type Codes

Code	Description
CA ^a	Calibrator
CO	Control
EX ^a	Extinction
RO	Routine
SC	STAT control
ST	STAT

^a Not programmable by the Host.

Sample Type Codes

Sample Type Codes

Table 4.2 Sample Type Codes

Code	Description
SE	Serum
SF	Spinal fluid
UR	Urine (Random)
TU	Timed urine (Can not be used with control or STAT Control)
PL	Plasma
OT	Other (Refer to Note after this table)
AM	Amniotic (UniCel DxC 600i or DxC 800i only)
BL	Blood
CE	Cervical (UniCel DxC 600i or DxC 800i only)
SA	Saliva (UniCel DxC 600i or DxC 800i only)
SY	Synovial (UniCel DxC 600i or DxC 800i only)
UT	Urethral (UniCel DxC 600i or DxC 800i only)

Age Unit Codes

Table 4.3 Age Unit Codes

Code	Description
1	Hours
2	Days
3	Weeks
4	Months
5	Years

Chemistry Codes

Chemistry Codes

Table 4.4 Chemistry (Chem) Codes

Code	Chemistry
01A	Sodium (NA)
01B	Potassium (K)
02A	Carbon Dioxide (CO2)
03A	Creatinine (CREA)
03E	Creatinine (CREm)
04A	Chloride (CL)
05A	Urea Nitrogen (BUN)
05D	Urea Nitrogen (BUNm)
06A	Glucose (GLU)
06D	Glucose (GLUCm)
07A	Total Protein (TP)
07B	Micro Protein (M-TP)
07D	Total Protein (TPm)
08A	Albumin (ALB)
08D	Albumin (ALBm)
08E	Prealbumin (PAB)
08M	Microalbumin (MA)
09D	ISE Calcium (CALC)
10A	Amylase (AMY)
10C	P-Amylase (PAM)
11A	Total Bilirubin (TBIL)
12A	Direct Bilirubin (DBIL)
13B	Infinity Lithium ^a (LI)
14A	Hemoglobin, Total (Hb)
14B	Hemoglobin, A1c (A1c)
14I	Hemoglobin Total (Hb2)
14J	Hemoglobin, A1c (A1c2)

(1 of 4)

Table 4.4 Chemistry (Chem) Codes, continued

Code	Chemistry
15A	Homocysteine (HCY)
24B	Apolipoprotein B (ApoB)
24C	Apolipoprotein A (ApoA)
30A	Aspartate Aminotransferase (AST)
30B	Aspartate Aminotransferase (Pyridoxal-5'-phosphate) (AST-)
31A	Alanine Aminotransferase (ALT)
31B	Alanine Aminotransferase (Pyridoxal-5'-phosphate) (ALT-)
32A	Creatine Kinase (CK)
32B	Creatine Kinase (N-Acetyl-L-cysteine) (CK-)
33A	Lactate Dehydrogenase (LD)
34A	Lactate Dehydrogenase (LD-P)
34B	LDL Cholesterol (LDLD)
35A	Alkaline Phosphatase (ALP)
36A	Gamma Glutamyl Transferase (GGT)
40A	Creatine Kinase MB (CKMB)
41A	Uric Acid (URIC)
42B	Triglycerides GPO (TG)
42C	Triglycerides GPO Blanked (TG-B)
43D	Phosphorus (PHOSm)
43E	Phosphorus (PHS)
44A	Cholesterol (CHOL)
46B	Iron (FE)
48A	Magnesium (MG)
50A	Lipase (LIPA)
50B	Lipase Wash (LIWA)
50M	Haptoglobin (HPT)
51B	Immunoglobulin G (Ig-G)
52B	Immunoglobulin A (Ig-A)
53B	Immunoglobulin M (Ig-M)
54A	Ammonia (AMM)

(2 of 4)

Table 4.4 Chemistry (Chem) Codes, continued

Code	Chemistry
55A	Lactate (LAC)
56A	Ethanol Alcohol (ETOH)
59A	Cholinesterase (CHE)
60A	Lipase (LIP)
61M	Complement (C3)
61R	Complement (C4)
62B	Digoxin (DIGN)
63A	Methaqualone (METQ)
64A	Methadone (METD)
65A	Propoxyphene (PROX)
66A	Phencyclidine (PCP)
67C	Phenobarbital (PHE)
68C	Phenytoin (PHY)
69C	Theophylline (THE)
70C	Tobramycin (TOB)
71B	Transferrin (TRFN)
72C	Gentamicin (GEN)
72M	Vancomycin (VANC)
73A	Urea Nitrogen (UREA)
73D	Urea Nitrogen (UREAm)
79B	Iron Binding Capacity, Total (IBCT)
83D	HDL Cholesterol (HDL D)
84A	Amphetamines (AMPH)
85A	Barbiturates (BARB)
86A	Benzodiazepine (BENZ)
86B	Benzodiazepine (BNZG)
87A	Cocaine Metabolite (COCM)
88A	Cannabinoid 100 ng (THC)
88B	Cannabinoid 20 ng (THC2)
88C	Cannabinoid 50 ng (THC5)

(3 of 4)

Table 4.4 Chemistry (Chem) Codes, continued

Code	Chemistry
89A	C-Reactive Protein (CRP)
89E	High Sensitivity C-Reactive Protein (CRPH)
89G	C-Reactive Protein (C-RP)
90A	T-Uptake (TU)
91A	Thyroxine (T4)
92A	Opiate 300 ng (OP)
92B	Opiate 2000 ng (OP2)
93B	Antistreptolysin (ASO-)
93C	Rheumatoid Factor (RF)
94B	Salicylate (SALY)
94M	Acetaminophen (ACTM)
95A	Valproic Acid (VPA)
98A	Carbamazepene (CAR)
99B	Icterus (ICTER)
99C	Lipemia (LIPEM)
99D	Hemolysis (HEMOL)
99E ^b	Cartridge Chemistry Wash Solution (CCWA)

(4 of 4)

^a All trademarks are the property of their respective owners.^b Do not use as a result.

Access Chemistry Codes

The following chemistry (assay) codes pertain to chemistries run on the Access 2 analyzer. The Access 2 analyzer is part of the UniCel DxC 600i or DxC 800i System.

Table 4.5 Access Chemistry Codes

Code	Chemistry
A02	Alpha-fetoprotein (AFP)
A03	OV Monitor (OV125)
A04	BR Monitor (BR153)
A05	GI Monitor (GI199)
A06	CEA (CEA)
A07	Chlamydia (ChlAg)
A08	Chlamydia Block (ChlBk)
A10	Cortisol (Cort)
A11	Digoxin (Dig-)
A12	Ostase ^a (Ostas)
A13	Toxo IgM II (Tox-M)
A14	Estradiol (ESTRD)
A15	Unconjugated Estriol (uE3)
A16	Testosterone (Testo)
A17	Total BhCG (TBhCG)
A18	Diluted Total BhCG (DlhCG)
A19	Ferritin (Ferr)
A22	Ultrasensitive hGH (hGH)
A23	Thyroglobulin (TgAg)
A30	Myoglobin (MYO)
A36	hFSH (hFSH)
A39	hLH (hLH)
A40	Insulin (Insul)
A42	Prolactin (PRL)
A43	Progesterone (Prog)
A46	Rubella IgG (Rub-G)
A48	Theophylline (Theo)
A49	Creatine Kinase-MB (CK-MB)

(1 of 2)

Table 4.5 Access Chemistry Codes, continued

Code	Chemistry
A50	Total T3 (TotT3)
A51	Total IgE (IgE)
A52	Total T4 (TotT4)
A53	Toxo IgG (Tox-G)
A57	HYPERsensitive hTSH (TSH)
A59	Uptake (TU-)
A60	Vitamin B ₁₂ (B12)
A61	Free T4 (FT4)
A62	Hybritech ^a PSA (PSA)
A63	Hybritech ^a free PSA (fPSA)
A67	AccuTnI ^a (cTnI)
A68	Thyroglobulin Antibody (TgAb)
A72	Intrinsic Factor Antibody (IFAb)
A73	Dehydroepiandrosterone Sulfate (DHEAS)
A74	Triage BNP ^a (BNP)
A75	Folate (FOL2)
A76	RBC Folate (RBC2)
A77	Diluted Alpha-fetoprotein (dAFP)
A78	Fast hTSH (FTSH)
A80	Free T3 (FT3)

(2 of 2)

^a All trademarks are the property of their respective owners.

Unit Codes

Unit Codes

Table 4.6 Unit Codes

Code	Unit
00	mg/dL
01	mg/L
02	g/dL
03	g/L
04	mmol/L
05	μmol/L
06	mEq/L
07	nKat/L
08	μKat/L
09	IU/L
10	μg/mL
11	ng/mL
12	μg/dL
13	μg/L
14	nmol/L
15	Ku.u
16	U/L
17	Other
18	%
19	mA
20	mA/min
21	IU/mL
22	U/mL
23	Rate
24	ng/dL
25	μIU/mL
26	mIU/mL

(1 of 2)

Table 4.6 Unit Codes, continued

Code	Unit
27	KU/L
28	nIU/dL
29	mIU/L
30	Positive
31	Negative
32	pg/mL
33	pg/dL
34	mg/mL
35	ng/L
36	pmol/L
37	%Uptake
38	%GHb
39	%A1c
40	GPL
41	MPL
42	APL
43	RLU
44	S/CO (sample to cutoff ratio)
45	AU/mL (allergy units or arbitrary units)
46	AU/L (arbitrary units)
47	mAU/L (arbitrary units)
48	mAU/mL (arbitrary units)
49	(1) SI (T-uptake SI units)
50	%Supp

(2 of 2)

Instrument Codes

Instrument Codes

Table 4.7 Instrument Codes with Remarks

Instrument Report Printout	Printed in Remark Area of Instrument Report	Inst. Code	Description
"a printed result"			Result was calculated
"a printed result"	CHECK RESULT		DxC has processed the sample prior to Access 2. Results for specific Access 2 tests should be checked for potential sample-to-sample carryover. This error code applies to Access 2 assays where a small quantity of carryover could significantly increase the results. (DxC 600i or DxC 800i only)
"a printed result"	CRITICAL HIGH		Exceeds Critical range high
"a printed result"	CRITICAL LOW		Exceeds Critical range low
"a printed result"	HIGH		Exceeds reference range high
"a printed result"	LOW		Exceeds reference range low
"a printed result"	OIR LO		Result is outside of reportable range – low (DxC 600i or DxC 800i only)
"a printed result"	TEMP ERR	T	RA temperature
"a printed result"	TEMP ERR	T	PP temperature
"a printed result"	TEMP ERR	T	ISE Temperature
"a printed result"		C	Calibration overridden
"a printed result"		D	Reagent Days exceeded
"a printed result"		E	Calibration time extended
"a printed result"		J	Slope offset adjustment
"a printed result"		M	Set point modification
"a printed result"		O	ORDAC sample
"a printed result"		P	Default Panel

(1 of 6)

Table 4.7 Instrument Codes with Remarks, continued

Instrument Report Printout	Printed in Remark Area of Instrument Report	Inst. Code	Description
"a printed result"		R	Reagent Expired
"a printed result"		V	Rerun result
"a printed result"		S	Substrate expired
< X (X = lower inst range)			Less than X - OIR LO response (no suppressed result) from the CAU
< X (X = lower reportable range)			Less than X - ORR LO response from console but result is not suppressed
> Y (Y = upper reportable range)			Greater than Y - ORR HI response from console but result is not suppressed
>Y (Y = upper inst range)			Greater than Y - OIR HI response from CAU but result is not suppressed
BACK TO BACK	Calibration failed		Back to back
CAL FAILED	Calibration failed		Calibration failed
CAL TIMED OUT			Calibration timed out
CHEM ABORTED			Too many retries
CHEM BYPASSED			Channel bypassed
CHEM DB ERROR			Chemistry DB Error
CHEM NEEDS CAL			Chemistry needs calibration
CHEM NOT RUN			Not run
DAC ERROR			ISE DAC failure
DIL1 LEV SENSE			Level sense error in diluent reagent
DISABLED			Module Disabled
DIVISION ERROR			0 denominator
HI REAGENT			High reagent level
INCOMPLETE			Not all programmed tests have a valid result

(2 of 6)

Table 4.7 Instrument Codes with Remarks, continued

Instrument Report Printout	Printed in Remark Area of Instrument Report	Inst. Code	Description
INCOMPLETE	CHEM NEEDS CAL		No curve – calibration required (DxC 600i or DxC 800i only)
	CHEM NOT RUN		Cannot run test because: <ul style="list-style-type: none"> • Sample is not detected by CTA • Insufficient supply at CTA/ Access 2 • Loss of communication between CTA and Access 2 • Access 2 is not ready • CTA - No chemistry definition at Access 2 • CTA/Access 2 waste bag is full • AV expired • Access 2 cover is open (DxC 600i or DxC 800i only)
	DEVICE ERROR		Mechanical device error detected (DxC 600i or DxC 800i only)
	NO SAMPLE DETECT		Insufficient sample quantity at Access, or no sample detected - CTA level sense (DxC 600i or DxC 800i only)
	NO SAMPLE ID		Cannot run Access test because sample ID is not programmed (DxC 600i or DxC 800i only)
INSUFF DIL1			Insufficient diluent reagent
INSUFF REAG			Insufficient reagent
INVALID CHEM			Invalid chemistry
LEV SENSE ERR			Low reagent fluid
LO REAGENT			Low reagent level

(3 of 6)

Table 4.7 Instrument Codes with Remarks, continued

Instrument Report Printout	Printed in Remark Area of Instrument Report	Inst. Code	Description
MATH ERR	Calibration required		Math error
MOTION ERR			Motion error
NO DIL1 ON SYS			No diluent reagent on the reagent carousel
NO REAG ON SYS			Reagent not on board
NO SAMP DETECT			Low sample fluid
NO SAMP DETECT			No sample detected
NO SAMP ON SYS			Sample not on board
NOT ACCEPTED			Wrong state for sample program message
PENDING TEST			Not ready
PROBE OBSTRUCT			Probe obstruction detected
RACK CONFLICT			Test cancelled because sample is in wrong rack (HDL, IBCT, HbA1c chemistries)
RANGE	Calibration failed		Range error
RECOVERY	Calibration failed		Recovery error
RESULT ERR			Result error
Results Suppressed	AG EXCESS		Antigen Excess
Results Suppressed	BL ABS HI		Blank absorbance high
Results Suppressed	BL ABS LO		Blank absorbance low
Results Suppressed	BL MAX DEV		Blank outlier (maximum deviation)
Results Suppressed	BL MEAN DEV		Blank noise (mean deviation)
Results Suppressed	BL RATE HI		Blank rate high
Results Suppressed	BL RATE LO		Blank rate low
Results Suppressed	CAL REF DRIFT		Excessive reference drift - sample to calibration
Results Suppressed	INDETERMINATE		Unable to determine results - near cutoff value (DxC 600i or DxC 800i only)

(4 of 6)

Table 4.7 Instrument Codes with Remarks, continued

Instrument Report Printout	Printed in Remark Area of Instrument Report	Inst. Code	Description
Results Suppressed	INIT ABS HI		Initial Absorbance high
Results Suppressed	INIT ABS LO		Initial Absorbance too low
Results Suppressed	INIT ADC HI		Initial ADC error high
Results Suppressed	INIT ADC LO		Initial ADC error low
Results Suppressed	INIT COND HI		Initial conductance high
Results Suppressed	INIT COND LO		Initial conductance low
Results Suppressed	INIT RATE HI		Initial rate high
Results Suppressed	INIT RATE LO		Initial rate low
Results Suppressed	K REQ		Bad K value
Results Suppressed	NA/K REQ		Bad Na value
Results Suppressed	OIR HI		Out of Instrument Range high
Results Suppressed	OIR LO		Out of Instrument Range low
Results Suppressed	ORDAC HI		Out of Instrument range ORDAC high
Results Suppressed	ORDAC LO		Out of Instrument range ORDAC low
Results Suppressed	URDAC HI		URDAC high
Results Suppressed	URDAC LO		URDAC low
Results Suppressed	ORR HI		Out of reportable range high
Results Suppressed	ORR LO		Out of reportable range low
Results Suppressed	ORR O HI		Out of ORDAC reportable range high
Results Suppressed	OVERFLOW		Number overflow error
Results Suppressed	REF NOISE		Reference signal noise
Results Suppressed	RLU RANGE		Relative light units (RLUs) exceeds measuring range of Access luminometer (DxC 600i or DxC 800i only)
Results Suppressed	RX ABS HI		Reaction absorbance high
Results Suppressed	RX ABS LO		Reaction absorbance low

(5 of 6)

Table 4.7 Instrument Codes with Remarks, continued

Instrument Report Printout	Printed in Remark Area of Instrument Report	Inst. Code	Description
Results Suppressed	RX ERR		Reaction error
Results Suppressed	RX MAX DEV		Reaction outlier (maximum deviation)
Results Suppressed	RX MEAN DEV		Reaction noise (mean deviation)
Results Suppressed	RX NOISE	A	Erratic ADC (cup)
Results Suppressed	RX RATE HI		Reaction rate high
Results Suppressed	RX RATE LO		Reaction rate low
Results Suppressed	SAMP REF DRIFT		Excessive reference drift - sample to sample
Results Suppressed	SAMPLE NOISE		Sample signal noise
Results Suppressed	SUB DEPL		Substrate depleted
Results Suppressed	INCUBATOR TEMP		Access temperature error (DxC 600i or DxC 800i only)
Results Suppressed	SUBSTRATE TEMP		Access temperature error (DxC 600i or DxC 800i only)
Results Suppressed	CAROUSEL TEMP		Access temperature error (DxC 600i or DxC 800i only)
SENSITIVITY	Calibration failed		Sensitivity error
SEV RECOVERY	Calibration required		Severe recovery error
SEV SENSITIVITY	Calibration required		Severe Sensitivity error
SPAN	Calibration failed		Span error
TEST STOPPED			Test was aborted by console - recovery, cleanup, etc.
TOO MANY TESTS			Too many tests on CC

(6 of 6)

Result Error Codes

Result Error Codes

Table 4.8 Result Error Codes

Host Code	Description
AB	Not all programmed tests have a valid result
AX	Antigen excess
BH	Blank absorbance high
BL	Blank absorbance low
BN	Blank noise (mean deviation)
BO	Blank outlier (maximum deviation)
CH	Initial conductance high
CL	Initial conductance low
CR	Check result Access result may have been contaminated by Dx C processing of sample
DH	Out of instrument range high
DL	Out of instrument range low
DR	Reference signal noise
DS	Sample signal noise
EA	Erratic ADC
EC	Excessive reference drift - sample to calibration
ES	Excessive reference drift - sample to sample
GH	URDAC high
GL	URDAC low
GT	Greater than upper instrument or reportable range
HI	Initial ADC error high
HR	Reaction absorbance high
IA	Initial absorbance either too high or too low
IK	Bad K value
IL	Initial rate too low
IN	Bad NA value
IR	Initial rate too high

(1 of 2)

Table 4.8 Result Error Codes, continued

LI	Initial ADC error low
LO	Result is lower than minimum reportable result value defined in APF (Assess LOW flag)
LR	Reaction absorbance low
LT	Less than lower instrument or reportable range
OF	Number overflow error
OH	Out of instrument range ORDAC high
OK	Result was calculated
OL	Out of instrument range ORDAC low
OR	Relative light units (RLUs) are outside acceptable luminometer measuring range
RE	Reaction error
RH	Reaction rate high
RL	Reaction rate low
RN	Reaction noise (mean deviation)
RO	Reaction outlier (maximum deviation)
SD	Substrate depleted
SH	Blank rate high
SL	Blank rate low
TM	Temperature error
UH	Out of reportable range high
UL	Out of reportable range low
UO	Out of ORDAC reportable range high
ZD	Zero denominator

(2 of 2)

Result Error Code Descriptions

The following table provides the Dx C result error code descriptions.

Table 4.9 Error Code Descriptions

Error	Description
Antigen Excess	Failed antigen excess test for immunoglobulin reagents.
Blank Absorbance (Hi/Lo)	The mean absorbance measured during the reagent blank spin cycle (read window). The units are in absorbance and will characterize reagent quality.
Blank Outlier Maximum Deviation	A single blank absorbance data point obtained during the reagent blank read window deviates more than allowed from the line of regression. The Blank Deviation units are measured in delta absorbance. Blank Max Deviation will check for a constant rate during the blank spin cycles.
Blank Mean Deviation	The average difference between the absorbance readings and the line of regression is greater than allowed. This will check for a constant rate during the blank spin cycles.
Calibrator Order Error	Occurs in rack mode when not all calibrator levels are in ascending order (low to high) and/or not all levels reside in the same rack.
Erratic ADC's	The difference between the high and low value of the four reference electrode readings (taken in milliseconds apart) has exceeded the limits allowed. Units are in ADC's and are a measurement of noise in the ISE system.
Excessive Reference Drift	An ISE reference electrode measurement for a sample drifted above the reference electrode measurement from the calibration and exceeded the limit. The units are in ADC's and are a measurement of ISE reference electrode drift.
Initial Absorbance (Hi/Lo)	The reaction absorbance data taken from the first spin cycle after sample inject has exceeded specifications. This is a measure of sample integrity.
Initial Rate (Hi/Lo)	The reaction absorbance data obtained between 2 and 17 seconds after sample inject exceeds specifications.
Insufficient Reagent Strength	GLUCm and BUNm reagent conductivity is below specifications.
Iteration Error	Too many iterations to fit nonlinear serum to multilevel calibration.
Noise	During a particular spin cycle, the average difference between an absorbance reading and the line of regression exceeded specifications.

(1 of 4)

Table 4.9 Error Code Descriptions, continued

Error	Description
Outlier	During a particular spin cycle, the deviation of a single absorbance reading with respect to the line of regression exceeded specifications.
Out of Instrument Range (Hi/Lo)	The recovered value exceeds the value that the instrument reports.
Out of ORDAC Range (Hi/Lo)	Reaction exceeded the range of ORDAC value that the instrument reports.
Reaction Absorbance (Hi/Lo)	The mean absorbance measured during the reaction spin cycles.
Reaction Maximum Deviation	A single reaction absorbance data point obtained during the reaction read window deviates more than allowed from the line of regression. The units are in delta absorbance and will check for a constant rate during the reaction read window.
Reaction Mean Deviation	The average difference between the reaction absorbance readings and the line of regression is greater than allowed. This will check for a constant rate during the reaction spin cycle.
Reaction Rate (Hi/Lo)	The rate calculated during the reaction read window. Units are in delta absorbance.
Substrate Depletion	The difference (delta) between the initial absorbance taken after sample inject and the final absorbance data point within the reaction read window exceeds specifications.
Blank Rate (Hi/Lo)	Rate during blank timing period is not within specifications.
Temperature	The operating temperature of the system is beyond 0.1°C from the set point value. All chemistries will be reported with a flag.
Ig-A ORDAC (Hi/Lo)	Result is outside the specifications for the special below normal range Ig-A ORDAC range.
Back to Back (Linear Cal)	Difference between usable calibrator replicates. Units are dependent on the reagent. Measure of system precision during calibration.

(2 of 4)

Table 4.9 Error Code Descriptions, continued

Error	Description
Back to Back (ISE Cal)	Measure of system precision during calibration. The difference between the two middle Sample-Reference replicate values within a calibration level. If the back-to-back limit is exceeded, the calibration fails. Only the two consecutive replicates that meet back-to-back are used; all other replicates are discarded. Select Calibrator Acceptance Limits under options in the Rgts/Cal screen to view current calibration limits.
DAC	Verifies proper operation of electrodes. The Digital-to-Analog Conversion (DAC) check mimics expected ranges of electrode output. Failure to detect an output within a predetermined range for each channel will trigger a DAC INITIALIZATION ERROR HAS OCCURRED... flag.
Erratic ADC's	The difference between the high and low value of the four electrode readings (taken milliseconds apart) must be within the delta limit. Measure of noise for modular chemistries.
Math Error 1	Multi-point calibration error. Algorithm produces curve parameters that are very large or small and/or rounding errors will give an invalid standard curve. Curve parameters are not reported and calibration may not be overridden.
Math Error 2	Computer is not able to calculate the parameters needed to determine the other math errors. Curve parameters are not reported and calibration may not be overridden.
Math Error 3	The range of data points is on a very small percent of the overall calculated curve. Curve is not sensitive enough to get valid results. Curve parameters are not reported and calibration may not be overridden.
Math Error 4	Future use.
Math Error 5	Calculated curve has multiple inflection points within the dynamic range indicating more than one result is possible for the same data point. Curve parameters are not reportable and calibration may not be overridden.
Math Error 6	Data does not properly define the curve. Only one point is near the inflection point of the curve. Curve parameters are not reported and calibration may not be overridden.
Math Error 7	Future use.

(3 of 4)

Table 4.9 Error Code Descriptions, continued

Error	Description
Math Error 8	Data does not properly define the curve. Only one point is near the inflection point of the curve. Results need to be evaluated by the user. Curve parameters are reported and calibration MAY be overridden. A less severe MATH ERROR 6.
Math Error 9-15	Future use.
Math Error 16	Occurs when calibrations include division by 0.
Range (ISE cal)	Measure of reagent and electrode performance. Determines if the calibrator set points fall within the acceptable range of ADC values.
Range (linear cal)	Absorbance or rate range that the usable calibrator replicates must fall between for a valid calibration.
Recovery (Multi-point)	Determines if recovery errors of each calibration level are too large for a calibration to be useful. The difference between the recovered concentration and the actual concentration exceeds specifications. It is overrideable.
Recovery Severe (Multi-point)	Non-overrideable extreme recovery error. Refer to Recovery (Multi-point) above.
Sensitivity (Multi-point)	A check of the calibration slope for each calibrator level and several intermediate levels to check if the response is too high or low. Is overrideable.
Sensitivity Severe (Multi-point)	Non-overrideable extreme sensitivity error. Refer to Sensitivity (Multi-point) above.
Span (ISE cal)	Difference between the average ADC value of calibrator 1 and the average ADC value of calibrator 2. Measure of sensitivity.
Span (Multi-point)	Difference between consecutive calibrator rates or absorbances. There can be up to 6 multi-point spans for level 0-1, 1-2, 2-3, 3-4, 4-5, 5-0. This is a measure of sensitivity over the entire standard curve.
ADC Electrolyte Error	More than 1 second elapsed between the start and end of the ADC conversion process.

(4 of 4)

Special Calculation Status Codes

Special Calculation Status Codes

Table 4.10 Special Calculation Status Codes

Code	Description
AB	One of the involved chemistries was not run
OK	Valid result
UN	Units for the involved chemistries are not compatible
ZD	Denominator of a ratio is zero

Instrument Status Codes

Instrument Status Codes

Table 4.11 Instrument Status Codes

Code	Cartridge Chemistry Instrument Status	Code	Modular Chemistry Instrument Status
01	Stopped	01	Stopped
02	Pause initiated	02	Standby
		03	System home
		04	Reagent load
05	Running	05	Prime
06	Idle shutdown in progress (future use)		
07	Idle (future use)	07	Running
08	Initializing (future use)	08	Maintenance
09	Reagent load		
10	Homing		
11	Priming		
12	PVT	12	Boot (future use)
		13	Pause initiated
		14	Waiting
15	Standby	15	System idle (future use)
18	Procedure in progress		
21	Waiting		

Access Codes Converted to DxC Codes for DxC 600i Reporting

Access Codes Converted to DxC Codes (DxC 600i only)

Table 4.12 Access Codes Converted to DxC Codes (DxC 600i only)

Access 2			DxC			
Description	Value	Flags	Remarks	Inst Code	Result	Host Code
Normal result	result	none	none	none	print result	none
Indeterminate	"No value"	IND	INDETERMINATE	none	suppress	RE
Out of range high	>highest	OVR	none	none	print: >X	GT
Curve expired	result	CEX	none	E	print result	none
Substrate expired	result	EXS	none	S	print result	none
Lot expired	result	LEX	none	R	print result	none
Lower than reportable result	result	LOW	OIR LO	none	print result	LO
	< lowest	none	none	none	print: < X	LT
No curve – cal required	Raw RLUs	NCR	CHEM NEEDS CAL	NA ^a	incomplete	NA
Pack expired	result	PEX	none	D	print result	none
Control outside QC limits	NA	QCF	NA	NA	NA	NA
Control lot expired	NA	QEX	NA	NA	NA	NA
Quantity not sufficient	"No value"	QNS	NO SAMP DETECT	NA	incomplete	NA
Reflex test	NA	RFX	NA	NA	NA	NA
RLU is out of measuring range	"Cancelled"	RLU	RLU RANGE	none	suppress	OR
Instrument error	NA	SYS	DEVICE ERROR	NA	incomplete	NA
Incubator temperature error	"No value"	TRI	INCUBATOR TEMP	none	suppress	TM
Substrate temperature error	"No value"	TRS	SUBSTRATE TEMP	none	suppress	TM
Carousel temperature error	"No value"	TRW	CAROUSEL TEMP	none	suppress	TM
Can't run test (Undetected by CTA)	"Cancelled"	other	CHEM NOT RUN	NA	incomplete	NA

^a NA = Not applicable.

Key Code Conversion for Local Languages

Key Code Conversion

For communications between the host computer and the UniCel Dx_C, local language character handling is accomplished using the extended ASCII character set as defined by IBM for the IBM PC and IBM compatible units. The correct interpretation of foreign language characters on the UniCel Dx_C will require the use of 8-bit communication between the system and the host computer.

The defined foreign language characters and their extended ASCII character codes are shown below. No other extended ASCII characters are supported by the UniCel Dx_C system for display, printing and host communication.

DECIMAL VALUE	HEXA DECIMAL VALUE	0	16	32	48	64	80	96	112
0	0	BLACK (NULL)	BLANK (SPACE)	0	@	P	'	p	
1	1	☺	☹	!	1	A	Q	a	q
2	2	☺	☹	"	2	B	R	b	r
3	3	♥	!!	#	3	C	S	c	s
4	4	♦	♠	\$	4	D	T	d	t
5	5	♣	§	%	5	E	U	e	u
6	6	♠	■	&	6	F	V	f	v
7	7	•	↑	'	7	G	W	g	w
8	8	•	↑	(8	H	X	h	x
9	9	○	↓)	9	I	Y	i	y
10	A	○	→	*	:	J	Z	j	z
11	B	♂	←	+	;	K	l	k	{
12	C	♀	└	,	<	L	\	l	:
13	D	♫	↔	-	=	M	l	m	}
14	E	♫	▲	•	>	N	^	n	~
15	F	☼	▼	/	?	O	—	o	△

DECIMAL VALUE	HEXA DECIMAL VALUE	128	144	160	176	192	208	224	240
0	0	Ç	É	á	⋮	⋮	⋮	α	≡
1	1	ü	æ	í	⋮	⋮	⋮	β	±
2	2	é	Æ	ó	⋮	⋮	⋮	Γ	≥
3	3	â	ô	ú	⋮	⋮	⋮	π	≤
4	4	ä	ö	ñ	⋮	⋮	⋮	Σ	∫
5	5	à	ò	Ñ	⋮	⋮	⋮	σ	∫
6	6	å	û	ä	⋮	⋮	⋮	μ	÷
7	7	ç	ù	ó	⋮	⋮	⋮	τ	≈
8	8	ê	ÿ	¿	⋮	⋮	⋮	ϕ	°
9	9	ë	Ö	⋮	⋮	⋮	⋮	Θ	•
10	A	è	Ü	⋮	⋮	⋮	⋮	Ω	•
11	B	ï	ç	½	⋮	⋮	⋮	δ	√
12	C	î	£	¼	⋮	⋮	⋮	∞	n
13	D	ì	¥	ì	⋮	⋮	⋮	φ	²
14	E	Ä	ŕ	«	⋮	⋮	⋮	€	■
15	F	Å	ƒ	»	⋮	⋮	⋮	∩	BLACK FF

A009699L.EPS

Figure 4.1 Key Code Conversion Chart

Index

A

Access chemistry codes 4-8
Access codes converted to DxC codes for DxC 600i
reporting 4-26

B

Bidirectional protocol 1-11
Control characters 1-11
Data transfer 1-12
Downloading sample programming 1-18
Line bidding 1-11
Solicited messages 1-17
Unsolicited messages 1-15

C

Cable length 1-6
Chemistry codes 4-4
Connector pin assignments 1-6

D

Definitions 1-4

I

Instrument codes 4-12
Instrument status codes 4-25
Interface protocol options 1-5
Invalid/valid characters 2-2, 3-2

K

Key code conversion for local languages 4-27

M

Message format 2-2, 3-2
Message streams 2-3, 3-3
Stream 700-Special functions
Function 1: Are you there? 3-4
Function 2: Host setup 3-4
Stream 701-Sample/cup program
Function 1: Sample/cup program 3-6
Function 2: Sample/cup return status 3-9
Function 3: Clear rack/sample IDs 3-11
Function 4: Clear rack/sample IDs status 3-13
Function 6: Host query rack/sample IDs 3-14
Stream 702-Results
Function 1: Cup header 3-15
Function 3: Test results 3-17

Function 5: End of cup 3-21
Function 11: Special calculations results 3-22
Function 13: Timed urine results 3-23
Function 81: Expanded results 3-24
Stream 703-Instrument status
Function 1: Power up 3-29
Function 2: Bidirectional on 3-31
Function 3: Request instrument state 3-31
Function 4: Instrument state 3-32
Function 17: End of run 3-33
Stream 704-Setup status
Function 7: Request installed chemistries 3-34
Function 8: Installed chemistries 3-34
Stream 800-Special functions
Function 1: Are you there? 2-4
Function 2: Host setup 2-4
Stream 801-Sample/cup program
Function 1: Sample/cup program 2-6
Function 2: Sample/cup return status 2-10
Function 3: Clear rack/sample IDs 2-12
Function 4: Clear rack/sample IDs status 2-13
Function 6: Host query rack/sample IDs 2-14
Stream 802-Results
Function 1: Cup header 2-15
Function 3: Test results 2-17
Function 5: End of cup 2-21
Function 11: Special calculations results 2-21
Function 13: Timed urine results 2-23
Function 81: Expanded results 2-24
Stream 803-Instrument status
Function 1: Power up 2-29
Function 2: Bidirectional on 2-31
Function 3: Request instrument state 2-32
Function 4: Instrument state 2-32
Function 17: End of run 2-33
Stream 804-Setup status
Function 7: Request installed chemistries 2-34
Function 8: Installed chemistries 2-34

R

Result error codes 4-18
Descriptions 4-20

S

Sample type codes 4-3
Special calculation status codes 4-24
Stream 700-Special functions 3-4
Stream 800-Special functions 2-4

T

Test type codes 4-2

Transmission format defaults	1-7
Transmission format options	1-7

U

Unidirectional protocol	1-8
Hardware control	1-9
Protocol control characters	1-9
Software control	1-9
Transmitted information	1-8
Unit codes	4-10

V

Valid/invalid characters	2-2, 3-2
------------------------------------	----------