



# **SYNCHRON CX<sup>®</sup>** **Clinical Systems**

## **HOST COMPUTER INTERFACE SPECIFICATIONS**

**This manual is intended for use with**  
**SYNCHRON CX<sup>®</sup> 4/CX<sup>®</sup> 7 DELTA**  
**SYNCHRON CX4CE/CX7**  
**SYNCHRON CX7 RTS/CX9 ALX**  
**SYNCHRON CX<sup>®</sup> 4 PRO/CX<sup>®</sup> 5 PRO/CX<sup>®</sup> 7 PRO/CX<sup>®</sup> 9 PRO**

### **NOTICE OF DISCLAIMER**

Beckman Coulter makes no warranty as to the compatibility of any particular host computer software to the SYNCHRON CX Systems. From time to time, Beckman Coulter may change SYNCHRON CX to host interface specifications and/or SYNCHRON CX software. Beckman Coulter disclaims any and all liability and the host computer owner and/or user assumes all risk and responsibility for any and all losses, expenses and/or damages alleged to have been caused by connection or use of a host computer to or with a SYNCHRON CX System, and/or Beckman Coulter's changing the interface specifications or SYNCHRON CX software.

SYNCHRON CX System operating and performance characteristics and the instructions contained in this interface specification are based on the computer system hardware configuration specified by Beckman Coulter. Beckman Coulter disclaims all liability, including all warranties, whether express or implied, for SYNCHRON CX System performance when operated other than the specified hardware configuration.

## CX4CE/CX5CE/CX7 DELTA Operational Notes For Host Interface Specifications

The SYNCHRON DELTA System host interface is basically equivalent to the SYNCHRON CX4CE/CX5CE/CX7 interface, except for the following additions:

- There are host query timeout options of 2.5, 5, 7.5, 10 minutes and OFF. Refer to Section 4 of Host Computer Interface Specifications.
- The Define Reportable Ranges feature uses three new Result Error Codes: UH, reportable range high; UL, reportable range low; and UO, ORDAC reportable range high. These are transmitted in Stream 702, Function 3. Refer to Section 4 and Appendix F of Host Computer Interface Specifications.
- The system software will prevent the operator from programming a chemistry combination for CSF samples that will exceed the sampling capacity of the ratio pump. If the programming is done at the instrument, a pop-up window will alert the operator to this condition. If the programming is done at the host system a Return Code will be transmitted in Stream 701, Function 2. Refer to Section 4 of Host Computer Interface Specifications.
- If the system is set to program controls by reagent cartridge positions, the host will not be able to differentiate which control results apply to each reagent cartridge position if either of the following scenarios exist:
  1. Multiple reagent cartridges of the same reagent are manually loaded and the serial number is not entered; or
  2. Multiple User Defined Reagent cartridges are loaded, all with identical chemistry names.
- Alternate method of downloading sample programming that is more efficient and easier for the Host. Instead of downloading one sample programming and then releasing the line, as described in SECTION 3, *Bidirectional Protocol*, Table 3.4, the Host may download up to 7 sample programming at one time. The Host then releases the line and waits until all sample program return status messages are issued by SYNCHRON. Note that these return status messages might not be sent immediately after the line is released by the Host. Also, they might not be sent sequentially nor corresponding to the same order as the sample programming download. The bidirectional handshaking protocol and time-outs described in Section 3.1 - 3.3 still apply in this case.

Example

**SYNCHRON CX**

**HOST**  
<EOT><SOH>

|       |  |
|-------|--|
| <ACK> | < _____ >  |
|       | < _____ >  |
| <ETX> | < _____[00,701,01,Sample Programming 1...]CS<CR><LF> > |
|       | < _____ >  |
| <ACK> | < _____[00,701,01,Sample Programming 2...]CS<CR><LF> > |
|       | < _____ >  |
| <ETX> | < _____[00,701,01,Sample Programming 3...]CS<CR><LF> > |
|       | < _____ >  |
| <ACK> | < _____[00,701,01,Sample Programming 4...]CS<CR><LF> > |
|       | < _____ >  |
| <ETX> | < _____[00,701,01,Sample Programming 5...]CS<CR><LF> > |
|       | < _____ >  |
| <ACK> | < _____[00,701,01,Sample Programming 6...]CS<CR><LF> > |
|       | < _____ >  |
| <EOT> | < _____[00,701,01,Sample Programming 7...]CS<CR><LF> > |
|       | < _____ >  |

# Table of Contents

---

|                  |   |      |
|------------------|---|------|
| <b>SECTION 1</b> | Description of Interface .....  | 1-1  |
| 1.1              | Introduction .....  | 1-1  |
| 1.2              | Interface Protocol Options .....  | 1-1  |
| 1.3              | Hardware Interface .....  | 1-1  |
| 1.4              | Transmission .....  | 1-2  |
| <b>SECTION 2</b> | Unidirectional Protocol .....   | 2-1  |
| 2.1              | Description .....   | 2-1  |
| 2.2              | Software Control .....  | 2-2  |
| 2.3              | Hardware Control .....  | 2-2  |
| 2.4              | Transmission Examples .....   | 2-3  |
| <b>SECTION 3</b> | Bidirectional Protocol .....  | 3-1  |
| 3.1              | Description .....   | 3-1  |
| 3.2              | Line Bidding .....  | 3-1  |
| 3.3              | Data Transfer .....   | 3-2  |
| 3.4              | Unsolicited Messages .....  | 3-5  |
| 3.5              | Solicited Messages .....  | 3-8  |
| 3.6              | Downloading Sample Programming from Host to SYNCHRON CX .....           | 3-8  |
| <b>SECTION 4</b> | Description of Messages .....   | 4-1  |
| 4.1              | Introduction .....  | 4-1  |
| 4.2              | Message Format .....  | 4-1  |
| 4.3              | Stream 700 - Special Function .....                                     | 4-3  |
| 4.4              | Stream 701 - Sample/Cup Program .....                                   | 4-6  |
| 4.5              | Stream 702 - Results .....  | 4-13 |
| 4.6              | Stream 703 - Instrument Status .....                                    | 4-37 |
| 4.7              | Stream 704 - Setup Status .....   | 4-43 |
| <b>SECTION 5</b> | Results and Sample Programming Sequence: .....                          | 5-1  |
| 5.1              | Sample Programming .....  | 5-1  |
| 5.2              | Host Query and Sample Programming .....                                 | 5-1  |
| 5.3              | Results .....   | 5-2  |
| 5.4              | Results Recalled When Running .....                                     | 5-3  |
| 5.5              | Options for Sending CX3 Results (on CX7) And ISE Results (on CX5) ..... | 5-4  |
| 5.6              | Reagent Pack .....  | 5-5  |
| <b>SECTION 6</b> | Operator Interface - Setting Host Communications Parameters .....       | 6-1  |
| <b>SECTION 7</b> | Appendices .....  | 7-1  |

# List of Tables

| Table | Title  | Page |
|-------|--|------|
| 1.1   | Definitions . . . . .  | 1-1  |
| 1.2   | Connector Pin Assignments . . . . .                                | 1-2  |
| 2.1   | Protocol Control Characters . . . . .                              | 2-1  |
| 3.1   | Bidirectional Protocol Control Characters . . . . .                | 3-1  |
| 3.2   | Unsolicited Message Transmission Example . . . . .                 | 3-7  |
| 3.3   | Solicited Message Transmission Example . . . . .                   | 3-8  |
| 3.4   | Host Downloading Sample Programming Example . . . . .              | 3-9  |
| 4.1   | SYNCHRON CX Streams and Functions . . . . .                        | 4-2  |
| 4.2   | Stream 700 - Function 1 Are you there? . . . . .                   | 4-3  |
| 4.3   | Stream 700 - Function 2 Host Setup . . . . .                       | 4-4  |
| 4.4   | Stream 700 - Function 7 Clear Queue . . . . .                      | 4-5  |
| 4.5   | Stream 701 - Function 1 Sample/Cup Program . . . . .               | 4-6  |
| 4.6   | Stream 701 - Function 2 Sample/Cup Return Status . . . . .         | 4-9  |
| 4.7   | Stream 701 - Function 3 Clear Sector/Sample IDs . . . . .          | 4-10 |
| 4.8   | Stream 701 - Function 4 Clear Sector/Sample IDs Status . . . . .   | 4-11 |
| 4.9   | Stream 701 - Function 6 HOST QUERY Sector/Sample IDs . . . . .     | 4-12 |
| 4.10  | Stream 702 - Function 1 Cup Header . . . . .                       | 4-13 |
| 4.11  | Stream 702 - Function 3 Test Results . . . . .                     | 4-15 |
| 4.12  | Stream 702 - Function 5 End of Cup . . . . .                       | 4-18 |
| 4.13  | Stream 702 - Function 7 Linear Calibration Result . . . . .        | 4-19 |
| 4.14  | Stream 702 - Function 9 Multipoint Calibration Result . . . . .    | 4-21 |
| 4.15  | Stream 702 - Function 11 Special Calculation Result . . . . .      | 4-23 |
| 4.16  | Stream 702 - Function 13 Timed Urine Result . . . . .              | 4-24 |
| 4.17  | Stream 702 - Function 21 Reagent Pack Header . . . . .             | 4-25 |
| 4.18  | Stream 702 - Function 23 Calibration Results . . . . .             | 4-26 |
| 4.19  | Stream 702 - Function 25 Expanded Result for Calibration . . . . . | 4-28 |
| 4.20  | Stream 702 - Function 27 End of Reagent Pack . . . . .             | 4-30 |
| 4.21  | Stream 702 - Function 81 . . . . .                                 | 4-31 |
| 4.22  | Stream 702 - Function 83 Expanded Multipoint Calibration . . . . . | 4-34 |
| 4.23  | Stream 703 - Function 1 Power Up . . . . .                         | 4-37 |
| 4.24  | Stream 703 - Function 2 Bidirectional On . . . . .                 | 4-39 |
| 4.25  | Stream 703 - Function 3 Request Instrument State . . . . .         | 4-40 |
| 4.26  | Stream 703 - Function 4 Instrument State . . . . .                 | 4-40 |
| 4.27  | Stream 703 - Function 5 Instrument Exception . . . . .             | 4-41 |
| 4.28  | Stream 703 - Function 7 Chemistry Configuration Change . . . . .   | 4-41 |
| 4.29  | Stream 703 - Function 13 Range Change . . . . .                    | 4-42 |
| 4.30  | Stream 703 - Function 17 End of Run . . . . .                      | 4-42 |
| 4.31  | Stream 704 - Function 7 Request Installed Chemistries . . . . .    | 4-43 |
| 4.32  | Stream 704 - Function 8 Installed Chemistries . . . . .            | 4-43 |
| 4.33  | Stream 704 - Function 9 Request Chemistry Ranges . . . . .         | 4-44 |
| 4.34  | Stream 704 - Function 10 Chemistry Ranges . . . . .                | 4-45 |
| 6.1   | Host Communication Parameter Setup Options . . . . .               | 6-3  |

# SECTION 1 Description of Interface

## 1.1 Introduction

These instructions provide the necessary information to interface the SYNCHRON CX3, CX4, CX5, CX7 DELTA, CX4CE, CX5CE, CX7 Systems, and CX PRO Systems to a laboratory computer system.

Table 1.1 lists several definitions that are useful in understanding this interface specification.

Table 1.1 Definitions

|                              |   |
|------------------------------|---|
| <b>SYNCHRON CX</b>           | SYNCHRON CX3, CX4, CX5, CX7 DELTA, and CX4CE, CX5CE, CX7 Systems, and CX PRO Systems.   |
| <b>Host</b>                  | Clinical laboratory computer system.  |
| <b>Download<sup>a</sup></b>  | The process of a host sending sample/cup programming to the SYNCHRON CX System.   |
| <b>Handshake<sup>a</sup></b> | The process of the SYNCHRON CX System and the host communicating with each other by requesting permission to transmit, granting permission, and acknowledging receipt of transmission (ACK - NAK protocol). |

<sup>a</sup> Bidirectional only.

## 1.2 Interface Protocol Options

SYNCHRON CX Systems provide two interface protocol options which are selected using the Define Host Screen:

- UNIDIRECTIONAL** Information is transmitted from the SYNCHRON CX System to the host computer. An X-ON, X-OFF protocol can be used by the host to control the information flow. (Section 2)
- BIDIRECTIONAL** Information is transmitted from the SYNCHRON CX System to the host computer and from the host computer to the SYNCHRON CX System. An X-ON, X-OFF protocol can be used by the host to control the information flow. The two systems communicate with each other using a software handshake, requesting permission to send information and acknowledging receipt of information (ACK - NAK protocol). (Section 3)

## 1.3 Hardware Interface

The hardware interface to the SYNCHRON CX Systems external communication port uses a standard RS-232-C 9-pin or 25-pin D-connector (male). Communication is done with asynchronous byte transmission. The host communication port is located on the back of the computer console.

### 1.3.1 Connector Pin Assignments

Connector pin assignments are shown in Table 1.2.

Table 1.2 Connector Pin Assignments

| 25-PIN          | 9-PIN          | Signal                                   | Direction |
|-----------------|----------------|--|-----------|
| 1               | ---            | Protected ground - used on one side only |           |
| 2               | 3              | Transmit data                            | Out       |
| 3               | 2              | Receive Data                             | In        |
| 4 <sup>a</sup>  | 8 <sup>a</sup> | Clear to send (CTS)                      | In        |
| 5 <sup>b</sup>  | 7 <sup>b</sup> | Request to send (RTS)                    | Out       |
| 6 <sup>c</sup>  | 4 <sup>c</sup> | Data terminal ready                      | Out       |
| 7               | 5              | Signal ground                            |           |
| 20 <sup>c</sup> | 6 <sup>c</sup> | Data set ready                           | In        |
| ---             | 1              | Data Carrier Detect                      | In        |
| ---             | 9              | Ring Indicator                           | In        |

<sup>a</sup> CTS is required before data will be transmitted, when using hardware flow control.

<sup>b</sup> RTS is turned on when operational, when using hardware flow control.

<sup>c</sup> Pins 6 and 20 (25-pin), and 4 and 6 (9-pin) are only needed when a modem is used.

### 1.3.2 Cable Length

The recommended cable length when using the RS-232-C interface should not exceed 40 feet (12 m).

## 1.4 Transmission

The SYNCHRON CX Systems use 7 bit ASCII (provides the 128 character ASCII set when 7 data bits are selected). For local languages, 8 bit ASCII must be used (provides the 256 character ASCII set).

### 1.4.1 Transmission Format

The serial transmission format is configured from the SYNCHRON CX System console using the Define Host Screen. The following options are available:

Interface Options: Unidirectional, Bidirectional or None

Data Bits: 7 or 8\*

Parity: Even, Odd or None

Stop Bits: 1 or 2

Baud Rates: 600, 1200, 2400, 4800 or 9600

\* For use with non-English keyboards, select 8-bit ASCII for the host interface in order to send and receive special characters (for key code conversion information, refer to Appendix J).

### 1.4.2 Data Record Format

The data record format, in both the unidirectional and bidirectional modes, used for all transmissions to and from the SYNCHRON CX System is:

' [ ' <MESSAGE> ' ] ' <CS> <CR> <LF>

where:

- ' [ '            Is the start of message character (HEX = 5B).
- <MESSAGE>    Is the data to be transmitted (refer to Section 4).
- ' ] '            Is the end of message character (HEX = 5D).
- <CS>            <CS> is the checksum. It is two hexadecimal digits.  
                 256 - ((Sum of ASCII value of all characters in <MESSAGE>,  
                 including spaces, commas, ' [ ' and ' ] ' ) MODULE 256) = CS
- <CR>            Is an ASCII carriage return.
- <LF>



## SECTION 2 Unidirectional Protocol

### 2.1 Description

In unidirectional transmission, information is transmitted in one direction only - from the SYNCHRON CX System to the host computer. The information flow can be controlled by the host through software protocols with software versions prior to 1.7. The information flow can be controlled by the host through software or hardware protocols with software version 1.7 and above (Sections 2.2 and 2.3). This information is transmitted as soon as it is available. The information transmitted includes:

|                            |   |
|----------------------------|---|
| <b>CUP HEADER</b>          | When a cup has finished and results are ready to send to the host, the cup header information is transmitted. (Stream 702 - Function 1)   |
| <b>REAGENT PACK HEADER</b> | The reagent pack header contains reagent and chemistry information and is sent before sending any calibration results. (Stream 702, Function 21)  |
| <b>TEST RESULTS</b>        | Test results are sent as a group when all the tests in the cup are completed. CX3 results (on CX7) and ISE results (on CX5) may be sent before all tests are completed if CX3/ISE immediate output is enabled (refer to Paragraph 6.5.6). An expanded results record is transmitted when the expanded result option is enabled.<br>(Stream 702 - Functions 3, 7, 9, 11, 13, 23, 25, 81, 83) |
| <b>END OF CUP</b>          | When all tests programmed for a sample are completed and sent, an end of cup message is transmitted. (Stream 702 - Function 5)  |
| <b>END OF REAGENT PACK</b> | A reagent pack end is sent after all the calibration results<br>(Stream 702, Function 27).  |
| <b>END OF RUN</b>          | When all tests programmed have been completed and the SYNCHRON CX System has gone into the standby mode an end of run message is transmitted. (Stream 703 - Function 17)  |
| <b>HOST SETUP CHANGES</b>  | When any options in the Define Host Screen are changed a host setup change message is transmitted. (Stream 700 - Function 2)  |
| <b>POWER UP</b>            | Transmitted when the SYNCHRON CX System is booted.<br>(Stream 703 - Function 1)   |

#### NOTICE

Output field width is dependent on the units and decimal precision selected. This should not exceed the fixed field width specified in this document.

The same transmission format is used in both the unidirectional and bidirectional modes. Refer to Section 4 for the details of the transmitted messages.

Table 2.1 Protocol Control Characters

| Character | Name                 | ASCII | HEX |
|-----------|----------------------|-------|-----|
| X-ON      | Resume transmission  | DL1   | 11  |
| X-OFF     | Suspend transmission | DL3   | 13  |

## 2.2 Software Control

The flow of information from the SYNCHRON CX Systems can be controlled by the host through an X-ON, X-OFF protocol (Table 2.1). The host transmits the character X-OFF prior to its buffer overflowing. The SYNCHRON CX suspends transmission after a few characters. When the host's buffer is ready to accept information again, the character X-ON is transmitted. The SYNCHRON CX then resumes transmission. Also, the host should transmit the X-ON character after power up to indicate it is ready to receive data.



### **WARNING - Lost Data**

**After the host transmits X-OFF, the SYNCHRON CX stores the data to be transmitted in a buffer. If the host does not transmit X-ON before the SYNCHRON CX buffer overflows, the overflow data will be lost. The size of the buffer is 600 samples.**

## 2.3 Hardware Control

The flow of information from the SYNCHRON CX Systems can be controlled by the host through the CTS hardware signal (refer to Paragraph 1.3.1). When the CTS signal is brought negative, the SYNCHRON CX will suspend transmission. The SYNCHRON CX will resume transmission when the CTS signal is brought positive.



### **WARNING - Lost Data**

**When the CTS signal is brought negative, the SYNCHRON CX stores the data to be transmitted in a buffer. If the CTS signal is not brought positive before the SYNCHRON CX buffer overflows, the overflow data will be lost. The size of the buffer is 600 samples.**

## 2.4 Transmission Examples

### EXAMPLE 1:

Unidirectional mode is selected and the host is not using software control of transmission. The accession number is 1100, sector 1, cup 3, and the selected chemistries are creatinine (03A), chloride (04A), potassium (01B), and sodium (01A). The special calculation anion gap with potassium (K) is also transmitted for the sample.

```
[- 0,702,01,270291,113121,- 1100,RG,- 1,- 3,RO,#####,TU,SAMPLE1.01- ,
-----,
Lee-----,John-----,- ,39485-----,Nelson-----,270291,
1044,304-----,- 39,5-----,M,-----,1400.00,
24.0,2.30,1.7300,- - 4,03A- ,04A- ,01B- ,01A- J35<CR><LF>
[- 0,702,03,270291,114132,- 1100,- - - - 3442,- 1,- 3,SAMPLE1.01- ,03A- ,2sk,911150,
19,- 1,#####,#####2,1,- 0,NA,NR,NR,0,NA,#####,- - - - - ,DL,NO,NO,
NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,1.0000,#####]
9D<CR><LF>
[- 0,702,03,270291,113741,- 1100,- - - - 3430,- 1,- 3,SAMPLE1.01- ,04A- ,###,#####,
##,- 1,- - - - 123.9,#####2,0,- 4,NA,NR,NA,0,NA,123.91887,- - - - - ,NO,NO,NO,
NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,1.0000,#####]
E9<CR><LF>
[- 0,702,03,270291,113741,- 1100,- - - - 3429,- 1,- 3,SAMPLE1.01- ,01B- ,###,#####,
##,- 1,- - - - 3.60,#####2,0,- 4,NA,NR,NA,0,NA,3.6003411,- - - - - ,NO,NO,NO,
NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,1.0000,#####]
0E<CR><LF>
[- 0,702,03,270291,113741,- 1100,- - - - 3428,- 1,- 3,SAMPLE1.01- ,01A- ,###,#####,
##,- 1,- - - - 174.3,#####2,0,- 4,NA,NR,NA,0,NA,174.28595,- - - - - ,NO,NO,NO,
NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,1.0000,#####]
E3<CR><LF>
[- 0,702,11,270291,114148,-1100,-1,- 3,SAMPLE1.01- ,-1,USER- SPL- CALC- - - - - ,OK,
120.31853,UN/UN- - - J4B<CR><LF>
[- 0,702,11,270291,114148,-1100,-1,- 3,SAMPLE1.01- ,-1,USER- SPL- CALC2- - - - - ,OK,
50.367081,UNIT- - - J47<CR><LF>
[- 0,702,13,270291,114148,- 1100,-1,- 3,SAMPLE1.01- ,-1,CREA- - - - - ,AB,
0.0000000,- - - - - J92<CR><LF>
[- 0,702,13,270291,114148,- 1100,-1,- 3,SAMPLE1.01- ,-1,CL- - - - - ,OK,
173.48641,mmol/24.J2D<CR><LF>
[- 0,702,13,270291,114148,- 1100,-1,- 3,SAMPLE1.01- ,-1,K- - - - - ,OK,
5.0404774,mmol/24.J54<CR><LF>
[- 0,702,13,270291,114148,- 1100,-1,- 3,SAMPLE1.01- ,-1,NA- - - - - ,OK,
244.00033,mmol/24.J3F<CR><LF>
[- 0,702,05,270291,114148,- 1100,SAMPLE1.01- ,-1,- 3J0A<CR><LF>
[- 0,703,17,190291,103626]DE
```

(- indicates space holder)

The same information is transmitted as in EXAMPLE 1, but software transmission control is implemented by the host.

## HOST

```

===== X - ON
[- 0,702,01,270291,113121,- 1100,RG,- 1,- 3,RO,#####,TU,SAMPLE1.01- ,
-----,-----,
-----,Lee-----,John-----,- ,
39485- -----,Nelson-----,270291,1044,304- -----,
- 39,5,- -----,M,- -----,1400.00,24.0,2.30,
1.7300,- - 4,03A- ,04A- ,01B- ,01A- ]35<CR><LF>

<===== X - OFF
<===== X - ON
[- 0,702,03,270291,114132,- 1100,- - - - 3442,- 1,- 3,SAMPLE1.01- ,03A- ,
2sk,911150,19,- 1,#####,#####,2,1,- 0,NA,NR,NR,0,NA,
#####,- - - - - ,DL,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,
NO,NO,1.0000,#####]9D<CR><LF>
[- 0,702,03,270291,113741,- 1100,- - - - 3430,- 1,- 3,SAMPLE1.01- ,04A- ,
###,#####,##,- 1,- - - 123.9,#####,2,0,- 4,NA,NR,NA,0,NA,
123.91887,- - - - - ,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,
NO,NO,1.0000,#####]E9<CR><LF>
[- 0,702,03,270291,113741,- 1100,- - - - 3429,- 1,- 3,SAMPLE1.01- ,01B- ,
###,#####,##,- 1,- - - - 3.60,#####,2,0,- 4,NA,NR,NA,0,NA,
3.6003411,- - - - - ,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,
NO,NO,1.0000,#####]0E<CR><LF>
[- 0,702,03,270291,113741,- 1100,- - - - 3428,- 1,- 3,SAMPLE1.01- ,01A- ,
###,#####,##,- 1,- - - 174.3,#####,2,0,- 4,NA,NR,NA,0,NA,
174.28595,- - - - - ,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,
NO,NO,1.0000,#####]E3<CR><LF>
[- 0,702,11,270291,114148,- 1100,- 1,- 3,SAMPLE1.01- ,- 1,
USER- SPL- CALC- - - - - ,OK,120.31853,UN/UN- - - ]4B<CR><LF>
[- 0,702,11,270291,114148,- 1100,- 1,- 3,SAMPLE1.01- ,- 1,
USER- SPL- CALC2- - - - - ,OK,50.367081,UNIT- - - ]47<CR><LF>
[- 0,702,13,270291,114148,- 1100,- 1,- 3,SAMPLE1.01- ,- 1,
CREA- - - - - ,AB,0.0000000,- - - - - ]92<CR><LF>
[- 0,702,13,270291,114148,- 1100,- 1,- 3,SAMPLE1.01- ,- 1,
CL- - - - - ,OK,173.48641,mmol/24.]2D<CR><LF>
[- 0,702,13,270291,114148,- 1100,- 1,- 3,SAMPLE1.01- ,- 1,
K- - - - - ,OK,5.0404774,mmol/24.]54<CR><LF>
[- 0,702,13,270291,114148,- 1100,- 1,- 3,SAMPLE1.01- ,- 1,
NA- - - - - ,OK,244.00033,mmol/24.]3F<CR><LF>
[- 0,702,05,270291,114148,- 1100,SAMPLE1.01- ,- 1,- 3]0A<CR><LF>
[- 0,703,17,190291,103626]DE

```

(- indicates space holder)

## SECTION 3 Bidirectional Protocol

### 3.1 Description

In bidirectional transmission, information is transmitted in two directions - from the SYNCHRON CX System to the host and from the host to the SYNCHRON CX System. The two systems communicate with each other using a software handshake, requesting permission to send information, granting permission, and acknowledging receipt of information (ACK - NAK protocol). Table 3.1 describes the bidirectional protocol control characters that are referred to throughout this section.

Table 3.1 Bidirectional Protocol Control Characters


| Acronym | Name   | ASCII | HEX |
|---------|--|-------|-----|
| LB      | Line Bid                                     | SOH   | 01  |
| LBO     | Line Bid Override Request                    | STX   | 02  |
| LF      | Line Feed                                    | LF    | 0A  |
| ENQ     | Enquiry (used to reestablish communications) | ENQ   | 05  |
| EOT     | End Of Transmission                          | EOT   | 04  |
| ACK-0   | Even Acknowledgement                         | ACK   | 06  |
| ACK-1   | Odd Acknowledgement                          | ETX   | 03  |
| NAK     | Negative Acknowledgement                     | NAK   | 15  |

### 3.2 Line Bidding


Before either the SYNCHRON CX System or the host transmit a message, they must bid for the communication line. Line bidding uses the LB, LBO, EOT, ACK-0, and NAK characters.

To bid for the line, <EOT><SOH> is transmitted when the line is idle. The EOT clears the line. The receiving system responds by transmitting ACK to acknowledge the line bid and allow message transmission or NAK to deny the line bid. If the line bid is denied, the originator should wait a short period of time (e.g. one second) before bidding for the line again.

#### SUCCESSFUL LINE BID



#### UNSUCCESSFUL LINE BID



The examples above are also correct if the labels SYNCHRON CX and HOST are reversed.

LINE BID - When the SYNCHRON CX is very busy, the response to a LINE BID may take up to a maximum of 15 seconds. The host may have to wait up to 15 seconds for the LINE BID response.

When the SYNCHRON CX has some messages in the output queue, it does a LINE BID. When the LINE BID is accepted by the host, the instrument sends one set of messages at a time, then gives up the line and waits for 2 seconds, then repeats the process if other messages are to be sent.

If the host has messages to transmit, it can do a LINE BID (EOT SOH) even if the instrument is busy transmitting messages. After sending a message, the SYNCHRON CX looks for a NAK, ACK, or SOH. When SYNCHRON CX sends a message, host can respond with a line bid <EOT> <SOH>. The message in progress will be saved by SYNCHRON CX and the line will be granted to the host by sending an <ACK>. This is not recommended because it will stop the collation of results.

LINE BID TIME OUT occurs if the receiving system does not respond to the EOT LB within fifteen (15) seconds. After seven (7) consecutive time outs or unrecognizable responses (i.e. not ACK or NAK), the originator considers the line "DOWN". The originator waits twenty (20) seconds and tries again.

CONTENTION occurs when both systems bid for the line at the same time. The SYNCHRON CX will be considered the master and the host should respond with ACK. However, the host may override the SYNCHRON CX line bid by transmitting LBO in response. The SYNCHRON CX will respond by transmitting ACK to acknowledge the line bid override and allow message transmission or will ignore the message if SYNCHRON CX was not bidding for the line when the message was sent.

### 3.3 Data Transfer

After successfully bidding for the line (Paragraph 3.2), the originator or sending system transmits its message. The message format is described in Paragraph 1.4.2 and Section 4.

A successful data transfer consists of:

[ <MESSAGE> ]<CS><CR><LF>

If the checksum is incorrect or any element is missing, the data transfer is unsuccessful.

After data transmission, the receiving system acknowledges transmission as follows:


1. If data transfer was successful, then the receiving system alternately returns ETX and ACK after each message. Since the receiving system responded ACK to the line bid, ETX is the correct acknowledgement to the first message, ACK to the second message, then ETX again, etc. When data transfer is complete, the sending system transmits EOT.
2. If data transfer was unsuccessful due to a bad checksum or other problem, then the receiving system responds NAK. The sending system retransmits the message up to seven times. If unsuccessful, the sending system rebids for the line and retransmits the last message.
3. If the sending system does not receive an acknowledgement within fifteen (15) seconds after data transfer or an incorrect acknowledgement is received, it transmits ENQ. The receiving system retransmits its last acknowledgement (e.g., ACK, ETX, NAK, or EOT). The sending system responds to the acknowledgement in one of the following manners:
  - If the correct acknowledgement is transmitted, the sending system transmits the next message or EOT.
  - If an incorrect acknowledgement or NAK is transmitted, the sending system retransmits the last message.
  - If the receiving system transmits an EOT, the sending system rebids for the line and retransmits the last message.
4. TIME OUT occurs:
  - If the sending system does not transmit data in twenty (20) seconds after successfully bidding for the line (Paragraph 3.2) the receiving system times out and returns to idle.
  - If the sending system does not receive an acknowledgement within fifteen (15) seconds after data transfer or an incorrect acknowledgement is received, it transmits ENQ. If an acknowledgement is not received after seven (7) consecutive ENQ transmissions, the SYNCHRON CX waits twenty (20) seconds and attempts a line bid and will continue to do so until an acknowledgement is received.

5. IDLE state occurs:


- The receiving system requests an IDLE state by transmitting an EOT. The sending system will retransmit the interrupted message.

Examples of data transfer (These examples are also correct if the labels SYNCHRON CX and HOST are reversed.):


#### NAK Example




#### No Response



### No Response Return to Idle



### Request for Idle State





### 3.4 Unsolicited Messages

Unsolicited messages are those messages which are automatically transmitted by the SYNCHRON CX System as the information becomes available. The following unsolicited messages are transmitted in the bidirectional mode:

|                            |   |
|----------------------------|---|
| <b>REAGENT PACK HEADER</b> | The reagent pack header contains reagent and chemistry information, and is sent before sending any calibration results.<br>(Stream 702 - Function 21).                              |
| <b>REAGENT RESULTS</b>     | Calibration results are sent as a group when all tests in the cup are completed.  |
| <b>END OF REAGENT</b>      | A reagent pack end is sent after all calibration results and indicating that all information for that reagent pack is complete.<br>(Stream 702 - Function 7, 9, 21, 23, 25, 27, 83) |

#### NOTICE

Output field width is dependent on the units and decimal precision selected. This should not exceed the fixed field width specified in this document.

|  |   |
|--|---|
| <b>CUP HEADER</b>  | When a cup has finished and results are ready to send to the host, the cup header information is transmitted.<br>(Stream 702 - Function 1)  |
| <b>TEST RESULTS</b>  | Test results are sent as a group when all the tests in the cup are completed. CX3 results (on CX7) and ISE results (on CX5) may be sent before all tests are completed if CX3/ISE immediate output is enabled (refer to Paragraph 6.5.6). An expanded results record is transmitted when the expanded result option is enabled.<br>(Stream 702 - Functions 3, 7, 9, 11, 13, 23, 25, 81, 83) |
| <b>END OF CUP</b>  | When all tests programmed for a sample are completed and sent, an end of cup message is transmitted.<br>(Stream 702 - Function 1, 3, 5, 11, 81)   |
| <b>END OF RUN</b>  | When all tests programmed have been completed and the SYNCHRON CX has gone into the idle mode an end of run message is transmitted. (Stream 703 - Function 17)  |
| <b>NORMAL/CRITICAL RANGES CHANGED</b><br><b>OR</b><br><b>CHEMISTRY CONFIGURATION CHANGED</b> | When any normal and/or critical range is changed, or the chemistry configuration is changed the appropriate message is transmitted.<br>(Stream 703 - Functions 7, 13)   |
| <b>HOST SETUP CHANGES</b>  | When any options in the Define Host Screen are changed a host setup change message is transmitted. (Stream 700 - Function 2)  |
| <b>POWER UP</b>  | Transmitted when the SYNCHRON CX is booted.<br>(Stream 703 - Function 1)  |

#### NOTICE

Undefined characters may be transmitted during system boot-up.

**BIDIRECTIONAL STARTUP**

Transmitted when the bidirectional interface option is enabled.  
(Stream 703 - Function 2)

Refer to Section 4 for the details of the transmitted messages and Table 3.2 for an example of SYNCHRON CX unsolicited message transmission.

Table 3.2 Unsolicited Message Transmission Example

| SYNCHRON CX  | HOST |
|--|------|
| EOT SOH =====>   |      |
| <=====ACK  |      |
| [-0,702,01,270291,113121,-1100,RG,-1,-3,RO,#####,TU,SAMPLE1.01-, |      |
| -----,-----,   |      |
| -----,Lee-----,John-----,-,                                      |      |
| 39485-----,Nelson-----,270291,1044,304-----,                     |      |
| -39,5,-----,M,-----,1400.00,24.0,2.30,1.7300,                    |      |
| --4,03A-,04A-,01B-,01A-]35<CR><LF>                               |      |
| <=====ETX  |      |
| -0,702,03,270291,114132,-1100,----3442,-1,-3,SAMPLE1.01-,03A-,   |      |
| 2sk,911150,19,-1,#####,#####2,1,-0,NA,NR,NR,0,NA,                |      |
| #####-----,DL,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,               |      |
| NO,NO,1.0000,#####]9D<CR><LF>                                    |      |
| <=====ACK  |      |
| [-0,702,03,270291,113741,-1100,----3430,-1,-3,SAMPLE1.01-,04A-   |      |
| ###,#####,-1,----123.9,#####2,0,-4,NA,NR,NA,0,NA,                |      |
| 123.91887,-----,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,          |      |
| NO,NO,1.0000,#####]E9<CR><LF>                                    |      |
| <=====ETX  |      |
| [-0,702,03,270291,113741,-1100,----3429,-1,-3,SAMPLE1.01-,01B-,  |      |
| ###,#####,-1,----3.60,#####2,0,-4,NA,NR,NA,0,NA,                 |      |
| 3.6003411,-----,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,          |      |
| NO,NO,1.0000,#####]0E<CR><LF>                                    |      |
| <=====ACK  |      |
| [-0,702,03,270291,113741,-1100,----3428,-1,-3,SAMPLE1.01-,01A-,  |      |
| ###,#####,-1,----174.3,#####2,0,-4,NA,NR,NA,0,NA,                |      |
| 174.28595,-----,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,          |      |
| NO,NO,1.0000,#####]E3<CR><LF>                                    |      |
| <=====ETX  |      |
| [-0,702,11,270291,114148,-1100,-1,-3,SAMPLE1.01-,-1,             |      |
| USER-SPL-CALC-----,OK,120.31853,UN/UN---]4B<CR><LF>              |      |
| <=====ACK  |      |
| [-0,702,11,270291,114148,-1100,-1,-3,SAMPLE1.01-,-1,             |      |
| USER-SPL-CALC2-----,OK,50.367081,UNIT---]47<CR><LF>              |      |
| <=====ETX  |      |
| [-0,702,13,270291,114148,-1100,-1,-3,SAMPLE1.01-,-1,             |      |
| CREA-----,AB,0.0000000,-----]92<CR><LF>                          |      |
| <=====ACK  |      |
| [-0,702,13,270291,114148,-1100,-1,-3,SAMPLE1.01-,-1,             |      |
| CL-----,OK,173.48641,mmol24.]2D<CR><LF>                          |      |
| <=====ETX  |      |
| [-0,702,13,270291,114148,-1100,-1,-3,SAMPLE1.01-,-1,             |      |
| K-----,OK,5.0404774,mmol24.]54<CR><LF>                           |      |
| <=====ACK  |      |
| [-0,702,13,270291,114148,-1100,-1,-3,SAMPLE1.01-,-1,             |      |
| NA-----,OK,244.00033,mmol24.]3F<CR><LF>                          |      |
| <=====ETX  |      |
| [-0,702,05,270291,114148,-1100,SAMPLE1.01-,-1,-3]0A<CR><LF>      |      |
| <=====ACK  |      |
| [-0,703,17,190291,103626]DE                                      |      |
| <=====ETX  |      |
| EOT=====>  |      |

(-indicates space holder)

### 3.5 Solicited Messages

Solicited messages are transmitted by the host to request information from the SYNCHRON CX System. Solicited messages are used by the host to:

- Identify that a bidirectional system is attached to the host serial port. (Stream 700 - Function 1)
- Obtain SYNCHRON CX instrument status information. (Stream 703 - Function 3)
- Obtain a list of the SYNCHRON CX installed chemistries and their corresponding defined normal and critical ranges. (Stream 704 - Functions 7, 9)

Refer to Section 4 for the details of the transmitted messages and Table 3.3 for an example of host solicited message transmission.

Table 3.3 Solicited Message Transmission Example

| SYNCHRON CX |  | HOST    |
|-------------|--|---------|
|             | < _____  | EOT SOH |
| ACK         | < _____>   |         |
|             | < _____ [00,703,03] <CS><CR><LF>                               |         |
| ETX         | < _____  |         |
|             | < _____  | EOT     |
| EOT SOH     | < _____>   |         |
|             | < _____  | ACK     |
|             | [00,703,04,120888,093044,CX7--,01,01,##,01]<CS><CR><LF> _____> |         |
|             | < _____  | ETX     |
| EOT         | _____>   |         |

(<CS> indicates checksum.)

(-indicates space holder)

### 3.6 Downloading Sample Programming from Host to SYNCHRON CX

Sample programming may be transmitted directly from the host to the SYNCHRON CX System. The information which can be programmed includes sector and cup number, sample ID, test type, fluid type, demographics, and chemistry requests (Stream 701 - Function 1). The dilution factor cannot be downloaded from the host. It must be edited at the SYNCHRON CX System before the sample is recognized on the sample carousel. The response by the SYNCHRON CX is used to notify the host whether or not the program was accepted (e.g. a BUSY response may indicate programming is being done at the SYNCHRON CX console). There is an interlock prohibiting simultaneous programming of a single sample from the host and the SYNCHRON CX console. After the SYNCHRON CX releases the line (sends an EOT), the host may then repeat the process for as many sample cups as necessary for completion of the download.

In addition, the host can clear previous sector programming before transmitting the new sample programming (Stream 701 - Function 3). In response to the clear sector message transmitted by the host, the SYNCHRON CX transmits an OK, BUSY or SYNTAX ERROR message (Stream 701 - Function 4). This response by the SYNCHRON CX is used to notify the host whether or not the sector was cleared (e.g. A BUSY response may indicate programming is being done at the SYNCHRON CX console or the sector is on the sample wheel.) There is an interlock prohibiting simultaneous clearing from the host and clearing or programming from the SYNCHRON CX console.

Refer to Section 4 for the details of the transmitted messages and Table 3.4 for an example of host downloading.

Table 3.4 Host Downloading Sample Programming Example

| SYNCHRON CX |   | HOST    |
|-------------|---|---------|
| CK          | <=====  | EOT SOH |
|             | =====   |         |
|             | <===== [00,701,01,-2,-1,1,ST,SE,235-----,               |         |
|             | -----, FASTING-----,                                    |         |
|             | LIPEMIC-----,JOHNSON-----,HAROLD-----,-,                |         |
|             | -39-44-6207-,WASHINGTON-----,110391,1130,CCU-2-----,    |         |
|             | 035,5,170852,M,PATIENT-IS-DIABETIC-----,-----,-----,    |         |
|             | 004,01A-,0,01B-,0,04A-,0,02A-,0]<CS><CR><LF>            |         |
| ETX         | =====   |         |
|             | <=====  | EOT     |
| EOT SOH     | =====   |         |
|             | <=====  | ACK     |
|             | [-0,701,02,-0,-230,-2,-1,235-----]<CS><CR><LF>          |         |
|             | =====   |         |
|             | <=====  | ETX     |
| EOT         | =====   |         |
|             | <=====  | EOT SOH |
| ACK         | =====   |         |
|             | <===== [00,701,01,02,02,0,RO,SE,238-----,               |         |
|             | -----,  |         |
|             | -----,BROWN-----,DOROTHY-----,                          |         |
|             | -,480-22-0894-,JONES-----,110391,1115,                  |         |
|             | 202-----,060,5,-----,F,-----,                           |         |
|             | -----,-----,-----,003,03A-,0,05A-,0,06A-,1]<CS><CR><LF> |         |
| ETX         | =====   |         |
|             | <=====  | EOT     |
| EOT SOH     | =====   |         |
|             | <=====  | ACK     |
|             | [-0,701,02,-0,-231,-2,-2,238-----]<CS><CR><LF>          |         |
|             | =====   |         |
|             | <=====  | ETX     |
| EOT         | =====   |         |
|             | <=====  | EOT SOH |
| ACK         | =====   |         |
|             | <===== [00,701,01,02,03,0,RO,UR,239-----,               |         |
|             | -----,  |         |
|             | -----,SMITH-----,SAM-----,                              |         |
|             | -,123-45-6789-,WASHINGTON-----,110891,1000,             |         |
|             | 445-1-----,028,5,-----,M,-----,                         |         |
|             | -----,-----,-----,004,03A-,0,01B-,0,04A-,0,06A-,0]      |         |
|             | <CS><CR><LF>  |         |
| ETX         | =====   |         |
|             | <=====  | EOT     |
| EOT SOH     | =====   |         |
|             | <=====  | ACK     |
|             | [-0,701,02,-0,-232,-2,-3,239-----]<CS><CR><LF>=====     |         |
|             | <=====  | ETX     |
| EOT         | =====   |         |

(<CS> indicates checksum, -indicates space holder)

## SECTION 4 Description of Messages

### 4.1 Introduction

The data record format for all messages transmitted to and from the SYNCHRON CX System is:

' [ '<MESSAGE>' ] '<CS><CR><LF>

This Section describes the information contained in the <MESSAGE> portion. Refer to Paragraph 1.4.2 for a description of ' [ ', ' ] ', <CS>, <CR> and <LF>.

### 4.2 Message Format

The format of the <MESSAGE> is:

<DEVICE ID>,<STREAM>,<FUNCTION>,<FIELD1>,...,<FIELDN>

where:

- <DEVICE ID> is a number between 0 and 99. The default SYNCHRON CX ID is 0 but can be changed using the Define Host Screen. All messages which have a different Device ID from the Define Host Screen will not be processed.
- <STREAM> is a number between 700 and 799 (refer to Table 4.1).
- <FUNCTION> is a number between 1 and 99 (refer to Table 4.1).
- <FIELD1>,...,<FIELDN> are data fields associated with the <STREAM> and <FUNCTION>.

#### NOTICE

All data fields (<FIELDN>) are fixed length. Numeric fields, excluding function numbers, are right justified and blank filled to the maximum length. Character and string fields are left justified and right blank filled to the maximum length. If numeric data exceeds the maximum field length, the field will be filled with asterisks (\*). If a field does not apply in a record it will be filled with pound signs (#).

Alphanumeric fields received from the host cannot contain commas(.). Commas received in an alphanumeric field will generate an HCP error. A comma entered in an alphanumeric field at the instrument is transmitted as a semicolon (;). A semicolon received from the host is converted to a comma. Acceptable characters for Sample ID ASCII characters are 33-126 (except , ; \$ \* ? [ ] \ and ^). For other entries, ASCII characters 32-175 are acceptable except 127 and 44.

Messages sent to and from the SYNCHRON CX are divided into streams. Each stream corresponds to one group of related operations within the interface. Within each stream are one or more functions. For each stream only one function can be active at a time. However, multiple streams can be active at the same time. The SYNCHRON CX's streams and associated functions are described in Table 4.1. Note that only the items with an asterisk (\*) are used by the unidirectional interface.

Stream transmission options are selected through the Special Functions: 4. System Setup, 10. Host Communications Parameters.

## SYNCHRON CX Streams and Functions

Table 4.1 SYNCHRON CX Streams and Functions

| Operation         | Stream | Function        | Information                                      | Sent By                  |
|-------------------|--------|-----------------|--|--------------------------|
| Special Functions | 700    | 1               | Are you there?                                   | Host                     |
|                   |        | 2               | Host Setup                                       | SYNCHRON CX (U)          |
|                   |        | 7               | Protocol control message (Clear Queue)           | Host                     |
| SampleCup Program | 701    | 1               | Sample/cup program                               | Host                     |
|                   |        | 2               | Sample/cup return status                         | SYNCHRON CX(S)           |
|                   |        | 3               | Clear sector/sample IDs                          | Host                     |
|                   |        | 4               | Clear sector/sample IDs                          | SYNCHRON CX(S)<br>status |
|                   |        | 6               | Host Query sector/sample IDs<br>Auto Clear Queue | SYNCHRON CX (U)          |
| Results           | 702    | 1 <sup>a</sup>  | Cup header                                       | SYNCHRON CX(U)           |
|                   |        | 3 <sup>a</sup>  | Test result                                      | SYNCHRON CX(U)           |
|                   |        | 5 <sup>a</sup>  | End of cup.                                      | SYNCHRON CX(U)           |
|                   |        | 7 <sup>a</sup>  | Linear calibration result.                       | SYNCHRON CX(U)           |
|                   |        | 9 <sup>a</sup>  | Multipoint calibration result.                   | SYNCHRON CX(U)           |
|                   |        | 11 <sup>a</sup> | Special calculation result.                      | SYNCHRON CX(U)           |
|                   |        | 13 <sup>a</sup> | Timed Urine Result                               | SYNCHRON CX(U)           |
|                   |        | 21 <sup>a</sup> | Reagent pack header.                             | SYNCHRON CX(U)           |
|                   |        | 23 <sup>a</sup> | Results for calibration.                         | SYNCHRON CX(U)           |
|                   |        | 25 <sup>a</sup> | Expanded result for calibration.                 | SYNCHRON CX(U)           |
|                   |        | 27 <sup>a</sup> | End of Reagent Pack                              | SYNCHRON CX(U)           |
|                   |        | 81 <sup>a</sup> | Expanded result.                                 | SYNCHRON CX(U)           |
|                   |        | 83 <sup>a</sup> | Expanded multipoint calibration<br>result.       | SYNCHRON CX(U)           |
| Instrument Status | 703    | 1 <sup>a</sup>  | Power up.  | SYNCHRON CX(U)           |
|                   |        | 2               | Bidirectional interface.                         | SYNCHRON CX(U)           |
|                   |        | 3               | Request instrument status                        | Host                     |
|                   |        | 4               | Instrument status.                               | SYNCHRON CX(S)           |
|                   |        | 5 <sup>a</sup>  | Instrument exception.                            | SYNCHRON CX(U)           |
|                   |        | 7               | Chemistry configuration change.                  | SYNCHRON CX(U)           |
|                   |        | 13              | Normal/critical range change.                    | SYNCHRON CX(U)           |
|                   |        | 17 <sup>a</sup> | End of run.                                      | SYNCHRON CX(U)           |
| Setup Status      | 704    | 7               | Request for installed<br>chemistries.            | Host                     |
|                   |        | 8               | Installed chemistries.                           | SYNCHRON CX(S)           |

1 of 2

Table 4.1 SYNCHRON CX Streams and Functions, continued

| Operation | Stream | Function | Information                         | Sent By        |
|-----------|--------|----------|-------------------------------------|----------------|
|           |        | 9        | Request for normal/critical ranges. | Host           |
|           |        | 10       | Normal/critical ranges.             | SYNCHRON CX(S) |

2 of 2

\* Supported by unidirectional interface.  
(S) indicates solicited message.  
(U) indicates unsolicited message.

Paragraphs 4.3 through 4.7 describe the transmission format details of each stream and its associated functions.

### 4.3 Stream 700 - Special Function

#### Stream 700 - Function 1 Are you there?

SENT BY: Host

PURPOSE: Identify that a bidirectional system is attached to the host port; also clears the status power-up flag.

Table 4.2 Stream 700 - Function 1 Are you there?

| Field                          | Length | Format  | Description       |
|--------------------------------|--------|---------|-------------------|
| Device I.D.                    | 2      | 0 -- 99 | Operator assigned |
| Stream                         | 3      | 700     | Always 700        |
| Function                       | 2      | 01      | Always 01         |
| EXAMPLE: [00,700,01]98<CR><LF> |        |         |                   |



## Stream 700 - Function 2 Host Setup

SENT BY: Host

PURPOSE: Identify that a bidirectional system is attached to the host port; also clears the status power-up flag.

Table 4.3 Stream 700 - Function 2 Host Setup

| Field                               | Length | Format       | Description                      |
|-------------------------------------|--------|--------------|----------------------------------|
| Device ID                           | 2      | 0 - 99       | Operator assigned                |
| Stream                              | 3      | 700          | Always 700                       |
| Function                            | 2      | 02           | Always 02                        |
| Date                                | 6      | ddmmyy       | Always ddmmyy                    |
| Time                                | 6      | hhmmss       | Always hhmmss                    |
| Software Revision                   | 10     | Alphanumeric | REV NNN.N                        |
| Driver Mode                         | 1      | B or U       | Bidirectional or Unidirectional  |
| STREAM 700 FUNCTION 2               | 1      | 0 or 1       | 0 = off, 1 = on                  |
| STREAM 701 FUNCTION 6               | 1      | 0 or 1       | 0 = off, 1 = on                  |
| STREAM 702 FUNCTION 3 <sup>1</sup>  | 1      | 0 or 1       | 0 = off, 1 = on                  |
| STREAM 702 FUNCTION 7 <sup>4</sup>  | 1      | 0 or 1       | 0 = off, 1 = on                  |
| STREAM 702 FUNCTION 9 <sup>4</sup>  | 1      | 0 or 1       | 0 = off, 1 = on                  |
| STREAM 702 FUNCTION 11              | 1      | 0 or 1       | 0 = off, 1 = on                  |
| STREAM 702 FUNCTION 23              | 1      | 0 or 1       | 0 = off, 1 = on                  |
| STREAM 702 FUNCTION 25 <sup>5</sup> | 1      | 0 or 1       | 0 = off, 1 = on                  |
| STREAM 702 FUNCTION 81 <sup>1</sup> | 1      | 0 or 1       | 0 = off, 1 = on                  |
| STREAM 702 FUNCTION 83              | 1      | 0 or 1       | 0 = off, 1 = on                  |
| STREAM 703 FUNCTION 1               | 1      | 0 or 1       | 0 = off, 1 = on                  |
| STREAM 703 FUNCTION 2               | 1      | 0 or 1       | 0 = off, 1 = on                  |
| STREAM 703 FUNCTION 5               | 1      | 0 or 1       | 0 = off, 1 = on                  |
| STREAM 703 FUNCTION 7               | 1      | 0 or 1       | 0 = off, 1 = on                  |
| STREAM 703 FUNCTION 13              | 1      | 0 or 1       | 0 = off, 1 = on                  |
| STREAM 703 FUNCTION 17              | 1      | 0 or 1       | 0 = off, 1 = on                  |
| STREAM 702 FUNCTION 13              | 1      | 0 or 1       | 0 = off, 1 = on                  |
| Unused Stream and Function          | 1      | 0 or 1       | 0 = off, 1 = on, Refer to Note 3 |
| Unused Stream and Function          | 1      | 0 or 1       | 0 = off, 1 = on, Refer to Note 3 |
| Unused Stream and Function          | 1      | 0 or 1       | 0 = off, 1 = on, Refer to Note 3 |
| Unused Stream and Function          | 1      | 0 or 1       | 0 = off, 1 = on, Refer to Note 3 |
| Unused Stream and Function          | 1      | 0 or 1       | 0 = off, 1 = on, Refer to Note 3 |
| Unused Stream and Function          | 1      | 0 or 1       | 0 = off, 1 = on, Refer to Note 3 |
| Unused Stream and Function          | 1      | 0 or 1       | 0 = off, 1 = on, Refer to Note 3 |

1 of 2

Table 4.3 Stream 700 - Function 2 Host Setup, continued

| Field                             | Length   | Format | Description |
|-----------------------------------|--|--------|-------------|
| EXAMPLE:                          | [ -0,700,02,270291,163322,v0.0.04.07,B,1,0,1,1,1,1,<br>1,1,0,1,1,1,1,1,1,1,1,1,1,1,1,1] 6E<CR><LF>   |        |             |
| <b>(- Indicates space holder)</b> |  |        |             |
| Note 1:                           | Stream 702 function 3, or Stream 702 function 81, or both Stream 702 functions 3 and 81 will be turned ON. It is not permitted to turn OFF both function 3 and 81. |        |             |
| Note 2:                           | Stream 702 function 1 and Stream 702 function 5 will always be sent.   |        |             |
| Note 3:                           | Unused stream and function fields are for future use as new streams and functions are defined.   |        |             |
| Note 4:                           | Complete reagent pack messages can be turned off/on.   |        |             |
| Note 5:                           | Should be turned off/on separately.  |        |             |

2 of 2

### Stream 700 - Function 7 Clear Queue

SENT BY: Host

**PURPOSE:** Controls the transmission of unsolicited messages (refer to Section 3.4) from the SYNCHRON CX. In response to this message the SYNCHRON CX clears the unsolicited and query queues.

## NOTICE

When the SYNCHRON CX transmission is enabled, the SYNCHRON CX transmission queue is cleared to ensure transmission of current information. All items in unsolicited and query queue are deleted and are not sent to the host.

Table 4.4 Stream 700 - Function 7 Clear Queue

| Field                              | Length | Format | Description       |
|------------------------------------|--------|--------|-------------------|
| Device ID                          | 2      | 0 - 99 | Operator assigned |
| Stream                             | 3      | 700    | Always 700        |
| Function                           | 2      | 07     | Always 07         |
| EXAMPLE:     [00,700,07]92<CR><LF> |        |        |                   |

#### 4.4 Stream 701 - Sample/Cup Program

##### Stream 701 - Function 1 Sample/Cup Program

SENT BY: Host

PURPOSE: Sample programming information, sent by the host, required by the SYNCHRON CX to process the sample in the specified cup.

Table 4.5 Stream 701 - Function 1 Sample/Cup Program

| Field                 | Length | Format       | Description  |
|-----------------------|--------|--------------|--|
| Device ID             | 2      | 0 - 99       | Operator assigned  |
| Stream                | 3      | 701          | Always 701   |
| Function              | 2      | 01           | Always 01  |
| Sector Number         | 2      | 0 - 60       | Numeric; 0 = no sector number in barcode mode  |
| Cup Number            | 2      | 0 - 7        | Numeric; 0 = no cup number in barcode mode   |
| Update Flag           | 1      | Flag         | 0 = Replace all programming<br>1 = Only add new chemistries if sample program exists, else program new sample; Refer to Note 4 |
| Test Type             | 2      | AA           | Refer to Appendix A  |
| Sample Type           | 2      | AA           | Refer to Appendix B  |
| Sample ID             | 11     | Alphanumeric | Refer to Note 2  |
| Control Name          | 20     | Alphanumeric |  |
| Sample Comment Code 1 | 25     | Alphanumeric | Refer to Note 1  |
| Sample Comment Code 2 | 25     | Alphanumeric |  |
| Last Name             | 18     | Alphanumeric |  |
| First Name            | 15     | Alphanumeric |  |
| Middle Initial        | 1      | Alphanumeric |  |
| Patient ID            | 12     | Alphanumeric |  |
| Doctor                | 18     | Alphanumeric |  |
| Draw Date             | 6      | ddmmyy       | Day Month Year<br>(e.g. 240795 = July 24, 1995)  |
| Draw Time             | 4      | hhmm         | Hours Minutes (military)*  |
| Location              | 20     | Alphanumeric |  |
| Age                   | 3      | 0-999        | in hours, days, weeks, months, years   |
| Age Units             | 1      | Numeric      | hours, days, weeks, months, years;<br>Refer to Appendix K  |
| Birth Date            | 6      | ddmmyy       | Day Month Year   |
| Sex                   | 1      | A            | M = Male; F = Female   |
| Patient Comments      | 25     | Alphanumeric |  |
| Timed Urine Volume    | 7      | Alphanumeric |  |

1 of 3

Table 4.5 Stream 701 - Function 1 Sample/Cup Program, continued

| Field  | Length | Format    | Description                   |
|--|--------|-----------|-------------------------------|
| Timed Urine Period   | 4      | Numeric   | Collection Period             |
| Timed Urine Creatinine   | 4      | Numeric   | Serum Creatinine (Note 6)     |
| Timed Urine Area   | 6      | Numeric   | Surface Area                  |
| Number of Tests  | 3      | 1 - 999   | Numeric                       |
| Programmed for Sample  |        |           |                               |
| Chemistry 1  | 4      | Chem Code |                               |
| ORDAC Chemistry 1  | 1      | Flag      | 1 = ORDAC<br>0 = Do not ORDAC |
| Chemistry 2  | 4      | Chem Code |                               |
| ORDAC Chemistry 2  | 1      | Flag      | 1 = ORDAC<br>0 = Do not ORDAC |
| :  |        |           |                               |
| Chemistry N<br>(Refer to Note 5)   | 4      | Chem Code | (N ≤ 72)                      |
| ORDAC Chemistry N  | 1      | Flag      | 1 = ORDAC<br>0 = Do not ORDAC |
| <p><b>EXAMPLE:</b>     [00,701,01,01,03,1,RO,SE,samp3-----,CONTROL_NAME_____,<br/> -----,SAMPLE_COMMENT_TWO_____,<br/> -----,PATIENT_F_NAME_,M,-----,<br/> -----,-----,-----,-----,<br/> 000,4,-----,M,-----,-----,<br/> ----,----,-----,001,01B-,0]2F&lt;CR&gt;&lt;LF&gt;</p> |        |           |                               |
| <b>(-Indicates space holder)</b>   |        |           |                               |

2 of 3

Table 4.5 Stream 701 - Function 1 Sample/Cup Program, continued

| Field   | Length   | Format            | Description |
|---------|--|-------------------|-------------|
| Note 1: | Control name is the unique identifier indicating which control this result is for. Only the selected fluid type for the control can be used. Any other fluid type is a syntax error. All lower case letters are converted to upper case.   |                   |             |
| Note 2: | The Sample ID must be filled in both modes. It is converted to upper case when received in lower case. Invalid characters:   |                   |             |
|         | 33 – 126   | (decimal) valid   |             |
|         | 44, 59, 36   | (decimal) invalid |             |
|         | * ? ] [ \ ^ , ;  | invalid           |             |
| Note 3: | The cup number and sector number must be filled in when the instrument is in sector and cup mode.  |                   |             |
| Note 4: | When update flag is on, chemistries that are programmed will not be added again to the chemistry list. When the update flag is sent from the host, the flag applies only to chemistries added to the sample program. This flag cannot be used for updating any other sample programming parameter. |                   |             |
| Note 5: | The number in the number of tests field must equal N. If Q C is programmed by cartridge then $N \leq 33$   |                   |             |
| Note 6: | Serum creatinine values must be received from the host in default units (mg/dL) for creatinine clearance calculations to be performed correctly, and to avoid erroneous answers.   |                   |             |

3 of 3

### Stream 701 - Function 2 Sample/Cup Return Status

SENT BY: SYNCHRON CX

PURPOSE: Notify the host whether or not the downloaded sample/cup program (Stream 701 - Function 1) was accepted. An OK message indicates the program was accepted. A BUSY, SYNTAX ERROR or INV. CHEM COMB. message indicates the program was not accepted (e.g. Programming was being done at the SYNCHRON CX console). An INVALID CHEMISTRY REQUESTED or INVALID ORDAC REQUESTED message indicates that a request for a particular chemistry was not accepted, and, consequently, remaining chemistries for that cup were not accepted and will not be run.

Table 4.6 Stream 701 - Function 2 Sample/Cup Return Status

| Field            | Length | Format    | Description   |
|------------------|--------|-----------|---|
| Device ID        | 2      | 0 - 99    | Operator assigned   |
| Stream           | 3      | 701       | Always 701  |
| Function         | 2      | 02        | Always 02   |
| Return Code      | 2      | NN        | 00 = OK<br>-1 = SYNTAX ERROR (Incorrect message syntax)<br>-2 = BUSY (The cup is running or being programmed at the SYNCHRON CX)<br>-3 = INVALID CHEMISTRY REQUESTED<br>-4 = INVALID ORDAC REQUESTED<br>-5 = INVALID CHEMISTRY COMBINATION PROGRAMMED (one or more chemistries require dilution or sample volume exceeds hardware limitations)<br>-6 = CONTROL NOT CONFIGURED<br>-7 = CALIBRATOR SECTOR ONLY<br>-8 = MODE MISMATCH (barcode/sector mode is not the same for SYNCHRON CX and host)<br>-9 = SYNCHRON CX Error (some internal logical error has occurred)<br>10 = COMPLETED SAMPLE (program update only)<br>11 = Incompatible Fluid Types. Used for program update or control sample program.<br>12 = Incompatible Test Types. Used only for program update.<br>13 = Incompatible Patient Name. Used only for program update (first, middle, last).<br>14-20 = Future Use. |
| Accession Number | 5      | 1 - 65535 | Unique ID assigned by SYNCHRON CX to each sample programmed; if any change to a programmed sample is made, a new accession number will be assigned.   |

1 of 2

Table 4.6 Stream 701 - Function 2 Sample/Cup Return Status, continued

| Field  | Length | Format | Description                           |
|--|--------|--------|---------------------------------------|
| Sector Number  | 2      | 0 - 60 | 0 = NO SECTOR NUMBER IN BAR CODE MODE |
| Cup Number   | 2      | 0 - 7  | 0 = NO CUP NUMBER IN BAR CODE MODE    |
| Sample ID  | 11     | Alpha  |                                       |
| <p>EXAMPLE:     [-0,701,02,-0,--230,-1,-3,samp3-----]5D&lt;CR&gt;&lt;LF&gt;</p> <p>(-Indicates space holder)</p> |        |        |                                       |

2 of 2

**Stream 701 - Function 3 Clear Sector/Sample IDs**

SENT BY: Host

PURPOSE: Clears previous programming from the indicated sector when SYNCHRON CX is in sector mode, or will clear a list of sample IDs when in bar code mode. If not sent, the new programming will cause the old sample program to be reassigned.

Table 4.7 Stream 701 - Function 3 Clear Sector/Sample IDs

| Field  | Length | Format | Description                    |
|--|--------|--------|--------------------------------|
| Device ID  | 2      | 0 - 99 | Operator assigned              |
| Stream   | 3      | 701    | Always 701                     |
| Function   | 2      | 03     | Always 03                      |
| Sector Number  | 2      | 0 - 60 | 0 = NO SECTOR IN BAR CODE MODE |
| Sample ID 1  | 11     | Alpha  | Sample ID for a cup            |
| Sample ID 2  | 11     | Alpha  | Sample ID for a cup            |
| Sample ID 3  | 11     | Alpha  | Sample ID for a cup            |
| Sample ID 4  | 11     | Alpha  | Sample ID for a cup            |
| Sample ID 5  | 11     | Alpha  | Sample ID for a cup            |
| Sample ID 6  | 11     | Alpha  | Sample ID for a cup            |
| Sample ID 7  | 11     | Alpha  | Sample ID for a cup            |
| <p>EXAMPLE:     [00,701,03,00,SAMP1-----,SAMP2-----,SAMP3-----,<br/> SAMP4-----,SAMP5-----,SAMP6-----,SAMP7-----]<br/> D2&lt;CR&gt;&lt;LF&gt;</p> <p>(-Indicates space holder)</p> |        |        |                                |

### Stream 701 - Function 4 Clear Sector/Sample IDs Status

SENT BY: SYNCHRON CX

PURPOSE: Notify the host whether or not the clear sector/sample ID command (Stream 701 - Function 3) was accepted. An OK message indicates the sector was cleared; a BUSY or BAD MESSAGE indicates the sector or not all of the sample IDs were cleared (e.g. Programming was being done at the SYNCHRON CX console).

Table 4.8 Stream 701 - Function 4 Clear Sector/Sample IDs Status

| Field   | Length | Format | Description  |
|---|--------|--------|--|
| Device ID   | 2      | 0 - 99 | Operator assigned  |
| Stream  | 3      | 701    | Always 701   |
| Function  | 2      | 04     | Always 04  |
| Sector Number   | 2      | 0 - 60 | 0 = NO SECTOR IN BAR CODE MODE   |
| Return Code   | 2      | 0 - 99 | -0 = COMPLETE SECTOR OR ALL SAMPLE IDs CLEARED<br>-1 = BAD MESSAGE (Incorrect message syntax)<br>-2 = BUSY (The sector or one of the sample IDs is running)<br>-3 = SYNCHRON CX ERROR<br>-4 = NON EXISTENT ERROR |
| Return Codes:   |        |        |  |
| Sample ID 1   | 2      | 0 - 99 | -0 = CLEARED   |
| Sample ID 2   | 2      | 0 - 99 | -2 = BUSY - NOT CLEARED  |
| Sample ID 3   | 2      | 0 - 99 | -3 = SYNCHRON CX ERROR   |
| Sample ID 4   | 2      | 0 - 99 | -4 = NON EXISTENT ERROR  |
| Sample ID 5   | 2      | 0 - 99 |  |
| Sample ID 6   | 2      | 0 - 99 |  |
| Sample ID 7   | 2      | 0 - 99 |  |
| EXAMPLE: [-0,701,04,-0,-4,-4,-4,-4,-4,-4,-4,-4]28<CR><LF> |        |        |  |
| (-Indicates space holder)                                 |        |        |  |



## Stream 701 - Function 6 HOST QUERY Sector/Sample IDs

SENT BY: SYNCHRON CX to Host

PURPOSE: To request from the host the sample program for the sample ID's specified. This message can be turned off in the host setup screen. After receiving this message, the host has a period of time to respond with the first sample program. For the CX DELTA Systems this time period is selected by the operator in the host setup screen. Options are OFF, 2.5, 5, 7.5, and 10 minutes. The CX DELTA System also has an AUTO CLEAR QUEUE function. This option, if ON, will clear all programming for the sample ID as soon as the sample is identified a second time, causing the instrument to requery the host. The host can respond by resending the same sample programming (rerun tests,), send new programming (including rerunning specific tests), or send nothing (no sample programming available).

Table 4.9 Stream 701 - Function 6 HOST QUERY Sector/Sample IDs

| Field   | Length | Format | Description       |
|---|--------|--------|-------------------|
| Device ID   | 2      | 0-99   | Operator assigned |
| Stream  | 3      | 701    |                   |
| Function  | 2      | 06     |                   |
| Sample ID 1   | 11     | Alpha  |                   |
| —   |        |        |                   |
| —   |        |        |                   |
| Sample ID 7   | 11     | Alpha  |                   |
| <br>EXAMPLE: [-0,701,06,samp1-----,samp2-----,samp3-----,<br>samp4-----,samp5-----,samp6-----,samp7-----]<br>EB<CR><LF> |        |        |                   |
| (-Indicates space holder)   |        |        |                   |

## 4.5 Stream 702 - Results

### Stream 702 - Function 1 Cup Header

SENT BY: SYNCHRON CX

PURPOSE: The cup header (sample ID and demographic information) is transmitted prior to the transmission of the test results.

Table 4.10 Stream 702 - Function 1 Cup Header

| Field                 | Length | Format       | Description   |
|-----------------------|--------|--------------|---|
| Device ID             | 2      | 0 - 99       | Operator assigned   |
| Stream                | 3      | 702          | Always 702  |
| Function              | 2      | 01           | Always 01   |
| Date Start            | 6      | ddmmyy       | Day Month Year<br>(e.g. 240795 = July 24, 1995)                             |
| Time Start            | 6      | hhmmss       | Hour Minute Second  |
| Accession Number      | 5      | number       | 1 - 65535 Unique ID assigned by<br>SYNCHRON CX to each sample<br>programmed |
| Print Type            | 2      | AA           | RG = Regular <sup>a</sup><br>RE = Recall <sup>a</sup><br>IN = INCOMPLETE    |
| Sector Number         | 2      | 1 - 60       | Numeric   |
| Cup Number            | 2      | 1 - 7        | Numeric   |
| Test Type             | 2      | Test Code    | Refer to Appendix A   |
| Future Use Space      | 9      | Alphanumeric |   |
| Sample Type           | 2      | Fluid Code   | Refer to Appendix B   |
| Sample ID             | 11     |              | Alphanumeric  |
| Control Name          | 20     |              | Alphanumeric (Note 1)   |
| Sample Comment Code 1 | 25     |              | Alphanumeric  |
| Sample Comment Code 2 | 25     |              | Alphanumeric  |
| Last Name             | 18     |              | Alphanumeric  |
| First Name            | 15     |              | Alphanumeric  |
| Middle Initial        | 1      |              | Alphanumeric  |
| Patient ID            | 12     |              | Alphanumeric  |
| Doctor                | 18     |              | Alphanumeric  |
| Draw Date             | 6      | ddmmyy       | Day Month Year  |
| Draw Time             | 4      | hhmm         | Alphanumeric  |
| Location              | 20     |              | Alphanumeric  |
| Age                   | 3      | 0 - 999      | Hours Days Weeks Months Years   |
| Age Units             | 1      |              | Number; Refer to Appendix K   |
| Birth Date            | 6      | ddmmyy       | Day Month Year  |

1 of 2

Table 4.10 Stream 702 - Function 1 Cup Header, continued

| Field   | Length | Format    | Description                           |
|---|--------|-----------|---------------------------------------|
| Sex   | 1      | A         | M = Male; F = Female                  |
| Patient Comments  | 25     |           | Alphanumeric                          |
| Timed Urine Volume  | 7      | Real      | 0-99999.0 (Note 2)                    |
| Timed Urine Period  | 4      | Real      | Collection Period                     |
| Timed Urine Creatinine  | 4      | Real      | Serum Creatinine 0-99.9 (Note 2)      |
| Timed Urine Area  | 6      | Real      | Surface Area 0-99.999 (Note 2)        |
| Number of Tests<br>Programmed for Sample  | 3      | 1 - 999   | Numeric (does not include replicates) |
| Chemistry 1   | 4      | Chem Code | Refer to Appendix C                   |
| Chemistry 2   | 4      | Chem Code | Refer to Appendix C                   |
| Chemistry N   | 4      | Chem Code | Refer to Appendix C (N= 0 - 999)      |
| <p><b>EXAMPLE:</b>    [-0,702,01,270291,113120,--212,RE,-1,-1,RO,#####,<br/> SE,28903T-----,-----,<br/> -----,<br/> Wilson-----,Joe-----,-,23895-----,<br/> Kildare-----,270291,1039,220-----,<br/> -34,5,-----,M,-----,#####,####,<br/> ####,#####,--5,44A-,04A-,02A-,01B-,01A-]53&lt;CR&gt;&lt;LF&gt;</p> <p><b>(-Indicates space holder)</b></p>                       |        |           |                                       |
| <p>Note 1: Control name is the unique identifier indicating which control this result is for.</p> <p>Note 2: The number of decimal places is adjusted so that the number can fit in this field. The urine volume field length in Sample Programming is 5 characters long. When results are transmitted to the host, all 7 characters for timed urine are transmitted.</p> |        |           |                                       |

2 of 2

<sup>a</sup> All automatic transmissions are identified as "RG"; transmission initiated through the SYNCHRON CX recall function are identified as "RE".

**Stream 702 - Function 3 Test Results**

SENT BY: SYNCHRON CX

PURPOSE: Transmission of test results when cup is complete.

Table 4.11 Stream 702 - Function 3 Test Results

| Field                     | Length | Format    | Description  |
|---------------------------|--------|-----------|--|
| Device ID                 | 2      | 0 - 99    | Operator assigned  |
| Stream                    | 3      | 702       | Always 702   |
| Function                  | 2      | 03        | Always 03  |
| Date Complete             | 6      | ddmmyy    | Day Month Year<br>(e.g. 240795 = July 24, 1995)  |
| Time Complete             | 6      | hhmmss    | Hour Minute Second   |
| Accession Number          | 5      | 1 - 65535 | Unique ID assigned by SYNCHRON CX to each sample programmed  |
| Result Record Number      | 9      | 1 - 10000 | Numeric  |
| Sector Number             | 2      | 1 - 60    | Numeric  |
| Cup Number                | 2      | 1 - 7     | Numeric  |
| Sample ID                 | 11     |           | Alphanumeric   |
| Chemistry                 | 4      | Chem Code | Refer to Appendix C  |
| Reagent Serial No.        | 3      |           | Alphanumeric   |
| Reagent Lot No.           | 6      |           | Alphanumeric   |
| Cuvette No.               | 2      |           | Numeric  |
| Replicate No.             | 2      |           | Numeric  |
| Result in Selected Units* | 9      | Real      | Numeric or ***** = invalid data or<br>##### = result suppressed. Result in user-selected units which is reported on the SYNCHRON CX. |
| Calibration Rate          | 9      | Real      |  |
| Positive or Negative      | 1      | Code      | Positive or Negative<br>0 = Negative <sup>a</sup><br>1 = Positive <sup>a</sup><br>2 = Not Applicable<br>3 = No Decision              |
| Suppress Result           | 1      | Flag      | 0 = Do not suppress result<br>1 = Suppress result  |
| Units                     | 2      | Code      | Refer to Appendix D  |
| Normal Range Flag         | 2      | Code      | **Refer to Range Codes below   |
| Instrument Range Flag     | 2      | Code      | **Refer to Range Codes below   |
| Critical Range Flag       | 2      | AA        | **Refer to Range Codes below   |

1 of 3

Table 4.11 Stream 702 - Function 3 Test Results, continued

| Field              | Length | Format | Description  |
|--------------------|--------|--------|--|
| ORDAC Result       | 1      | N      | 0 = ORDAC not used<br>1 = ORDAC used   |
| Control Range Flag | 2      | AA     | NA = Not applicable<br>NR = Within 2 SD of mean<br>H2 = Between 2 SD and 3 SD above mean<br>H3 = Greater than 3 SD above mean<br>H4 = Greater than 4SD above mean<br>L2 = Between 2SD + 3SD below mean<br>L3 = Greater than 3SD below mean<br>L4 = Greater than 4SD below mean<br>IT = Invalid SD<br>IC = Incomplete Ranges<br>CD = Control or chemistry deleted |
| Calculated Result  | 9      | Real   | Numeric or ***** = invalid data or<br>##### = result suppressed. Result<br>calculated using the SYNCHRON CX<br>default units.  |
| Instrument Codes   | 9      |        | Refer to Appendix E  |
| Result Error 1     | 2      | Code   | Refer to Appendix F  |
| Result Error 2     | 2      | Code   | Refer to Appendix F  |
| Result Error 3     | 2      | Code   | Refer to Appendix F  |
| Result Error 4     | 2      | Code   | Refer to Appendix F  |
| Result Error 5     | 2      | Code   | Refer to Appendix F  |
| Result Error 6     | 2      | Code   | Refer to Appendix F  |
| Result Error 7     | 2      | Code   | Refer to Appendix F  |
| Result Error 8     | 2      | Code   | Refer to Appendix F  |
| Result Error 9     | 2      | Code   | Refer to Appendix F  |
| Result Error 10    | 2      | Code   | Refer to Appendix F  |
| Result Error 11    | 2      | Code   | Refer to Appendix F  |
| Result Error 12    | 2      | Code   | Refer to Appendix F  |
| Result Error 13    | 2      | Code   | Refer to Appendix F  |
| Result Error 14    | 2      | Code   | Refer to Appendix F  |
| Result Error 15    | 2      | Code   | Refer to Appendix F  |
| Result Error 16    | 2      | Code   | Refer to Appendix F  |
| Dilution Factor    | 6      | Real   | 0.0 - 1000.0 (refer to Note 1)   |
| Reportable Range   | 2      | Code   | **Refer to range codes below   |
| Future Use Space   | 20     | Alpha  |  |

2 of 3

Table 4.11 Stream 702 - Function 3 Test Results, continued

| Field    | Length   | Format | Description |
|----------|--|--------|-------------|
| EXAMPLE: | [-0,702,03,270291,114739,--212,-----3468,-1,-1,28903T-----,44A-,2D6,008150,12,-1,-----150,#####,2,0,-0,NA,NR,NR,0,NA,149.51265,-----,NO, |        |             |

3 of 3

<sup>a</sup> 0 and 1 are used only for drugs of abuse. All other chemistries use 2, not applicable.

### Stream 702 - Function 5 End of Cup

SENT BY: SYNCHRON CX

PURPOSE: Notify host that all tests on a sample have been completed.

Table 4.12 Stream 702 - Function 5 End of Cup

| Field   | Length | Format    | Description   |
|---|--------|-----------|---|
| Device ID   | 2      | 0 - 99    | Operator assigned   |
| Stream  | 3      | 702       | Always 702  |
| Function  | 2      | 05        | Always 05   |
| Date Completed  | 6      | ddmmyy    | Day Month Year<br>(e.g. 240795 = July 24, 1995)                   |
| Time Completed  | 6      | hhmmss    | Hour Minute Second  |
| Accession Number  | 5      | 1 - 65535 | Unique<br>ID assigned by SYNCHRON CX to<br>each sample programmed |
| Sample ID   | 11     |           | Alphanumeric  |
| Sector Number   | 2      | 1 - 60    | Numeric   |
| Cup Number  | 2      | 1 - 7     | Numeric   |
| EXAMPLE: [-0,702,05,220291,155308,--105,1-----,58,-1]<br>2F<CR><LF> |        |           |   |
| (-Indicates space holder)   |        |           |   |

### Stream 702 - Function 7 Linear Calibration Result

SENT BY: SYNCHRON CX

PURPOSE: Transmission of linear calibration results. This message is only transmitted after all tests for the calibration are complete. This message can be turned off in the Host Setup screen.

Table 4.13 Stream 702 - Function 7 Linear Calibration Result

| Field                  | Length | Format    | Description                                     |
|------------------------|--------|-----------|---|
| Device ID              | 2      | 0 - 99    | Operator assigned                               |
| Stream                 | 3      | 702       | Always 702                                      |
| Function               | 2      | 07        | Always 07                                       |
| Date Completed         | 6      | ddmmyy    | Day Month Year<br>(e.g. 240795 = July 24, 1995) |
| Time Completed         | 6      | hhmmss    | Hour Minute Second                              |
| Chemistry              | 4      | Chem Code | Refer to Appendix C                             |
| Reagent Serial No.     | 3      |           | Alphanumeric                                    |
| Reagent Lot No.        | 6      |           | Alphanumeric                                    |
| Reag Record Number     | 9      |           | Numeric   |
| Set Point 1            | 9      | Real      | Numeric - Calibrator 1 or ne<br>(Note 1)        |
| Set Point 2            | 9      |           | Numeric - Calibrator 2 or lo<br>(Note 1)        |
| Set Point 3            | 9      | Real      | Calibrator 3 or mi (Note 1)                     |
| Cal Rate               | 9      | Real      | Calibrator ne (Note 2)                          |
| Cal Rate               | 9      | Real      | Calibrator lo (Note 2)                          |
| Cal Rate               | 9      | Real      | Calibrator mi (Note 2)                          |
| Calibration            | 1      | Flag      | 0 = Not bypassed<br>1 = Bypassed                |
| Calibration Override   | 1      | Flag      | 0 = Not overridden<br>1 = Overridden            |
| Calibration Timeout    | 1      | Flag      | 0 = Not timed out<br>1 = Timed out              |
| Calibration Calibrated | 1      | Flag      | 0 = Not calibrated<br>1 = Calibrated            |
| Set Points Modified    | 1      | Flag      | 0 = Not modified<br>1 = Modified                |
| Calibration Error 1    | 2      | Code      | Refer to Appendix G                             |
| Calibration Error 2    | 2      | Code      | Refer to Appendix G                             |
| Calibration Error 3    | 2      | Code      | Refer to Appendix G                             |
| Calibration Error 4    | 2      | Code      | Refer to Appendix G                             |
| Calibration Error 5    | 2      | Code      | Refer to Appendix G                             |

1 of 2



Table 4.13 Stream 702 - Function 7 Linear Calibration Result, continued

| Field  | Length | Format     | Description                  |
|--|--------|------------|------------------------------|
| Calibration Error 6  | 2      | Code       | Refer to Appendix G          |
| Calibration Error 7  | 2      | Code       | Refer to Appendix G          |
| Calibration Error 8  | 2      | Code       | Refer to Appendix G          |
| Calibration Error 9  | 2      | Code       | Refer to Appendix G          |
| Calibration Error 10   | 2      | Code       | Refer to Appendix G          |
| Calibration Error 11   | 2      | Code       | Refer to Appendix G          |
| Calibration Error 12   | 2      | Code       | Refer to Appendix G          |
| Calibration Error 13   | 2      | Code       | Refer to Appendix G          |
| Calibration Error 14   | 2      | Code       | Refer to Appendix G          |
| Calibration Error 15   | 2      | Code       | Refer to Appendix G          |
| Calibration Error 16   | 2      | Code       | Refer to Appendix G          |
| Calibration Slope  | 10     | NNNN.NNNNN | Real; ***** = Does not fit   |
| Calibration Offset   | 10     | NNNN.NNNNN | Real; ***** = Does not fit   |
| Slope Adjustment   | 10     | NNNN.NNNNN | Real; ***** = Does not fit   |
| Offset Adjustment  | 10     | NNNN.NNNNN | Real; ***** = Does not fit   |
| Cal Factor   | 9      |            | Real; (Note 2)               |
| Cal Span   | 9      |            | Real; calibrator ne (Note 2) |
| Cal Span   | 9      |            | Real; calibrator lo (Note 2) |
| Cal Span   | 9      |            | Real; calibrator mi (Note 2) |
| Future Use Space   | 30     |            | alpha                        |
| <p><b>EXAMPLE:</b> [-0,702,07,280291,120101,83A-,0RO,010226,-----34,<br/> 48.500000,#####,#####,#####,#####,<br/> #####,0,0,0,1,0,00,00,00,00,00,00,00,00,00,00,00,<br/> 00,00,00,00,-281.10751,---0.00000,---1.00000,---0.00000,<br/> #####,#####,#####,#####,<br/> #####] 2F&lt;CR&gt;&lt;LF&gt;</p> <p><b>(-Indicates space holder)</b></p> <p>Note 1: The number of decimal places in this field is the number in the selected precision for the default units field in the chem database plus an additional decimal place</p> <p>Note 2: 5 decimal places</p> |        |            |                              |

2 of 2

### Stream 702 - Function 9 Multipoint Calibration Result

SENT BY: SYNCHRON CX

PURPOSE: Transmission of multipoint calibration results. This message is only transmitted after all tests for the calibration are complete. This message can be turned off on the Host Setup screen.

Table 4.14 Stream 702 - Function 9 Multipoint Calibration Result

| Field                  | Length | Format    | Description                                     |
|------------------------|--------|-----------|---|
| Device ID              | 2      | 0 - 99    | Operator assigned                               |
| Stream                 | 3      | 702       | Always 702                                      |
| Function               | 2      | 09        | Always 09                                       |
| Date Completed         | 6      | ddmmyy    | Day Month Year<br>(e.g. 240795 = July 24, 1995) |
| Time Completed         | 6      | hhmmss    | Hour Minute Second                              |
| Chemistry              | 4      | Chem Code | Refer to Appendix C                             |
| Reagent Serial No.     | 3      |           | Alphanumeric                                    |
| Reagent Lot No.        | 6      |           | Alphanumeric                                    |
| Reagent Record Number  | 9      |           | Numeric   |
| Set Point 1            | 9      |           | Real - Calibrator 1 (Note 1)                    |
| Set Point 2            | 9      |           | Real - Calibrator 2 (Note 1)                    |
| Set Point 3            | 9      |           | Real - Calibrator 3 (Note 1)                    |
| Set Point 4            | 9      |           | Real - Calibrator 4 (Note 1)                    |
| Set Point 5            | 9      |           | Real - Calibrator 5 (Note 1)                    |
| Set Point 6            | 9      |           | Real - Calibrator 6 (Note 1)                    |
| Calibration Bypassed   | 1      | Flag      | 0 = Not bypassed<br>1 = Bypassed                |
| Calibration Override   | 1      | Flag      | 0 = Not overridden<br>1 = Overridden            |
| Calibration Timeout    | 1      | Flag      | 0 = Not timed out<br>1 = Timed out              |
| Calibration Calibrated | 1      | Flag      | 0 = Not calibrated<br>1 = Calibrated            |
| Set Points Modified    | 1      | Flag      | 0 = Not modified<br>1 = Modified                |
| Calibration Error 1    | 2      | Code      | Refer to Appendix G                             |
| Calibration Error 2    | 2      | Code      | Refer to Appendix G                             |
| Calibration Error 3    | 2      | Code      | Refer to Appendix G                             |
| Calibration Error 4    | 2      | Code      | Refer to Appendix G                             |
| Calibration Error 5    | 2      | Code      | Refer to Appendix G                             |
| Calibration Error 6    | 2      | Code      | Refer to Appendix G                             |
| Calibration Error 7    | 2      | Code      | Refer to Appendix G                             |

1 of 2

Table 4.14 Stream 702 - Function 9 Multipoint Calibration Result, continued

| Field  | Length | Format    | Description         |
|--|--------|-----------|---------------------|
| Calibration Error 8  | 2      | Code      | Refer to Appendix G |
| Calibration Error 9  | 2      | Code      | Refer to Appendix G |
| Calibration Error 10   | 2      | Code      | Refer to Appendix G |
| Calibration Error 11   | 2      | Code      | Refer to Appendix G |
| Calibration Error 12   | 2      | Code      | Refer to Appendix G |
| Calibration Error 13   | 2      | Code      | Refer to Appendix G |
| Calibration Error 14   | 2      | Code      | Refer to Appendix G |
| Calibration Error 15   | 2      | Code      | Refer to Appendix G |
| Calibration Error 16   | 2      | Code      | Refer to Appendix G |
| Model  | 2      | 1 - 5,8,9 | Model number        |
| RO   | 9      |           | Real (Note 2)       |
| Kc   | 9      |           | Real (Note 2)       |
| a  | 9      |           | Real (Note 2)       |
| b  | 9      |           | Real (Note 2)       |
| c  | 9      |           | Real (Note 2)       |
| d  | 9      |           | Real (Note 2)       |
| Slope Adjustment   | 10     |           | Real (Note 2)       |
| Offset Adjustment  | 10     |           | Real (Note 2)       |
| Future Use Space   | 30     |           | Alpha               |
| <p><b>EXAMPLE:</b> [-0,702,09,060391,103806,MPT6,###,#####,------10,<br/> 0.0000000,1.0000000,2.0000000,4.0000000,8.0000000,<br/> 16.000000,0,0,0,1,0,00,00,00,00,00,00,00,00,00,00,<br/> 00,00,00,00,00,00,-1,879.68268,-787.9960,0.1827289,<br/> 0.7196383,0.0000000,0.0000000,---1.00000,---0.00000,<br/> #####]90&lt;CR&gt;&lt;LF&gt;</p> <p><b>(-Indicates space holder)</b></p> <p>Note 1: The number of decimal places in this field is the number in the selected precision for the default units field in the chem database plus an additional decimal place</p> <p>Note 2: The number of decimal places is adjusted so that the number can fit in this field</p> |        |           |                     |

2 of 2

### Stream 702 - Function 11 Special Calculation Result

SENT BY: SYNCHRON CX

PURPOSE: Transmission of special calculation results. One message is transmitted for each special calculation result. This message cannot be transmitted until all test results are completed. This message can be turned off on the Host Setup screen.

Table 4.15 Stream 702 - Function 11 Special Calculation Result

| Field   | Length | Format    | Description   |
|---|--------|-----------|---|
| Device ID   | 2      | 0 - 99    | Operator assigned   |
| Stream  | 3      | 702       | Always 702  |
| Function  | 2      | 11        | Always 11   |
| Date Run  | 6      | ddmmyy    | Day Month Year<br>(e.g. 240795 = July 24, 1995)             |
| Time Run  | 6      | hhmmss    | Hour Minute Second  |
| Accession Number  | 5      | 1 - 65535 | Unique ID assigned by SYNCHRON CX to each sample programmed |
| Sector Number   | 2      | 1 - 60    | Numeric   |
| Cup Number  | 2      | 1 - 7     | Numeric   |
| Sample ID   | 11     |           | Alphanumeric  |
| Replicate No.   | 2      |           | Numeric   |
| Special Calculation Name  | 20     |           | Alphanumeric; Name given by operator                        |
| Special Calculation Status  | 2      | Code      | Refer to Appendix H   |
| Special Calculation Result  | 9      |           | Real (Note 1)   |
| Special Unit String   | 8      |           | Alpha   |
| <p>EXAMPLE: [-0,702,11,060391,123246,---19,18,-3,-----,-1,OSMOLALITY-(1)-----,OK,308.60977,-----]46&lt;CR&gt;&lt;LF&gt;</p> <p>(- Indicates space holder)</p> |        |           |   |
| Note 1: The number of decimal places is adjusted so that the number can fit in this field.  |        |           |   |

### Stream 702 - Function 13 Timed Urine Result

SENT BY: SYNCHRON CX

PURPOSE: This message is sent whenever all results are completed. There is one message for each timed urine result. This message can be turned off on the Host Setup screen.

Table 4.16 Stream 702 - Function 13 Timed Urine Result

| Field  | Length | Format | Description         |
|--|--------|--------|---------------------|
| Device ID  | 2      | 0 - 99 | Operator assigned   |
| Stream   | 3      | 702    |                     |
| Function   | 2      | 13     |                     |
| Date   | 6      | ddmmyy |                     |
| Time   | 6      | hhmmss |                     |
| Accession Number   | 5      | Number | 1 to 65535          |
| Sector Number  | 2      | 1 - 60 |                     |
| Cup Number   | 2      | 1 - 7  |                     |
| Sample ID  | 11     | Alpha  |                     |
| Replicate #  | 2      | Number |                     |
| TU Calc Name   | 20     | Alpha  |                     |
| TU Calc Status   | 2      | Code   | Refer to Appendix H |
| TU Calc Results  | 9      | Real   | Refer to Note 1     |
| TU Unit String   | 8      | Alpha  |                     |
| <p>EXAMPLE: [-0,702,13,270291,114148,--215,-1,-3,866463K----,-1,CL-----,OK,173.48641,mmol/24.]CD&lt;CR&gt;&lt;LF&gt;</p> <p>(- Indicates space holder)</p> |        |        |                     |
| <p>Note 1: The number of decimal places is adjusted so that the number can fit in this field.</p>  |        |        |                     |

**Stream 702 - Function 21 Reagent Pack Header**

SENT BY: SYNCHRON CX

PURPOSE: A reagent pack header is sent before sending any calibration results. The header contains reagent and chemistry information.

Table 4.17 Stream 702 - Function 21 Reagent Pack Header

| Field                | Length | Format    | Description                                   |
|----------------------|--------|-----------|---|
| Device ID            | 2      | 0 - 99    | Operator assigned                             |
| Stream               | 3      | 702       |   |
| Function             | 2      | 21        |   |
| Start Date           | 6      | ddmmyy    | Day Month Year                                |
| Start time           | 6      | hhmmss    | Hours Minutes Seconds                         |
| COM CHEM ID          | 4      | Chem Code | Refer to Appendix C or User Defined Chemistry |
| Print Type           | 2      | Code      | RG=Regular<br>RE=Recall<br>IN=Incomplete      |
| Serial reagent       | 3      | Alpha     |   |
| Lot reagent          | 6      | Alpha     |   |
| Calibrator Lot       | 7      | Alpha     |   |
| Number of Level      | 2      | Number    | 1-6   |
| Number of reps cup 1 | 2      | Number    | 1-5   |
| Number of reps cup 2 | 2      | Number    | 1-5   |
| Number of reps cup 3 | 2      | Number    | 1-5   |
| Number of reps cup 4 | 2      | Number    | 1-5   |
| Number of reps cup 5 | 2      | Number    | 1-5   |
| Number of reps cup 6 | 2      | Number    | 1-5   |
| Reagent position     | 2      | Number    | 25-33 = CX3/ISE Module<br>1-24 CX4 positions  |
| Sector number        | 2      | 1-60      |   |
| Cup number           | 2      | 1-7       | first cup                                     |
| Cup number           | 2      | 1-7       |   |
| Cup number           | 2      | 1-7       |   |
| Cup number           | 2      | 1-7       |   |
| Cup number           | 2      | 1-7       |   |
| Cup number           | 2      | 1-7       | sixth cup                                     |
| Future Use Space     | 30     | Alpha     |   |

1 of 2

Table 4.17 Stream 702 - Function 21 Reagent Pack Header, continued

| Field  | Length | Format | Description |
|--|--------|--------|-------------|
| <p>EXAMPLE:     [-0,702,21,270291,102849,01A-,RG,---,N/A---,N/A----,<br/>                        -2,-2,-2,00,00,00,00,25,-0,-0,-0,-0,-0,-0,-0,<br/>                        #####] FF&lt;CR&gt;&lt;LF&gt;</p> <p>(- Indicates space holder)</p> |        |        |             |

2 of 2

**Stream 702 - Function 23 Calibration Results**

SENT BY:     SYNCHRON CX

PURPOSE:    This message is sent when a chemistry is calibrated

Table 4.18 Stream 702 - Function 23 Calibration Results

| Field                | Length | Format    | Description                                   |
|----------------------|--------|-----------|---|
| Device ID            | 2      | 0 - 99    | Operator assigned                             |
| Stream               | 3      | 702       |   |
| Function             | 2      | 23        |   |
| Completion date      | 6      | ddmmyy    | Days Months Years                             |
| Completion time      | 6      | hhmmss    | Hours Minutes Seconds                         |
| Result record number | 9      | Number    | 1 - 10,000                                    |
| CHEM ID              | 4      | CHEM CODE | Refer to Appendix C or User Defined Chemistry |
| Cuvette              | 2      | Number    |   |
| Replicate            | 2      | Number    |   |
| Cal Level            | 2      | 1 - 6     |   |
| Units                | 2      | Code      | Refer to Appendix D                           |
| Calculated Result    | 9      | Real      | Refer to Note 1                               |
| Suppress Value       | 1      | Flag      | 1 = suppress value and Refer to Note 2        |
| Instrument Codes     | 9      | Alpha     | Refer to Appendix E                           |
| Result Error 1       | 2      | Alpha     | Refer to Appendix F                           |
| Result Error 2       | 2      | Alpha     | Refer to Appendix F                           |
| Result Error 3       | 2      | Alpha     | Refer to Appendix F                           |
| Result Error 4       | 2      | Alpha     | Refer to Appendix F                           |
| Result Error 5       | 2      | Alpha     | Refer to Appendix F                           |
| Result Error 6       | 2      | Alpha     | Refer to Appendix F                           |
| Result Error 7       | 2      | Alpha     | Refer to Appendix F                           |
| Result Error 8       | 2      | Alpha     | Refer to Appendix F                           |

1 of 2

Table 4.18 Stream 702 - Function 23 Calibration Results, continued

| Field  | Length | Format | Description         |
|--|--------|--------|---------------------|
| Result Error 9   | 2      | Alpha  | Refer to Appendix F |
| Result Error 10  | 2      | Alpha  | Refer to Appendix F |
| Result Error 11  | 2      | Alpha  | Refer to Appendix F |
| Result Error 12  | 2      | Alpha  | Refer to Appendix F |
| Result Error 13  | 2      | Alpha  | Refer to Appendix F |
| Result Error 14  | 2      | Alpha  | Refer to Appendix F |
| Result Error 15  | 2      | Alpha  | Refer to Appendix F |
| Result Error 16  | 2      | Alpha  | Refer to Appendix F |
| Future Use Space   | 30     | Alpha  |                     |
| <p><b>EXAMPLE:</b>     [-0,702,23,200291,124727,-----44,83A-,76,-1,-1,-0,<br/> -999.0000,1,-----,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,NO,<br/> NO,NO,NO,NO,NO,#####] F9&lt;CR&gt;&lt;LF&gt;</p> <p><b>(- Indicates space holder)</b></p> <p>Note 1: The number of decimal places is adjusted so that the number can fit in this field.</p> <p>Note 2: If answer suppressed = 1</p> |        |        |                     |

2 of 2



## Stream 702 - Function 25 Expanded Result for Calibration

SENT BY: SYNCHRON CX

PURPOSE: This message is sent when a calibration is complete and contains an image of the data base result record. It is used primarily for diagnostic and research purposes. This message can be turned off on the Host Setup screen.

Table 4.19 Stream 702 - Function 25 Expanded Result for Calibration

| Field                 | Length | Format    | Description                                   |
|-----------------------|--------|-----------|---|
| Device ID             | 2      | 0 - 99    | Operator assigned                             |
| Stream                | 3      | 702       |   |
| Function              | 2      | 25        |   |
| Date                  | 6      | ddmmyy    | Day Month Year                                |
| Time                  | 6      | hhmmss    | Hours Minutes Seconds                         |
| Result record number  | 9      | Number    | 1 - 10,000                                    |
| Reagent record number | 9      | Number    |   |
| Replicate number      | 2      | Number    |   |
| COM CHEM ID           | 4      | CHEM CODE | Refer to Appendix C or User Defined Chemistry |
| Cuvette number        | 2      | Number    |   |
| Cal level             | 2      | 1-6       |   |
| Noise Maximum         | 5      | Number    |   |
| Outlier Maximum       | 5      | Number    |   |
| Noise Thresh Total    | 3      | Number    |   |
| Outlier Thresh Total  | 3      | Number    |   |
| Initial Abs           | 9      | Real      | Refer to Note 1                               |
| Final Abs             | 9      | Real      | Refer to Note 1                               |
| Water Blank- Rate     | 9      | Real      | Refer to Note 1                               |
| - Abs                 | 9      | Real      | Refer to Note 1                               |
| - Mean Dev            | 9      | Real      | Refer to Note 1                               |
| - Max Dev             | 9      | Real      | Refer to Note 1                               |
| Blank - Rate          | 9      | Real      | Refer to Note 1                               |
| - Abs                 | 9      | Real      | Refer to Note 1                               |
| - Mean Dev            | 9      | Real      | Refer to Note 1                               |
| - Max Dev             | 9      | Real      | Refer to Note 1                               |
| Reaction - Rate       | 9      | Real      | Refer to Note 1                               |
| - Abs                 | 9      | Real      | Refer to Note 1                               |
| - Mean Dev            | 9      | Real      | Refer to Note 1                               |
| - Max Dev             | 9      | Real      | Refer to Note 1                               |
| Raw Result            | 9      | Real      | Refer to Note 1                               |
| Calculated Result     | 9      | Real      | Refer to Note 1                               |
| Units                 | 2      | Code      | Refer to Appendix D                           |
| Instrument Codes      | 9      | Alpha     | Refer to Appendix E                           |

1 of 2



### Stream 702 - Function 27 End of Reagent Pack

SENT BY: SYNCHRON CX

PURPOSE: A reagent pack end is sent after all the calibration results for a pack, indicating information for that reagent pack is completed.

Table 4.20 Stream 702 - Function 27 End of Reagent Pack

| Field  | Length | Format    | Description                                      |
|--|--------|-----------|--|
| Device ID  | 2      | 0 - 99    | Operator assigned                                |
| Stream   | 3      | 702       | Always 702                                       |
| Function   | 2      | 27        | Always 27  |
| Start Date   | 6      | ddmmyy    | Day Month Year<br>(e.g. 240795 = July 24, 1995)  |
| Start Time   | 6      | hhmmss    | Hour Minute Second                               |
| COM CHEM ID  | 4      | CHEM CODE | Refer to Appendix C or User<br>Defined Chemistry |
| Reagent Serial No.   | 3      | Alpha     |  |
| Reagent Lot No.  | 6      | Alpha     |  |
| Calibrator Lot No.   | 7      | Alpha     |  |
| Reagent Position   | 2      | Numeric   | 0=CX3/CX3 DELTAISE Module,<br>1-24 for CX4 chems |
| Sector Number  | 2      | 1 - 60    |  |
| Cup Number   | 2      | 1 - 7     | first cup  |
| Cup Number   | 2      | 1 - 7     |  |
| Cup Number   | 2      | 1 - 7     |  |
| Cup Number   | 2      | 1 - 7     |  |
| Cup Number   | 2      | 1 - 7     |  |
| Cup Number   | 2      | 1 - 7     | sixth cup  |
| EXAMPLE: [-0,702,27,200291,124727,83A-,0RO,010226,M911051,<br>-6,17,-2,-0,-0,-0,-0,-0]EB<CR><LF> |        |           |  |
| (-Indicates space holder)  |        |           |  |

**Stream 702 - Function 81**

SENT BY: SYNCHRON CX

PURPOSE: Transmission of expanded test result message. An expanded test result message is only transmitted when the expanded result option is enabled at the SYNCHRON CX. This message can be turned off on the Host Setup screen.

Table 4.21 Stream 702 - Function 81

| Field                 | Length | Format    | Description   |
|-----------------------|--------|-----------|---|
| Device ID             | 2      | 0 - 99    | Operator assigned   |
| Stream                | 3      | 702       | Always 702  |
| Function              | 2      | 81        | Always 81   |
| Date Run              | 6      | ddmmyy    | Day Month Year<br>(e.g. 240795 = July 24, 1995)                   |
| Time Run              | 6      | hhmmss    | Hour Minute Second  |
| Accession Number      | 5      | 1 - 65535 | Unique ID assigned by<br>SYNCHRON CX to each sample<br>programmed |
| Sample Record Number  | 9      |           | Numeric   |
| Result Record Number  | 9      |           | Numeric   |
| Sample ID             | 11     |           | Alphanumeric  |
| Sector Number         | 2      | 1 - 60    | Numeric   |
| Cup Number            | 2      | 1 - 7     | Numeric   |
| Reagent Record Number | 9      |           | Numeric   |
| Replicate No.         | 2      |           | Numeric   |
| COM CHEM ID           | 4      | Chem Code | Refer to Appendix C   |
| Cuvette Number        | 2      |           | Numeric   |
| Calibrator Level      | 1      | 0 - 6     | Not used  |
| Suppressed Result     | 1      | Flag      | 0 = Result not suppressed<br>1 = Suppressed result                |
| ORDAC Result          | 1      | Flag      | 0 = ORDAC not used<br>1 = ORDAC used                              |
| Noise Maximum         | 5      |           | Numeric   |
| Outlier Maximum       | 5      |           | Numeric   |
| Noise Thresh Total    | 3      |           | Numeric   |
| Outlier Thresh Total  | 3      |           | Numeric   |
| Initial Absorbance    | 9      |           | Real Refer to Note 1  |
| Final Absorbance      | 9      |           | Real Refer to Note 1  |
| Water Blank - Rate    | 9      |           | Real Refer to Note 1  |
| - Abs                 | 9      |           | Real Refer to Note 1  |
| - Mean Dev            | 9      |           | Real Refer to Note 1  |
| - Max Dev             | 9      |           | Real Refer to Note 1  |

1 of 3

Table 4.21 Stream 702 - Function 81, continued

| Field                | Length     | Format | Description  |
|----------------------|------------|--------|--|
| Blank                | - Rate     | 9      | Real Refer to Note 1   |
|                      | - Abs      | 9      | Real Refer to Note 1   |
|                      | - Mean Dev | 9      | Real Refer to Note 1   |
|                      | - Max Dev  | 9      | Real Refer to Note 1   |
| Reaction             | - Rate     | 9      | Real Refer to Note 1   |
|                      | - Abs      | 9      | Real Refer to Note 1   |
|                      | - Mean Dev | 9      | Real Refer to Note 1   |
|                      | - Max Dev  | 9      | Real Refer to Note 1   |
| Raw Result           |            | 9      | Real Refer to Note 1   |
| Calculated Result    |            | 9      | Real Refer to Note 1   |
| Calculated Rate      |            | 9      | Real Refer to Note 1   |
| Positive or Negative |            | 1      | Flag 0 = negative<br>1 = positive<br>2 = not applicable<br>3 = no decision |
| Units                |            | 2      | Code Refer to Appendix D   |
| Instrument Codes     |            | 9      | Alpha Refer to Appendix E  |
| Result Error 1       |            | 2      | Code Refer to Appendix F   |
| Result Error 2       |            | 2      | Code Refer to Appendix F   |
| Result Error 3       |            | 2      | Code Refer to Appendix F   |
| Result Error 4       |            | 2      | Code Refer to Appendix F   |
| Result Error 5       |            | 2      | Code Refer to Appendix F   |
| Result Error 6       |            | 2      | Code Refer to Appendix F   |
| Result Error 7       |            | 2      | Code Refer to Appendix F   |
| Result Error 8       |            | 2      | Code Refer to Appendix F   |
| Result Error 9       |            | 2      | Code Refer to Appendix F   |
| Result Error 10      |            | 2      | Code Refer to Appendix F   |
| Result Error 11      |            | 2      | Code Refer to Appendix F   |
| Result Error 12      |            | 2      | Code Refer to Appendix F   |
| Result Error 13      |            | 2      | Code Refer to Appendix F   |
| Result Error 14      |            | 2      | Code Refer to Appendix F   |
| Result Error 15      |            | 2      | Code Refer to Appendix F   |
| Result Error 16      |            | 2      | Code Refer to Appendix F   |
| Intermediate ADC 1   |            | 5      | Numeric  |
| Intermediate ADC 2   |            | 5      | Numeric  |
| Intermediate ADC 3   |            | 5      | Numeric  |
| Intermediate ADC 4   |            | 5      | Numeric  |
| Intermediate ADC 5   |            | 5      | Numeric  |
| Dilution Factor      |            | 6      | Real 0.0 - 1000 Refer to Note 2  |

2 of 3

Table 4.21 Stream 702 - Function 81, continued

[illegible]

3 of 3

**Stream 702 - Function 83 Expanded Multipoint Calibration**

SENT BY: SYNCHRON CX

PURPOSE: This message is the result of a multipoint calibration and is sent whenever all tests for a multipoint calibration are completed.

Table 4.22 Stream 702 - Function 83 Expanded Multipoint Calibration

| Field                 | Length | Format       | Description         |
|-----------------------|--------|--------------|---------------------|
| Device ID             | 2      | 0 - 99       | Operator assigned   |
| Stream                | 3      | 702          | Always 702          |
| Function              | 2      | 83           | Always 83           |
| Date                  | 6      | ddmmyy       |                     |
| Time                  | 6      | hhmmss       |                     |
| Comm Chem ID          | 4      | Chem Code    |                     |
| Reagent Serial No.    | 3      | Alphanumeric |                     |
| Reagent Lot No.       | 6      | Alphanumeric |                     |
| Reagent Record Number | 9      | Numeric      |                     |
| Set Point 1           | 9      | Numeric      | Refer to Note 1     |
| Set Point 2           | 9      | Numeric      | Refer to Note 1     |
| Set Point 3           | 9      | Numeric      | Refer to Note 1     |
| Set Point 4           | 9      | Numeric      | Refer to Note 1     |
| Set Point 5           | 9      | Numeric      | Refer to Note 1     |
| Set Point 6           | 9      | Numeric      | Refer to Note 1     |
| Cal Bypassed          | 1      | Flag         | 1 = bypassed        |
| Cal Override          | 1      | Flag         | 1 = bypassed        |
| Cal Time-out          | 1      | Flag         | 1 = bypassed        |
| Cal Calibrated        | 1      | Flag         | 1 = calibrated      |
| Set Points Modified   | 1      | Flag         | 1 = modified        |
| Cal Error - 1         | 2      | Code         | Refer to Appendix G |
| Cal Error - 2         | 2      | Code         | Refer to Appendix G |
| Cal Error - 3         | 2      | Code         | Refer to Appendix G |
| Cal Error - 4         | 2      | Code         | Refer to Appendix G |
| Cal Error - 5         | 2      | Code         | Refer to Appendix G |
| Cal Error - 6         | 2      | Code         | Refer to Appendix G |
| Cal Error - 7         | 2      | Code         | Refer to Appendix G |
| Cal Error - 8         | 2      | Code         | Refer to Appendix G |
| Cal Error - 9         | 2      | Code         | Refer to Appendix G |
| Cal Error - 10        | 2      | Code         | Refer to Appendix G |
| Cal Error - 11        | 2      | Code         | Refer to Appendix G |
| Cal Error - 12        | 2      | Code         | Refer to Appendix G |

1 of 3

Table 4.22 Stream 702 - Function 83 Expanded Multipoint Calibration, continued

| Field                  | Length | Format      | Description         |
|------------------------|--------|-------------|---------------------|
| Cal Error - 13         | 2      | Code        | Refer to Appendix G |
| Cal Error - 14         | 2      | Code        | Refer to Appendix G |
| Cal Error - 15         | 2      | Code        | Refer to Appendix G |
| Cal Error - 16         | 2      | Code        | Refer to Appendix G |
| Span - 1               | 9      | Real        | 5 decimal places    |
| Span - 2               | 9      | Real        | 5 decimal places    |
| Span - 3               | 9      | Real        | 5 decimal places    |
| Span - 5               | 9      | Real        | 5 decimal places    |
| Span - 6               | 9      | Real        | 5 decimal places    |
| Rate - 1               | 9      | Real        | 5 decimal places    |
| Rate - 2               | 9      | Real        | 5 decimal places    |
| Rate - 3               | 9      | Real        | 5 decimal places    |
| Rate - 5               | 9      | Real        | 5 decimal places    |
| Rate - 6               | 9      | Real        | 5 decimal places    |
| Conc_Hi                | 9      | Real        | Refer to Note 2     |
| Iteration Number       | 5      | Numeric     |                     |
| Iteration Tolerance    | 9      | Real        | Refer to Note 2     |
| Standard Deviation     | 9      | Real        | Refer to Note 2     |
| Rate Deviation - 1     | 9      | Real        | Refer to Note 2     |
| Rate Deviation - 2     | 9      | Real        | Refer to Note 2     |
| Rate Deviation - 3     | 9      | Real        | Refer to Note 2     |
| Rate Deviation - 4     | 9      | Real        | Refer to Note 2     |
| Rate Deviation - 5     | 9      | Real        | Refer to Note 2     |
| Rate Deviation - 6     | 9      | Real        | Refer to Note 2     |
| Recovery Deviation - 1 | 9      | Real        | Refer to Note 2     |
| Recovery Deviation - 2 | 9      | Real        | Refer to Note 2     |
| Recovery Deviation - 3 | 9      | Real        | Refer to Note 2     |
| Recovery Deviation - 4 | 9      | Real        | Refer to Note 2     |
| Recovery Deviation - 5 | 9      | Real        | Refer to Note 2     |
| Recovery Deviation - 6 | 9      | Real        | Refer to Note 2     |
| Model                  | 2      | 1 - 5, 8, 9 | model number        |
| RO                     | 9      | Real        | Refer to Note 2     |
| Kc                     | 9      | Real        | Refer to Note 2     |
| a                      | 9      | Real        | Refer to Note 2     |
| b                      | 9      | Real        | Refer to Note 2     |
| c                      | 9      | Real        | Refer to Note 2     |

2 of 3



Table 4.22 Stream 702 - Function 83 Expanded Multipoint Calibration, continued

| Field  | Length | Format       | Description                              |
|--|--------|--------------|--|
| d  | 9      | Real         | Refer to Note 2                          |
| Slope Adjustment   | 10     | Real         | 5 decimal places<br>***** = does not fit |
| Offset Adjustment  | 10     | Real         | 5 decimal places                         |
| Future Use Space   | 30     | Alphanumeric |  |
| <p><b>EXAMPLE:</b> [-0,702,83,060391,103825,MPT6,###,#####,------26,<br/> 0.0000000,1.0000000,2.0000000,4.0000000,8.0000000,<br/> 16.000000,0,0,0,1,0,00,00,00,00,00,00,00,00,00,00,<br/> 00,00,00,00,00,00,00,--0.42980,--0.10018,--0.08096,<br/> --0.05396,--0.04705,--0.71195,881.13220,451.33038,<br/> 351.15433,270.19046,216.23090,169.18050,16.000000,<br/> ---9,27.955742,0.0023646,-0.000037,0.0009348,<br/> -0.001393,-0.001537,0.0042340,-0.002202,0.0000020,<br/> 0.9937373,2.0214555,4.0612559,7.5468578,16.750963,<br/> -1,105.32443,775.84485,-0.221832,-0.776385,0.0000000,<br/> 0.0000000,--1.00000,-0.00000,<br/> #####] ED&lt;CR&gt;&lt;LF&gt;</p> <p><b>(-Indicates space holder)</b></p> <p>Note 1: The number of decimal places in this field is the number in the selected precision for the default units field in the chem database plus an additional decimal place.</p> <p>Note 2: The number of decimal places is adjusted so that the number can fit in this field.</p> |        |              |  |

3 of 3

## 4.6 Stream 703 - Instrument Status

### Stream 703 - Function 1 Power Up

SENT BY: SYNCHRON CX to host

PURPOSE: This message is sent whenever the SYNCHRON CX goes through power up or is reset. This message can be turned off on the host setup screen.

Table 4.23 Stream 703 - Function 1 Power Up

| Field   | Length | Format | Description       |
|---|--------|--------|-------------------|
| Device ID   | 2      | 0 - 99 | Operator assigned |
| Stream  | 3      | 703    |                   |
| Function  | 2      | 01     |                   |
| Date  | 6      | ddmmyy |                   |
| Time  | 6      | hhmmss |                   |
| S.W. REVISIONS from versions files<br>master version number | 10     | ALPHA  |                   |
| CX4 zos version   | 10     | ALPHA  |                   |
| CX4 cpu boot prom   | 10     | ALPHA  |                   |
| CX4 host version  | 10     | ALPHA  |                   |
| CX4 host diag version                                       | 10     | ALPHA  |                   |
| CX4 barcode   | 10     | ALPHA  |                   |
| CX4 mac version   | 10     | ALPHA  |                   |
| CX4 mac prom  | 10     | ALPHA  |                   |
| CX4 mac diag  | 10     | ALPHA  |                   |
| CX4 mmc version   | 10     | ALPHA  |                   |
| CX4 mmc prom  | 10     | ALPHA  |                   |
| CX4 mmc seqr  | 10     | ALPHA  |                   |
| CX4 mmc seqr diag   | 10     | ALPHA  |                   |
| CX4 mmc seqr tbl  | 10     | ALPHA  |                   |
| CX4 mmc dev func  | 10     | ALPHA  |                   |
| CX4 mmc motor func  | 10     | ALPHA  |                   |
| CX4 msc slave   | 10     | ALPHA  |                   |
| CX4 msc prom  | 10     | ALPHA  |                   |
| CX4 msc diag  | 10     | ALPHA  |                   |
| CX3/ISE zos version   | 10     | ALPHA  |                   |
| CX3/ISE cpu boot prom                                       | 10     | ALPHA  |                   |
| CX3/ISE host version  | 10     | ALPHA  |                   |
| CX3/ISE host diag version                                   | 10     | ALPHA  |                   |
| CX3/ISE mmc version   | 10     | ALPHA  |                   |
| CX3/ISE mmc prom  | 10     | ALPHA  |                   |

1 of 2

Table 4.23 Stream 703 - Function 1 Power Up, continued

| Field  | Length | Format | Description |
|--|--------|--------|-------------|
| PS2 QNX version  | 10     | ALPHA  |             |
| PS2 CPU BOOT PROM  | 10     | ALPHA  |             |
| PS2 CLOCK TASK   | 10     | ALPHA  |             |
| PS2 DBM TASK   | 10     | ALPHA  |             |
| PS2 DISPMGR TASK   | 10     | ALPHA  |             |
| PS2 EMULATE TASK   | 10     | ALPHA  |             |
| PS2 EVENT TASK   | 10     | ALPHA  |             |
| PS2 KEYDEV TASK  | 10     | ALPHA  |             |
| PS2 LINK_CRT TASK  | 10     | ALPHA  |             |
| PS2 MASTER_SCR TASK  | 10     | ALPHA  |             |
| PS2 NOTE TASK  | 10     | ALPHA  |             |
| PS2 REC_ISR TASK   | 10     | ALPHA  |             |
| PS2 REC_TASK   | 10     | ALPHA  |             |
| PS2 SAMP_PROG TASK   | 10     | ALPHA  |             |
| PS2 SCSITASK   | 10     | ALPHA  |             |
| PS2 SENDTASK   | 10     | ALPHA  |             |
| PS2 SETUP TASK   | 10     | ALPHA  |             |
| PS2 CXSERVER TASK  | 10     | ALPHA  |             |
| PS2 DOWNLOAD TASK  | 10     | ALPHA  |             |
| PS2 CONSOLE TASK   | 10     | ALPHA  |             |
| PS2 SHUTDOWN TASK  | 10     | ALPHA  |             |
| PS2 ALARM TASK   | 10     | ALPHA  |             |
| PS2 QUEUE TASK   | 10     | ALPHA  |             |
| PS2 HCP TASK   | 10     | ALPHA  |             |
| PS2 RESULT TASK  | 10     | ALPHA  |             |
| PS2 PRINTER TASK   | 10     | ALPHA  |             |
| <p><b>EXAMPLE:</b> [-0,703,01,280291,115648,v0.0.04.07,v0.0.05.00,<br/> v1.0.11.00,v0.0.04.07,v0.0.04.07,v0.0.02.00,v0.0.02.00,<br/> v0.0.02.00,v0.0.02.00,v0.0.02.00,v0.0.02.00,v0.0.03.00,<br/> v0.0.02.00,v0.0.03.00,v0.0.02.00,v0.0.02.00,v0.0.02.00,<br/> v0.0.02.00,v0.0.02.00,v0.0.05.00,v1.0.11.00,v0.0.04.07,<br/> v0.0.04.07,v0.0.02.00,v0.0.02.00,v0.2.15.00,v0.0.01.01,<br/> v0.0.04.07,v0.0.04.07,v0.0.04.07,v0.0.04.07,v0.0.04.07,<br/> v0.0.04.07,v0.0.04.07,v0.0.04.07,v0.0.04.07,v0.0.04.07,<br/> v0.0.04.07,v0.0.04.07,v0.0.04.07,v0.0.04.07,v0.0.04.07,<br/> v0.0.04.07,v0.0.04.07,v0.0.04.07,v0.0.04.07,v0.0.04.07,<br/> v0.0.04.07,v0.0.04.07,v0.0.04.07,v0.0.04.07] 3F&lt;CR&gt;&lt;LF&gt;</p> <p><b>(-Indicates space holder)</b></p> |        |        |             |

2 of 2

**Stream 703 - Function 2 Bidirectional On**

SENT BY: SYNCHRON CX

PURPOSE: Notify the host when the SYNCHRON CX bidirectional interface option is enabled.

Table 4.24 Stream 703 - Function 2 Bidirectional On

| Field   | Length | Format       | Description                                     |
|---|--------|--------------|---|
| Device ID   | 2      | 0 - 99       | Operator assigned                               |
| Stream  | 3      | 703          | Always 703                                      |
| Function  | 2      | 02           | Always 02                                       |
| Date  | 6      | ddmmyy       | Day Month Year<br>(e.g. 240795 = July 24, 1995) |
| Time  | 6      | hhmmss       | Hour Minute Second                              |
| Software Set  | 10     | Alphanumeric | PGM: SOFT. VER                                  |
| EXAMPLE:     [-0,703,02,270291,163322,v0.0.04.07]8F<CR><LF> |        |              |   |
| (-Indicates space holder)                                   |        |              |   |

### Stream 703 - Function 3 Request Instrument State

SENT BY: Host

PURPOSE: Request the SYNCHRON CX transmit its instrument status information.

Table 4.25 Stream 703 - Function 3 Request Instrument State

| Field                          | Length | Format | Description       |
|--------------------------------|--------|--------|-------------------|
| Device ID                      | 2      | 0 - 99 | Operator assigned |
| Stream                         | 3      | 703    | Always 703        |
| Function                       | 2      | 03     | Always 03         |
| EXAMPLE: [00,703,03]93<CR><LF> |        |        |                   |

### Stream 703 - Function 4 Instrument State

SENT BY: SYNCHRON CX

PURPOSE: Transmission of SYNCHRON CX instrument status information.

Table 4.26 Stream 703 - Function 4 Instrument State

| Field  | Length | Format | Description                                     |
|--|--------|--------|---|
| Device ID  | 2      | 0 - 99 | Operator assigned                               |
| Stream   | 3      | 703    | Always 703                                      |
| Function   | 2      | 04     | Always 04                                       |
| Date   | 6      | ddmmyy | Day Month Year<br>(e.g. 240795 = 24 July 1995)  |
| Time   | 6      | hhmmss | Hour Minute Second                              |
| Instrument Configuration   | 5      | CXNAA  | SYNCHRON CX instrument configuration (e.g. CX7) |
| CX3/ISE State  | 2      | Code   | Refer to Appendix I                             |
| CX4 State  | 2      | Code   | Refer to Appendix I                             |
| Future Use Space   | 2      |        |   |
| Instrument Mode  | 2      | Code   | -0 = cup and sector mode<br>-1 = barcode mode   |
| EXAMPLE: [-0,703,04,270291,163324,CX7--, -2,15,##, -1]A5<CR><LF> |        |        |   |
| (-Indicates space holder)  |        |        |   |

### Stream 703 - Function 5 Instrument Exception

SENT BY: SYNCHRON CX to host

PURPOSE: This message is sent by the SYNCHRON CX whenever it encounters an exception condition. This corresponds to the same messages seen on the operator terminal.

Table 4.27 Stream 703 - Function 5 Instrument Exception

| Field        | Length | Format  | Description                                    |
|--------------|--------|---------|--|
| Device ID    | 2      | 0 - 99  | Operator assigned                              |
| Stream       | 3      | 703     | Always 703                                     |
| Function     | 2      | 05      | Always 05                                      |
| Date         | 6      | ddmmyy  | Day Month Year<br>(e.g. 240795 = 24 July 1995) |
| Time         | 6      | hhmmss  | Hour Minute Second                             |
| Error Number | 5      | Numeric | This is the error number.                      |
| Data         | 10     | Alpha   | Optional data associated with the error        |

### Stream 703 - Function 7 Chemistry Configuration Change

SENT BY: SYNCHRON CX to host

PURPOSE: Notify the host when the installed chemistry list has been changed at the SYNCHRON CX. This message is transmitted when the SYNCHRON CX Installed Chemistries function is exited. This message can be turned off on the Host Setup screen.

Table 4.28 Stream 703 - Function 7 Chemistry Configuration Change

| Field  | Length | Format | Description                                     |
|--|--------|--------|---|
| Device ID                                    | 2      | 0 - 99 | Operator assigned                               |
| Stream                                       | 3      | 703    | Always 703                                      |
| Function                                     | 2      | 07     | Always 07                                       |
| Date   | 6      | ddmmyy | Day Month Year<br>(e.g. 240795 = July 24, 1995) |
| Time   | 6      | hhmmss | Hour Minute Second                              |
| EXAMPLE: [-0,703,07,280291,100914]E2<CR><LF> |        |        |   |
| (-Indicates space holder)                    |        |        |   |

### Stream 703 - Function 13 Range Change

SENT BY: SYNCHRON CX

PURPOSE: Notify the host when the defined normal and/or critical ranges have been changed at the SYNCHRON CX. This message is transmitted when the SYNCHRON CX Define Ranges function is exited. This message can be turned off on the Host Setup screen.

Table 4.29 Stream 703 - Function 13 Range Change

| Field   | Length | Format    | Description                                     |
|---|--------|-----------|---|
| Device ID   | 2      | 0 - 99    | Operator assigned                               |
| Stream  | 3      | 703       | Always 703                                      |
| Function  | 2      | 13        | Always 13                                       |
| Date  | 6      | ddmmyy    | Day Month Year<br>(e.g. 240795 = July 24, 1995) |
| Time  | 6      | hhmmss    | Hour Minute Second                              |
| Comm Chem ID                                      | 4      | Chem Code | Refer to Appendix C; 0000 = all                 |
| EXAMPLE: [-0,703,13,280291,101724,TH01]9D<CR><LF> |        |           |   |
| (-Indicates space holder)                         |        |           |   |

### Stream 703 - Function 17 End of Run

SENT BY: SYNCHRON CX

PURPOSE: Notify the host when all programmed tests have been completed. This message can be turned off on the Host Setup screen.

Table 4.30 Stream 703 - Function 17 End of Run

| Field  | Length | Format | Description                                     |
|--|--------|--------|---|
| Device ID                                    | 2      | 0 - 99 | Operator assigned                               |
| Stream                                       | 3      | 703    | Always 703                                      |
| Function                                     | 2      | 17     | Always 17                                       |
| Date Queued                                  | 6      | ddmmyy | Day Month Year<br>(e.g. 240795 = July 24, 1995) |
| Time Queued                                  | 6      | hhmmss | Hour Minute Second                              |
| EXAMPLE: [-0,703,17,190291,103626]DE<CR><LF> |        |        |   |
| (-Indicates space holder)                    |        |        |   |

#### 4.7 Stream 704 - Setup Status

##### Stream 704 - Function 7 Request Installed Chemistries

SENT BY: Host

PURPOSE: Request transmission of the list of installed chemistries from the SYNCHRON CX.

Table 4.31 Stream 704 - Function 7 Request Installed Chemistries

| Field                          | Length | Format | Description       |
|--------------------------------|--------|--------|-------------------|
| Device ID                      | 2      | 0 - 99 | Operator assigned |
| Stream                         | 3      | 704    | Always 704        |
| Function                       | 2      | 07     | Always 07         |
| EXAMPLE: [00,704,07]8E<CR><LF> |        |        |                   |

##### Stream 704 - Function 8 Installed Chemistries

SENT BY: SYNCHRON CX

PURPOSE: Transmit the list of installed chemistries on the SYNCHRON CX.

Table 4.32 Stream 704 - Function 8 Installed Chemistries

| Field     | Length | Format | Description  |
|-----------|--------|--------|--|
| Device ID | 2      | 0 - 99 | Operator assigned  |
| Stream    | 3      | 704    | Always 704   |
| Function  | 2      | 08     | Always 08  |
| Date      | 6      | ddmmyy | Day Month Year<br>(e.g. 240795 = July 24, 1995)  |
| Time      | 6      | hhmmss | Hour Minute Second   |
| Chem ID 1 | 4      | AAAA   | Refer to Appendix C. Up to 999 chemistry codes may be transmitted.   |
| Chem ID 2 | 4      | AAAA   |  |
| :         | :      | :      | The last chemistry code transmitted is '0000', indicating no more codes will follow. If there are no installed chemistries there will be one zero field. |
| Chem ID N | 4      | AAAA   |  |

1 of 2



Table 4.32 Stream 704 - Function 8 Installed Chemistries, continued

| Field  | Length | Format | Description |
|--|--------|--------|-------------|
| <b>EXAMPLE:</b> [-0,704,08,280291,103941,MIGA,MDLC,CAL1,CAL2,CAL3,CAL4,THO1,74A-,59A-,40A-,75A-,46A-,76A-,77A-,57A-,78A-,48A-,55A-,82A-,90A-,79A-,83A-,91A-,89A-,08A-,35A-,31A-,10A-,30A-,12A-,11A-,05A-,09A-,44A-,42B-,42C-,32A-,03A-,36A-,06A-,33A-,34A-,43A-,07A-,42A-,41A-,62A-,72A-,67A-,68A-,69A-,70A-,52A-,51A-,53A-,71A-,04A-,02A-,01B-,01A-,05C-,06C-,03C-,09C-,0000]<br>C6<CR><LF> |        |        |             |
| <b>(-Indicates space holder)</b>   |        |        |             |

2 of 2

**Stream 704 - Function 9 Request Chemistry Ranges**

SENT BY:    Host

PURPOSE:    Request transmission of the defined normal and critical ranges from the SYNCHRON CX.

Table 4.33 Stream 704 - Function 9 Request Chemistry Ranges

| Field                                      | Length | Format    | Description                              |
|--|--------|-----------|--|
| Device ID                                  | 2      | 0 - 99    | Operator assigned                        |
| Stream                                     | 3      | 704       | Always 704                               |
| Function                                   | 2      | 09        | Always 09                                |
| Chem ID                                    | 4      | Chem Code | Refer to Appendix C Chem Code 0000 = all |
| <b>EXAMPLE:</b> [00,704,09,THO1]44<CR><LF> |        |           |  |

### Stream 704 - Function 10 Chemistry Ranges

SENT BY: SYNCHRON CX

PURPOSE: Transmit the defined normal and critical ranges for all installed chemistries on the SYNCHRON CX.

Table 4.34 Stream 704 - Function 10 Chemistry Ranges

| Field             | Length | Format  | Description  |
|-------------------|--------|---------|--|
| Device ID         | 2      | 0 - 99  | Operator assigned  |
| Stream            | 3      | 704     | Always 704   |
| Function          | 2      | 10      | Always 10  |
| Date              | 6      | ddmmyy  | Day Month Year<br>(e.g. 240795 = July 24, 1995)          |
| Time              | 6      | hhmmss  | Hour Minute Second                                       |
| Last Message      | 1      | Flag    | 0 = last range message<br>1 = more range messages follow |
| COM CHEM ID       | 4      | CHEM    | zero if none, or refer to Appendix C                     |
| Fluid Type        | 2      | FLUID   | Refer to Appendix B                                      |
| Number of Ranges  | 3      | Numeric | 1-10, times 3 (up to 30)                                 |
| Selected Units    | 2      | Code    | Refer to Appendix D                                      |
| Sex               | 1      | Alpha   | M = Male, F = Female, B = Both                           |
| Age Units         | 1      | Numeric | Refer to Appendix K                                      |
| Age Low           | 3      | Real    |  |
| Age Units         | 1      | Number  | Refer to Appendix K                                      |
| Age High          | 3      | Real    |  |
| Age Normal Low    | 8      | Real    |  |
| Age Normal High   | 8      | Real    |  |
| Age Critical Low  | 8      | Real    |  |
| Age Critical High | 8      | Real    |  |
| Sex               | 1      | Alpha   | M = Male, F = Female, B = Both                           |
| Age Units         | 1      | Numeric | Refer to Appendix K                                      |
| Age Low           | 3      | Real    |  |
| Age Units         | 1      | Number  | Refer to Appendix K                                      |
| Age High          | 3      | Real    |  |
| Age Normal Low    | 8      | Real    |  |
| Age Normal High   | 8      | Real    |  |
| Age Critical Low  | 8      | Real    |  |
| Age Critical High | 8      | Real    |  |
| Sex               | 1      | Alpha   | M = Male, F = Female, B = Both                           |
| Age Units         | 1      | Numeric | Refer to Appendix K                                      |

1 of 2

Table 4.34 Stream 704 - Function 10 Chemistry Ranges, continued

| Field   | Length | Format  | Description                    |
|---|--------|---------|--------------------------------|
| Age High  | 3      | Real    |                                |
| Age Normal Low  | 8      | Real    |                                |
| Age Normal High   | 8      | Real    |                                |
| Age Critical Low  | 8      | Real    |                                |
| Age Critical High   | 8      | Real    |                                |
| :   | :      | :       |                                |
| Sex   | 1      | Alpha   | M = Male, F = Female, B = Both |
| Age Units   | 1      | Numeric | Refer to Appendix K            |
| Age High  | 3      | Real    |                                |
| Age Normal Low  | 8      | Real    |                                |
| Age Normal High   | 8      | Real    |                                |
| Age Critical Low  | 8      | Real    |                                |
| Age Critical High   | 8      | Real    |                                |
| <p>EXAMPLE: [-0,704,10,280291,103159,0,THO1,SE,--4,-0,B,2,--1,2,-10,-----12,-----15,-----11,-----16,B,3,-10,3,-15,-----30,-----35,-----29,-----37,B,4,-10,4,-12,-----35,-----40,-----30,-----41,B,5,-10,5,-81,-----50,-----60,-----49,-----65]1B&lt;CR&gt;&lt;LF&gt;</p> <p>(-Indicates space holder)</p> |        |         |                                |

2 of 2

## SECTION 5 Results and Sample Programming Sequence:

### 5.1 Sample Programming

For each sample program, the host needs to bid for the line and send 701-1. SYNCHRON CX will respond with program ok, 701-2 message. If the status return is different from "-0" then sample is not programmed.

| SYNCHRON CX | HOST  |
|-------------|---|
| <ACK>       | <EOT><SOH>  |
| <ETX>       | [00,701,01,51,01,0,RO,SE,A12345678.1,...]CS<CR><LF> |
| <EOT><SOH>  | <EOT>   |
|             | <ACK>   |
|             | [00,701,02,00,12345,51,01,A12345678.1..]CS<CR><LF>  |
| <EOT>       | <ETX>   |

### 5.2 Host Query and Sample Programming

Each time the SYNCHRON CX instrument loads a sector, it checks if those samples are programmed. If samples are not programmed and SYNCHRON CX is in barcode mode it will send a host query message, if that message is "turned on". The host has a period of time to respond with the first sample program for that sector. For the CX DELTA Systems this period of time is selected by the operator in the host setup screen. Options are OFF, 2.5, 5, 7.5, and 10 minutes. Other systems are set to 2.5 minutes. Host query will have 1 to 7 sample ids.

| SYNCHRON CX | HOST                                  |
|-------------|---------------------------------------|
| <EOT><SOH>  | <ACK>                                 |
|             | <ETX>                                 |
| <EOT>       | <EOT><SOH>                            |
| <ACK>       |                                       |
| <ETX>       |                                       |
| <EOT><LB>   | <EOT>                                 |
|             | <ACK>                                 |
|             | [00,701,02,00,12345678.1..]CS<CR><LF> |
| <EOT>       | <ETX>                                 |

### 5.3 Results

When a sample is completed, SYNCHRON CX will send the sample result collated, meaning no other messages will be sent when one set of result is being sent. Special calculation and timed urine results are sent just before end of cup.

#### SYNCHRON CX

#### HOST

|            |   |       |
|------------|---|-------|
| <EOT><SOH> | _____>  |       |
|            | <_____  | <ACK> |
|            | [-0,702,01,220191,123033,RG,51,-1, ...]CS<CR><LF>_____> |       |
|            | <_____  | <ETX> |
|            | [-0,702,03,220191,124013,12345,.....]CS<CR><LF> _____>  |       |
|            | <_____  | <ACK> |
|            | [-0,702,81,220191,124013,12345,.....]CS<CR><LF> _____>  |       |
|            | <_____  | <ETX> |
|            | [-0,702,03,220191,124013,12345,.....]CS<CR><LF> _____>  |       |
|            | <_____  | <ACK> |
|            | [-0,702,81,220191,124013,12345,.....]CS<CR><LF> _____>  |       |
|            | <_____  | <ETX> |
|            | [-0,702,11,220191,124059,12345,.....]CS<CR><LF> _____>  |       |
|            | <_____  | <ACK> |
|            | [-0,702,11,220191,124059,12345,.....]CS<CR><LF> _____>  |       |
|            | <_____  | <ETX> |
|            | [-0,702,13,220191,124059,12345,.....]CS<CR><LF> _____>  |       |
|            | <_____  | <ACK> |
|            | [-0,702,05,220191,124013,12345,.....]CS<CR><LF> _____>  |       |
|            | <_____  | <ETX> |
| <EOT>      | _____>  |       |

## 5.4 Results Recalled When Running

When the operator recalls results for a sample which are not complete, the SYNCHRON CX will send all the tests which are complete. When sample is completed, it will send the rest of the results. Complete special calculation and timed urine will be sent in both sets.

| SYNCHRON CX   | HOST  |
|---|-------|
| <EOT><SOH>  |       |
| <_____>   |       |
| <_____>   | <ACK> |
| [-0,702,01,220191,123033,RE,51,-1,..]CS<CR><LF>_____>       |       |
| <_____>   | <ETX> |
| [-0,702,03,220191,124013,12345,..]CS<CR><LF> (chem1) _____> |       |
| <_____>   | <ACK> |
| [-0,702,81,220191,124013,12345,..]CS<CR><LF> (chem1) _____> |       |
| <_____>   | <ETX> |
| [-0,702,11,220191,124059,12345,..]CS<CR><LF> _____>         |       |
| <_____>   | <ACK> |
| [-0,702,11,220191,124059,12345,..]CS<CR><LF> _____>         |       |
| <_____>   | <ETX> |
| [-0,702,13,220191,124059,12345,..]CS<CR><LF> _____>         |       |
| <_____>   | <ACK> |
| [-0,702,05,220191,124013,12345,..]CS<CR><LF> _____>         |       |
| <_____>   | <ETX> |
| <EOT>   |       |
| <EOT><SOH>  |       |
| <_____>   |       |
| <_____>   | <ACK> |
| [-0,702,01,220191,123033,RG,51,-1,..]CS<CR><LF>_____>       |       |
| <_____>   | <ETX> |
| [-0,702,03,220191,125113,12345,..]CS<CR><LF> (chem2) _____> |       |
| <_____>   | <ACK> |
| [-0,702,81,220191,125113,12345,..]CS<CR><LF> (chem2) _____> |       |
| <_____>   | <ETX> |
| [-0,702,11,220191,125113,12345,..]CS<CR><LF> _____>         |       |
| <_____>   | <ACK> |
| [-0,702,11,220191,125113,12345,..]CS<CR><LF> _____>         |       |
| <_____>   | <ETX> |
| [-0,702,13,220191,125113,12345,..]CS<CR><LF> _____>         |       |
| <_____>   | <ACK> |
| [-0,702,05,220191,125113,12345,..]CS<CR><LF> _____>         |       |
| <_____>   | <ETX> |
| <EOT>   |       |

## 5.5 Options for Sending CX3 Results (on CX7) And ISE Results (on CX5)

When a sample is completed, and it has CX3/ISE chemistries, the user can opt to send CX3/ISE chemistries three different ways. Users can have the system send results only when the sample is finished, or send STAT CX3/ISE results as soon as available, or send any CX3/ISE result as soon as available. If the user elects to send STAT CX3/ISE results or any/all routine CX3/ISE results prior to sample completion, two sets of results will be sent to the host. The first set will include the CX3/ISE results (STAT or routine) and the second set the remaining results.

|            |  |       |
|------------|--|-------|
| <EOT><SOH> | _____>   |       |
|            | <_____>  | <ACK> |
|            | [-0,702,01,220191,123033,RG,51,-1, ...]CS<CR><LF>_____>      |       |
|            | <_____>  | <ETX> |
|            | [-0,702,03,220191,124013,12345,.....]CS<CR><LF> (CX3) _____> |       |
|            | <_____>  | <ACK> |
|            | [-0,702,81,220191,124013,12345,.....]CS<CR><LF> (CX3) _____> |       |
|            | <_____>  | <ETX> |
|            | [-0,702,11,220191,124059,12345,.....]CS<CR><LF> _____>       |       |
|            | <_____>  | <ACK> |
|            | [-0,702,11,220191,124059,12345,.....]CS<CR><LF> _____>       |       |
|            | <_____>  | <ETX> |
|            | [-0,702,13,220191,124059,12345,.....]CS<CR><LF> _____>       |       |
|            | <_____>  | <ACK> |
|            | [-0,702,05,220191,124013,12345,.....]CS<CR><LF> _____>       |       |
|            | <_____>  | <ETX> |
| <EOT>      | _____>   |       |
| <EOT><SOH> | _____>   |       |
|            | <_____>  | <ACK> |
|            | [-0,702,-1,220191,123033,RG,51,-1, ...]CS<CR><LF> _____>     |       |
|            | <_____>  | <ETX> |
|            | [-0,702,03,220191,124013,12345,.....]CS<CR><LF> (CX4) _____> |       |
|            | <_____>  | <ACK> |
|            | [-0,702,81,220191,124013,12345,.....]CS<CR><LF> (CX4) _____> |       |
|            | <_____>  | <ETX> |
|            | [-0,702,11,220191,124059,12345,.....]CS<CR><LF> _____>       |       |
|            | <_____>  | <ACK> |
|            | [-0,702,11,220191,124059,12345,.....]CS<CR><LF> _____>       |       |
|            | <_____>  | <ETX> |
|            | [-0,702,13,220191,124059,12345,.....]CS<CR><LF> _____>       |       |
|            | <_____>  | <ACK> |
|            | [-0,702,05,220191,124013,12345,.....]CS<CR><LF> _____>       |       |
|            | <_____>  | <ETX> |
| <EOT>      | _____>   |       |

## 5.6 Reagent Pack

When a calibration is completed, SYNCHRON CX will send reagent result collated.

### SYNCHRON CX

### HOST

|   |       |
|---|-------|
| <EOT><SOH> _____ >  |       |
| < _____   | <ACK> |
| [-0,702,21,220191,123033,RG,51,-1, ...]CS<CR><LF> _____ > |       |
| < _____   | <ETX> |
| [-0,702,23,220191,124013,12345,...]CS<CR><LF> _____ >     |       |
| < _____   | <ACK> |
| [-0,702,25,220191,124013,12345,...]CS<CR><LF> _____ >     |       |
| < _____   | <ETX> |
| [-0,702,23,220191,124059,12345,...]CS<CR><LF> _____ >     |       |
| < _____   | <ACK> |
| [-0,702,25,220191,124059,12345,...]CS<CR><LF> _____ >     |       |
| < _____   | <ETX> |
| [-0,702,07,220191,124059,12345,...]CS<CR><LF> _____ >     |       |
| < _____   | <ACK> |
| or  |       |
| [-0,702,09,220191,124059,12345,...]CS<CR><LF> _____ >     |       |
| < _____   | <ETX> |
| [-0,702,83,220191,124059,12345,...]CS<CR><LF> _____ >     |       |
| < _____   | <ACK> |
| [-0,702,27,220191,124059,12345,...]CS<CR><LF> _____ >     |       |
| < _____   | <ETX> |
| <EOT> _____ >   |       |



## SECTION 6 Operator Interface - Setting Host Communications Parameters

This feature allows the operator to establish compatibility between the instrument and a host computer for data transmission. Host Communications Parameters may be viewed at any time, but modifications can be made only when the system is in Standby.

### NOTICE

The console displays a message that Host Communications is going to be disabled when **IDLE** is pressed (for Backup, Restore, Resume, or Rebuild Databases), when performing a Chemistry Update, or when switching from barcode mode to sector mode and vice versa. The system will reestablish Host Communications following Backup, Resume, Rebuild Databases, Chemistry Update and switching between sector and barcode mode. However, the operator must toggle Host Communications to the desired mode following the Restore function.

### Defining Host Communications Parameters

### NOTICE

Screens are for example only and may not reflect your specific instrument configuration.

1. From the MASTER Screen, press **F4 SPECIAL FUNCTION**.
2. Cursor and **SELECT 4**. System Setup, or type **4 ENTER**.
3. Cursor and **SELECT 10**. Host Communications Parameters, or type **10 ENTER**.

SYSTEM SETUP  
CONSOLE: CX3: STANDBY Jul 14 95  
CX4: STANDBY 9:56

Cursor and <SELECT> or enter a number to choose item.

SYSTEM SETUP OPTIONS

|                             |                                    |
|-----------------------------|------------------------------------|
| 1. Configure Chemistry Menu | 8. Enable Automatic ORDAC          |
| 2. Define Chemistry Panels  | 9. Set Date/Time/Temperature       |
| 3. Bar Code Setup           | 10. Host Communications Parameters |
| 4. Define Reference Ranges  | 11. Replicates/Statistics          |
| 5. Report Setup             | 12. Define Comments                |
| 6. Special Calculations     | 13. Patient Result Setup           |
| 7. Units Selection          | 14. Define Reportable Ranges       |

Item Number: 10

F1 F2 F3 F4 F5 F6 F7 F8

A\_06422s.eps

4. Refer to this Host Computer Interface Specifications and also to Table 6.1 for comprehensive documentation of the host parameters. Move cursor to the desired parameter and press **SELECT** to scroll through available options or enter the appropriate values.

**NOTICE**  
Automated Sample Handling is for use with Optional Robotics interface.

```
DEFINE HOST          CX3: STANDBY          Jul 14 95
CONSOLE:             CX4: STANDBY          9:56

Press <SELECT> to set desired option

HOST COMMUNICATION PARAMETERS

Mode: [UNIDIRECTIONAL]
Baud Rate: [4800]
Data Bits: [8]
Stop Bits: [1]
Parity: [NONE]
Device ID: 0
Flow Control: [XON/XOFF]
Automated Sample Handling: [OFF]

MESSAGES  PRINT
F1        F2        F3        F4        F5        F6        F7        F8
A_06423s.eps
```

5. Message options can be accessed by pressing the F1 messages key. Move cursor to desired option and press **SELECT** key to toggle OFF/ON. Use PAGE key to access additional messages.

```
DEFINE HOST          CX3: STANDBY          Dec  4 90
CONSOLE:             CX4: STANDBY          13:25

Press <SELECT> to set desired option

UNIDIRECTIONAL MESSAGE OPTIONS

Stream 700 SPECIAL FUNCTIONS
Func 2 Host Setup [ON ]
Stream 702 RESULTS
Func 3 Test Results [ON ]
Func 7 Linear Calibration [ON ]
Func 9 Multipoint Calibration [ON ]
Func 11 Special Calculation [ON ]
Func 23 Results for Calib [ON ]
Func 25 Expanded Results for Calib [ON ]
Func 81 Expanded Results [OFF]
Func 83 Expanded Calibration [ON ]

Press <PREV SCREEN> to exit

F1        F2        F3        F4        F5        F6        F7        F8
A_06655s.eps
```

6. Press **PREV SCREEN** to return to the SYSTEM SETUP Screen, or press **MASTER SCREEN** to exit.

Table 6.1 Host Communication Parameter Setup Options

| Field   | Available Options                                      |
|---|--|
| MODE:   | UNIDIRECTIONAL, BIDIRECTIONAL, OFF                     |
| BAUD RATE:  | 600, 1200, 2400, 4800, 9600                            |
| DATA BITS:  | 7, 8   |
| STOP BITS:  | 1, 2   |
| PARITY:   | NONE, ODD, EVEN  |
| DEVICE ID:  | 0 - 99   |
| FLOW CONTROL:   | XON/XOFF, RTS/CTS, NONE                                |
| <b>Status Messages: (U)nidirectional, (B)idirectional</b> |  |
| Stream 700/Func 2:  | Host Setup ON, OFF (U,B)                               |
| Stream 701/Func 6:  | Host Query ON, OFF (B)<br>Auto Clear Queue ON, OFF (B) |
| Stream 702/Func 3:  | Test Results ON, OFF (U,B)                             |
| Stream 702/Func 7:  | Linear Calibration ON, OFF (U,B)                       |
| Stream 702/Func 9:  | Multipoint Calibration ON, OFF (U,B)                   |
| Stream 702/Func 11:                                       | Special Calculations ON, OFF (U,B)                     |
| Stream 702/Func 13:                                       | Timed Urine Result ON, OFF (U,B)                       |
| Stream 702/Func 23:                                       | Results for Calibration ON, OFF (U,B)                  |
| Stream 702/Func 25:                                       | Expanded Calib Result ON, OFF (U,B)                    |
| Stream 702/Func 81:                                       | Expanded Results ON, OFF (U,B)                         |
| Stream 702/Func 83:                                       | Expanded Calibration ON, OFF (U,B)                     |
| Stream 703/Func 1:  | Power Up ON, OFF (U,B)                                 |
| Stream 703/Func 5:  | Instrument Exception ON, OFF (U,B)                     |
| Stream 703/Func 7:  | Chem Change ON, OFF (U,B)                              |
| Stream 703/Func 13:                                       | Range Change ON, OFF (U,B)                             |
| Stream 703/Func 17:                                       | End of Run ON, OFF (U,B)                               |

**IMPORTANT NOTE**

Regarding Stream 702, Functions 3 and 81 - at least one of these functions must be set ON. Both can be ON, but both cannot be set to OFF at the same time.

## **SECTION 7 Appendices**

### **APPENDIX A – Test Type Codes**

CA\* = Calibrator  
CO = Control  
EX\* = Extinction  
RO = Routine  
SC = STAT Control  
ST = STAT

\*Cannot be programmed by the host computer.

### **APPENDIX B – Sample Type Codes**

SE = Serum  
SF = Spinal Fluid  
UR = Urine  
TU = Timed Urine  
PL = Plasma  
Future Use  
Future Use

## APPENDIX C – Chemistry Codes

|     |   |  |
|-----|---|--|
| 01A | = | Sodium (NA)  |
| 01B | = | Potassium (K)  |
| 02A | = | Carbon Dioxide (CO <sub>2</sub> )                        |
| 03A | = | Creatinine (CRE)   |
| 03C | = | CX3/CX3 DELTA Creatinine (CRE3/CREA)                     |
| 03D | = | Creatinine Triggered (CR-T)                              |
| 04A | = | Chloride (CL)  |
| 05A | = | Blood Urea Nitrogen (BUN)                                |
| 05C | = | CX3/CX3 DELTA Blood Urea Nitrogen (BUN3/BUN)             |
| 06A | = | Glucose (GLU)  |
| 06C | = | CX3/CX3 DELTA Glucose (GLU3/GLU)                         |
| 07A | = | Total Protein (TP)                                       |
| 07B | = | Micro Protein (M-TP)                                     |
| 07C | = | CX3/CX3 DELTA Total Protein (TP3/TP)                     |
| 08A | = | Albumin (ALB)  |
| 08E | = | Prealbumin (PAB)   |
| 08M | = | Microalbumin (MA)  |
| 09A | = | Calcium (CA)   |
| 09C | = | CX3/CX7 Cup Calcium (CA/CA3)                             |
| 09D | = | CX3 DELTA ISE/CX5/CX7 DELTA ISE (Ca/CALC)                |
| 10A | = | Amylase (AMY)  |
| 10B | = | Pancreatic Amylase (PAMY)                                |
| 11A | = | Total Bilirubin (TBIL)                                   |
| 12A | = | Direct Bilirubin (DBIL)                                  |
| 14A | = | Hemoglobin (HB)  |
| 14B | = | Hemoglobin A1c (A1c)                                     |
| 30A | = | Aspartate Aminotransferase (AST)                         |
| 30B | = | Aspartate Aminotransferase-Pyridoxal-5'-Phosphate (AST-) |
| 31A | = | Alanine Aminotransferase (ALT)                           |
| 31B | = | Alanine Aminotransferase-Pyridoxal-5'-Phosphate (ALT-)   |
| 32A | = | Creatine Kinase (CK)                                     |
| 32B | = | Creatine Kinase-N-Acetyl-L-Cysteine (CK-)                |
| 33A | = | Lactate Dehydrogenase (LD-L)                             |
| 34A | = | Lactate Dehydrogenase (LD-P)                             |
| 35A | = | Alkaline Phosphatase (ALP)                               |
| 36A | = | Gamma Glutamyltransferase (GGT)                          |
| 40A | = | Creatine Kinase MB (CKMB)                                |
| 41A | = | Uric Acid (URIC)   |
| 42A | = | Triglycerides (TRIG)                                     |
| 42B | = | Triglycerides GPO (TG)                                   |
| 42C | = | Triglycerides (TG-B)                                     |
| 43A | = | Inorganic Phosphorus (PHOS)                              |
| 43B | = | Inorganic Phosphorus (PO <sub>4</sub> )                  |
| 44A | = | Cholesterol (CHOL)                                       |
| 46A | = | Iron (IRON)  |
| 46B | = | Iron (FE)  |
| 48A | = | Magnesium (MG)   |
| 49A | = | Acid Phosphatase (ACP)                                   |
| 50A | = | Lipase (LIPA)  |
| 50B | = | Lipase Wash (LIWA)                                       |

1 of 3

## APPENDIX C – Chemistry Codes, continued

|      |   |  |
|------|---|--|
| 51A  | = | Immunoglobulin G (IGG)                             |
| 52A  | = | Immunoglobulin A (IGA)                             |
| 53A  | = | Immunoglobulin M (IGM)                             |
| 54A  | = | Ammonia (AMM)                                      |
| 55A  | = | Lactate (LAC)                                      |
| 56A  | = | Ethyl Alcohol (ETOH)                               |
| 56B  | = | Alcohol (ALC)                                      |
| 57A* | = | a-Hydroxybutyrate Dehydrogenase (HBDH)             |
| 59A* | = | Cholinesterase (CHE)                               |
| 62A  | = | Digoxin (DIG)                                      |
| 62B  | = | Digoxin (DIGN)                                     |
| 63A  | = | Methaqualone (METQ)                                |
| 64A  | = | Methadone (METD)                                   |
| 65A  | = | Propoxyphene (PROX)                                |
| 66A  | = | Phencyclidine (PCP)                                |
| 67A  | = | Phenobarbital (PHNB)                               |
| 67B  | = | Phenobarbital (PBR)                                |
| 67C  | = | Phenobarbital (PHE)                                |
| 68A  | = | Phenytoin (PHNY)                                   |
| 68C  | = | Phenytoin (PHY)                                    |
| 69A  | = | Theophylline (THEO)                                |
| 69C  | = | Theophylline (THE)                                 |
| 70A  | = | Tobramycin (TOBR)                                  |
| 70C  | = | Tobramycin (TOB)                                   |
| 71A  | = | Transferrin (TRF)                                  |
| 72A  | = | Gentamicin (GENT)                                  |
| 72C  | = | Gentamicin (GEN)                                   |
| 73A  | = | Urea Nitrogen (UREA)                               |
| 73C  | = | CX3/CX3 DELTA Urea Nitrogen (URE3/UREA)            |
| 74A* | = | Alkaline Phosphatase (ALPd)                        |
| 75A* | = | Creatine Kinase NAC Buffer (CKNa)                  |
| 76A* | = | Glutamate Oxalacetate Transaminase (GOT)           |
| 77A* | = | Glutamate Pyruvate Transaminase (GPT)              |
| 78A* | = | Lactate Dehydrogenase (LDH)                        |
| 79A  | = | Total Iron Binding Capacity (TIBC)                 |
| 79B  | = | Total Iron Binding Capacity (IBCT)                 |
| 82A  | = | Leucine Aminopeptidase (LAP)                       |
| 83A  | = | High Density Lipoprotein Cholesterol (HDLc)        |
| 83D  | = | High Density Lipoprotein Cholesterol Direct (HDLd) |
| 84A  | = | Amphetamines (AMPH)                                |
| 85A  | = | Barbiturates (BARB)                                |
| 86A  | = | Benzodiazepine (BENZ)                              |
| 87A  | = | Cocaine Metabolites (COCM)                         |
| 88A  | = | Cannabinoids - 100ng (THC)                         |
| 88B  | = | Cannabinoids - 20ng (THC2)                         |
| 88C  | = | Cannabinoids - 50ng (THC5)                         |
| 89A  | = | C Reactive Protein (CRP)                           |
| 89B  | = | C Reactive Protein (CRP-)                          |
| 90A  | = | T Uptake (TU)                                      |
| 91A  | = | Thyroxine (T4)                                     |
| 92A  | = | Opiates (OP)                                       |
| 92B  | = | Opiate 2000 ng (OP2)                               |

### APPENDIX C – Chemistry Codes, continued

|     |   |                           |
|-----|---|---------------------------|
| 93A | = | Antistreptolysin O (ASO)  |
| 93B | = | Antistreptolysin O (ASO-) |
| 93C | = | Rheumatoid Factor (RF)    |
| 94A | = | Salicylate (SAL)          |
| 94M | = | Acetaminophen (ACTM)      |
| 95A | = | Valproic Acid (VPA)       |
| 98A | = | Carbamazepine (CAR)       |
| 99B | = | Icterus (ICTER)           |
| 99C | = | Lipemia (LIPEM)           |
| 99D | = | Hemolysis (HEMOL)         |
| 99G | = | Diluent 1 (DIL1)          |

3 of 3

\* Deutsche Gesellschaft für Klinische Chemie (German Clinical Chemistry Association) Formulations

#### NOTICE

Chemistry codes for user-defined chemistries correspond to the testname defined for the chemistry on the SYNCHRON CX USER-DEFINED SETUP Screen.

#### APPENDIX D – Unit Codes

|      |        |      |        |
|------|--------|------|--------|
| 00 = | mg/dL  | 14 = | nmol/L |
| 01 = | mg/L   | 15 = | Ku.u.  |
| 02 = | g/dL   | 16 = | U/L    |
| 03 = | g/L    | 17 = | Other  |
| 04 = | mmol/L | 18 = | %      |
| 05 = | μmol/L | 19 = | mA     |
| 06 = | mEq/L  | 20 = | mA/min |
| 07 = | nKat/L | 21 = | IU/mL  |
| 08 = | μKat/L | 22 = | U/mL   |
| 09 = | IU/L   | 23 = | Rate   |
| 10 = | μg/mL  | 24 = | ng/dL  |
| 11 = | ng/mL  | 25 = | μIU/mL |
| 12 = | μg/dL  | 26 = | mIU/mL |
| 13 = | μg/L   | 27 = | KU/L   |

#### APPENDIX E – Instrument Codes

##### CX4CX5

|     |                           |
|-----|---------------------------|
| C = | Calibration Overridden    |
| D = | Days Exceeded             |
| E = | Calibration Time Extended |
| H = | Temperature Error         |
| J = | Slope/Offset Adjustment   |
| M = | Set Point Modification    |
| O = | ORDAC Sample              |
| R = | Reagent Expired           |
| S = | Service Mode              |
| T = | Temperature Correction    |
| V = | Rerun Sample              |
| Z = | Edited Result             |

##### CX3CX3 DELTA

|     |                        |
|-----|------------------------|
| C = | Calibration Overridden |
| E = | Cal Time Extended      |
| T = | Cal Time Modified      |
| B*  | Channel Bypassed       |
| P = | Creatinine Bubble      |
| H*  | DAC Offset             |
| A*  | Erratic ADC            |
| I = | GLU IC                 |
| G = | GLU membrane           |
| D*  | No Sample              |
| N*  | Not Run                |
| O = | ORDAC Sample           |
| R = | Reagent Strength       |
| X*  | Reference Drift        |
| S*  | Service                |
| F = | ORDAC to follow        |
| V = | Rerun Sample           |
| Z = | Edited Result          |

\* Also applies to ISE Module on CX5



## APPENDIX F – Result Error Codes

|  |   |
|--|---|
| AE = ADC Error                                     | LR = Reaction Absorbance Low            |
| AH = Initial Absorbance High                       | NO = No Result Error                    |
| BH = Blank Absorbance High                         | NT = Noise Threshold                    |
| BL = Blank Absorbance Low                          | OH = ORDAC High                         |
| BN = Blank Mean Deviation                          | OL = ORDAC Low                          |
| BO = Blank Maximum Deviation                       | OT = Outliers Threshold                 |
| DE = Out of Instrument Electronic Range            | RH = Reaction Rate High                 |
| DH = Out of Instrument Range High                  | RL = Reaction Rate Low                  |
| DL = Out of Instrument Range Low                   | RN = Reaction Mean Deviation            |
| DR = Reference Drift (ISE)                         | RO = Reaction Maximum Deviation         |
| EA = Erratic ADC (ISE)                             | SD = Substrate Depleted                 |
| HI = Out of Instrument Range High (CX3 DELTA ONLY) | SH = Blank Rate High                    |
| HR = Reaction Absorbance High                      | SL = Blank Rate Low                     |
| IR = Initial Rate High                             | TM = Temperature                        |
| IT = Iteration Tolerance                           | UO = Out of ORDAC Reportable Range High |
| LO = Out of Instrument Range Low (CX3 DELTA ONLY)  | UH = Out of Reportable Range High       |
|  | UL = Out of Reportable Range Low        |

## APPENDIX G – Calibration Error Codes

|                            |                               |
|----------------------------|-------------------------------|
| 00 = No Error              | 19 = Math Error 12            |
| 01 = Back-To-Back Error    | 20 = Math Error 13            |
| 02 = Blank Absorbance High | 21 = Math Error 14            |
| 03 = Blank Absorbance Low  | 22 = Math Error 15            |
| 04 = DAC                   | 23 = Math Error 16            |
| 05 = Erratic ADCs          | 25 = Range Error              |
| 06 = Inconsistent Data     | 26 = Calibrator Range Hi      |
| 08 = Math Error 1          | 27 = Calibrator Range Lo      |
| 09 = Math Error 2          | 28 = Reaction Absorbance Hi   |
| 10 = Math Error 3          | 29 = Reaction Absorbance Lo   |
| 12 = Math Error 5          | 30 = Span                     |
| 13 = Math Error 6          | 31 = Sensitivity Error        |
| 14 = Math Error 7          | 32 = Severe Sensitivity Error |
| 15 = Math Error 8          | 33 = Recovery Error           |
| 16 = Math Error 9          | 34 = Severe Recovery Error    |
| 17 = Math Error 10         | 35 = ADC Electrolyte Error    |
| 18 = Math Error 11         |                               |

## APPENDIX H – Special Calculation Status Codes

|   |
|---|
| AB = One of the involved chemistries was not run.           |
| OK = Valid result.  |
| UN = Units for the involved chemistries are not compatible. |
| ZD = Denominator of a ratio is zero.                        |

## APPENDIX I – Instrument Status Codes

| CX4 Instrument<br>States Available     | CX3/ISE Module<br>States Available |
|--|------------------------------------|
| 00 = No state                          | 00 = No state                      |
| 01 = Stopped                           | 01 = Stopped                       |
| 02 = Pause initiated                   | 02 = Standby                       |
| 04 = Extinction-Coefficient            | 03 = System home                   |
| 05 = Running                           | 04 = Reagent load                  |
| 06 = Idle shutdown in progress         | 05 = Prime                         |
| 07 = Idle                              | 06 = Calibration                   |
| 08 = Initializing                      | 07 = Running                       |
| 09 = Reagent load                      | 08 = Maintenance                   |
| 10 = Homing                            | 09 = Autoprime                     |
| 11 = Priming                           | 10 = Calibration request           |
| 13 = Saving to disk                    | 11 = No state                      |
| 14 = reading from disk                 | 12 = Boot                          |
| 15 = standby                           | 13 = Pause-initiated               |
| 16 = Checking levels                   | 14 = Waiting                       |
| 18 = procedure in progress             | 15 = System-idle                   |
| 19 = procedure termination in progress | 16 = Loading*                      |
| 20 = Procedure complete                | 17 = System Stopping*              |
| 21 = Waiting                           |                                    |

\* Synchron CX3 DELTA only

## APPENDIX J – Key Code Conversion for Local Languages

For communications between the host computer and the SYNCHRON CX, local language character handling is accomplished via the extended ASCII character set as defined by IBM for the IBM PC and IBM compatible units. The correct interpretation of foreign language characters on the SYNCHRON CX will require the use of 8-bit communication between the system and the host computer.

The defined foreign language characters and their extended ASCII character codes are as follows:

| DECIMAL VALUE | HEXA DECIMAL VALUE | 0            | 16            | 32 | 48 | 64 | 80 | 96 | 112 |
|---------------|--------------------|--------------|---------------|----|----|----|----|----|-----|
| 0             | 0                  | BLACK (NULL) | BLANK (SPACE) | 0  | @  | P  | ‘  | p  |     |
| 1             | 1                  | ☺            | ◀             | !  | 1  | A  | Q  | a  | q   |
| 2             | 2                  | ☹            | ↕             | "  | 2  | B  | R  | b  | r   |
| 3             | 3                  | ♥            | !!            | #  | 3  | C  | S  | c  | s   |
| 4             | 4                  | ♦            | ⌘             | \$ | 4  | D  | T  | d  | t   |
| 5             | 5                  | ♣            | §             | %  | 5  | E  | U  | e  | u   |
| 6             | 6                  | ♠            | ■             | &  | 6  | F  | V  | f  | v   |
| 7             | 7                  | •            | ↕             | '  | 7  | G  | W  | g  | w   |
| 8             | 8                  | •            | ↑             | (  | 8  | H  | X  | h  | x   |
| 9             | 9                  | ○            | ↓             | )  | 9  | I  | Y  | i  | y   |
| 10            | A                  | ○            | →             | *  | :  | J  | Z  | j  | z   |
| 11            | B                  | ♂            | ←             | +  | ;  | K  | l  | k  | {   |
| 12            | C                  | ♀            | └             | ,  | <  | L  | \  | l  |     |
| 13            | D                  | ♪            | ↔             | –  | =  | M  | l  | m  | }   |
| 14            | E                  | ♪            | ▲             | ·  | >  | N  | ^  | n  | ~   |
| 15            | F                  | ☼            | ▼             | /  | ?  | O  | —  | o  | △   |

| DECIMAL VALUE | HEXA DECIMAL VALUE | 128 | 144 | 160 | 176 | 192 | 208 | 224 | 240      |
|---------------|--------------------|-----|-----|-----|-----|-----|-----|-----|----------|
| 0             | 0                  | Ç   | É   | á   | ⋮   |     |     | α   | ≡        |
| 1             | 1                  | ü   | æ   | í   | ⋮   |     |     | β   | ±        |
| 2             | 2                  | é   | Æ   | ó   | ⋮   |     |     | Γ   | ≥        |
| 3             | 3                  | â   | ô   | ú   |     |     |     | π   | ≤        |
| 4             | 4                  | ä   | ö   | ñ   |     |     |     | Σ   | ∫        |
| 5             | 5                  | à   | ò   | Ñ   |     |     |     | σ   | ∫        |
| 6             | 6                  | å   | û   | ä   |     |     |     | μ   | ÷        |
| 7             | 7                  | ç   | ù   | o   |     |     |     | τ   | ≈        |
| 8             | 8                  | ê   | ÿ   | ¿   |     |     |     | ø   | ◊        |
| 9             | 9                  | ë   | Ö   | Γ   |     |     |     | Θ   | •        |
| 10            | A                  | è   | Ü   | ¬   |     |     |     | Ω   | •        |
| 11            | B                  | ï   | ç   | ½   |     |     |     | δ   | √        |
| 12            | C                  | î   | £   | ¼   |     |     |     | ∞   | n        |
| 13            | D                  | ì   | ¥   | ì   |     |     |     | φ   | ²        |
| 14            | E                  | Ä   | Ð   | «   |     |     |     | €   | ■        |
| 15            | F                  | Å   | f   | »   |     |     |     | ∩   | BLACK FF |

A\_00071C.EPS

No other extended ASCII characters are supported by the SYNCHRON CX system for display, printing and host communication.

## **APPENDIX K – Age Units**

|   |   |        |
|---|---|--------|
| 1 | = | Hours  |
| 2 | = | Days   |
| 3 | = | Weeks  |
| 4 | = | Months |
| 5 | = | Years  |